

SOLUTIONS FOR
SAMPLE PREPARATION



SILICYCLE 


Greyhound
CHROMATOGRAPHY
AND ALLIED CHEMICALS



Founded in 1995, SiliCycle is specialized in the development, manufacturing and commercialization of high value silica gels and specialty products for chromatography, purification and synthesis.



Solutions for Sample Preparation

SiliaPrep & SiliaPrepX SPE Cartridges and Well Plates	2
Ordering Information	8
Applications.....	18
SiliaPrep Accessories	40
SiliaPrep Tips Micro-SPE Cartridges.....	46
SiliaQuick & SiliaFast Sample Preparation & Pesticide Analysis	50
Contact Us.....	52





SiliaPrep™ and SiliaPrepX™

SPE Cartridges and Well Plates

- Wide variety of sorbents
- Tight particle size distribution
- Very good packing (*no fines*)
- High recovery and yield

Silica-based and polymeric sorbents

Solid-phase extraction (SPE) is designed for rapid sample preparation and purification prior to chromatographic analysis.










Our SiliaPrep (*silica-based*) and SiliaPrepX (*polymeric*) families of SPE cartridges and well plates have been created to cover the entire spectrum of solid-phase extraction. This complete range of sorbents allows treatment of most common matrices:

- human and animal biological fluids
- waste waters
- petrochemical residues
- toxicological residues
- food and beverage

SiliaPrep and SiliaPrepX products are made using state-of-the-art technology, giving you the highest quality and the best lot-to-lot reproducibility. These advanced sorbents are providing you a clean extract, reducing ion suppression and increasing selectivity for LC/MS/MS applications. All our ultra pure SiliaFlash silica gels and functionalized SiliaBond silica gels are available in SPE formats. Just tell us what you need!

Cartridge sizes

We can provide a complete range of SPE cartridge lengths and diameters.

SiliaPrep Cartridge Sizes								
Tips Micro-SPE Cartridges	Mini-SiliaPrep SPE Cartridges		SPE Cartridges					
	500 mg	1 g	1 mL	3 mL	6 mL	10 mL	12 mL	25 mL
								

Bigger sizes (70 mL, 150 mL & 276 mL) are also available under SiliaSep OT branding.

Tips for your method development

Tips for Your Method Development		
Sorbent Type	Silica-Based (SiliaPrep)	Polymeric (SiliaPrepX)
Sorbent Capacity	Load up to 5 % of bed weight: 100 mg of silica-based sorbent will retain up to 5 mg of sample	Load up to 10 % of bed weight: 100 mg of polymeric sorbent will retain up to 10 mg of sample

Not enough sorbent: ANALYTE LOSS ► low recovery and reproducibility
 Too much sorbent: MORE EXPENSIVE ► more solvent used, taller SPE cartridges
 Concentrated samples: double the bed weight to avoid analyte loss

Product Selection Guide by Technical Characteristics

Product Selection Guide by Technical Characteristics								
SiliaPrep / SiliaPrepX	Sorbent Number	Particulate Size	Pore Size	Surface Area	Carbon Load*	Endcapping	Ionic Capacity	pH Stability
Silica-Based Non Polar Phases								
SiliaPrep C18 Plus	SPE-R00830B-xxx	40 - 63 µm	60 Å	500 m ² /g	17 %	Proprietary	-	2 - 10
SiliaPrep C18 nec	SPE-R35530B-xxx	40 - 63 µm	60 Å	500 m ² /g	17 %	No	-	2 - 10
SiliaPrep C18 WPD	SPE-R33229G-xxx	37 - 55 µm	125 Å	300 m ² /g	13 %	Yes	-	2 - 10
SiliaPrep C8	SPE-R31030B-xxx	40 - 63 µm	60 Å	500 m ² /g	12 %	Yes	-	2 - 10
SiliaPrep C8 nec	SPE-R31130B-xxx	40 - 63 µm	60 Å	500 m ² /g	12 %	No	-	2 - 10
SiliaPrep Phenyl (PH)	SPE-R34030B-xxx	40 - 63 µm	60 Å	500 m ² /g	9 %	Yes	-	2 - 10
SiliaPrep PFP	SPE-R67530B-xxx	40 - 63 µm	60 Å	500 m ² /g	11 %	Yes	-	2 - 10
Silica-Based Polar Phases								
SiliaPrep Cyano (CN)	SPE-R38030B-xxx	40 - 63 µm	60 Å	500 m ² /g	7 %	Yes	-	2 - 10
SiliaPrep Diol nec	SPE-R35030B-xxx	40 - 63 µm	60 Å	500 m ² /g	8 %	No	-	2 - 10
SiliaPrep Florisil	SPE-AUT-0014-xxx	40 - 75 µm	100 Å	250 m ² /g	-	-	-	3 - 8
SiliaPrep Florisil LP	SPE-AUT-0014LP-xxx	75 - 150 µm	100 Å	250 m ² /g	-	-	-	3 - 8
SiliaPrep Florisil PR	SPE-AUT-0015-xxx	150 - 250 µm	100 Å	200 m ² /g	-	-	-	3 - 8
SiliaPrep Silica	SPE-R10030B-xxx	40 - 63 µm	60 Å	500 m ² /g	-	-	-	2 - 9
SiliaPrep Silica WPD	SPE-R10029G-xxx	37 - 55 µm	125 Å	300 m ² /g	-	-	-	2 - 9
SiliaPrep Acidic Alumina	SPE-AUT-0053-xxx	75 - 150 µm	70 Å	150 - 320 m ² /g	-	-	-	3 - 8
SiliaPrep Neutral Alumina	SPE-AUT-0054-xxx	75 - 150 µm	70 Å	150 - 320 m ² /g	-	-	-	3 - 8
SiliaPrep Basic Alumina	SPE-AUT-0055-xxx	75 - 150 µm	70 Å	150 - 320 m ² /g	-	-	-	3 - 8
Silica-Based Ion Exchange Phases								
SiliaPrep SAX nec	SPE-R66530B-xxx	40 - 63 µm	60 Å	500 m ² /g	10 %	No	0.90 meq/g	2 - 10
SiliaPrep SAX-2 nec	SPE-R66430B-xxx	40 - 63 µm	60 Å	500 m ² /g	9 %	No	0.71 mmol/g	2 - 10
SiliaPrep Carbonate	SPE-R66030B-xxx	40 - 63 µm	60 Å	500 m ² /g	6 %	Yes	0.46 mmol/g	2 - 10
SiliaPrep Amine (WAX)	SPE-R52030B-xxx	40 - 63 µm	60 Å	500 m ² /g	7 %	Yes	1.2 mmol/g	2 - 10
SiliaPrep SCX	SPE-R60530B-xxx	40 - 63 µm	60 Å	500 m ² /g	9 %	Yes	0.54 meq/g	2 - 10
SiliaPrep SCX-2	SPE-R51230B-xxx	40 - 63 µm	60 Å	500 m ² /g	5 %	Yes	0.63 meq/g	2 - 10
SiliaPrep WCX	SPE-R70030B-xxx	40 - 63 µm	60 Å	500 m ² /g	7 %	Yes	0.92 mmol/g	2 - 10
Specialty Phases								
SiliaPrep PCB	SP2-R00650030B-xxx	40 - 63 µm	60 Å	500 m ² /g	3 %	Proprietary	-	2 - 10
SiliaPrep CleanDRUG	SPEC-R651230B-xxx	40 - 63 µm	60 Å	500 m ² /g	9 %	Proprietary	-	2 - 10
SiliaPrep CleanENVI	SPEC-R31930B-xxx	40 - 63 µm	60 Å	500 m ² /g	19 %	Proprietary	-	2 - 10
SiliaPrep PAH	SP2-R0610030B-xxx	40 - 63 µm	60 Å	500 m ² /g	13 %	Proprietary	-	2 - 10
Polymeric Phases								
SiliaPrepX DVB	SPE-P0001-xxx	85 µm	60 Å	1,000 m ² /g	90 %	-	-	1 - 14
SiliaPrepX HLB	SPE-P0002-xxx	40 µm	110 Å	850 m ² /g	88 %	-	-	1 - 14
SiliaPrepX SCX	SPE-P0005-xxx	85 µm	60 Å	800 m ² /g	80 %	-	0.80 meq/g	1 - 14
SiliaPrepX SAX	SPE-P0010-xxx	85 µm	60 Å	900 m ² /g	85 %	-	0.20 meq/g	1 - 14
SiliaPrepX WCX	SPE-P0015-xxx	85 µm	60 Å	800 m ² /g	85 %	-	0.70 meq/g	1 - 14
SiliaPrepX WAX	SPE-P0020-xxx	85 µm	60 Å	800 m ² /g	86 %	-	0.50 meq/g	1 - 14

* Typical values

Typical Applications - Reversed and Normal Phases

The table below will help you select the right media to purify your compounds of interest, in either reversed-phase or normal phase.

SPE Cartridges & Well Plates Portfolio (<i>Reversed and Normal Phases</i>)		
Mode	SiliaPrep Phases	Applications
Reversed-Phases: non-polar sorbents	SiliaPrep C18 (<i>Plus, WPD Widepore, nec</i>)	For organic compounds from water, drugs & metabolites from fluids
	SiliaPrep C8 (<i>endcapped & nec</i>)	For extremely non-polar and large compounds (<i>vitamin D, oils</i>)
	SiliaPrep Phenyl (<i>PH</i>) & Pentafluorophenyl (<i>PPF</i>)	For aromatic compounds, complex natural products
Polymeric Reversed-Phases	SiliaPrepX HLB & DVB	For drugs & metabolites from biological fluids, API from tablets and cream
Normal Phases: polar sorbents	SiliaPrep Cyano (<i>CN</i>)	For acidic, basic and neutral compounds from aqueous solutions
	SiliaPrep Diol <i>nec</i>	For polar compounds from non-polar solvents, structural isomers
	SiliaPrep Florisil & Florisil PR (<i>Pesticide Residues</i>)	For chlorinated pesticides, PCB's and polysaccharides
	SiliaPrep Silica & Silica WPD (<i>Widepore</i>)	For various compounds from non-polar solvents, structural isomers
	SiliaPrep Alumina (<i>Acidic, Neutral & Basic</i>)	For aromatic compounds and aliphatic amines

Experimental Procedures - Reversed and Normal Phases

Generic protocols are presented below, for reversed-phase and normal phase SPE, to help you develop your method depending on the sorbent used, the sample matrix and the analyte properties.

These are only convenient starting points for method development. Further optimization may be required to tailor the method to the application needs.

Reversed-Phases

Extraction of neutral, acidic & basic organic compounds

Extraction of neutral, acidic & basic organic compounds	
CONDITIONING STEP	1 x CV ⁽¹⁾ of Methanol
EQUILIBRATION STEP	1 x CV of water
LOADING STEP	Aqueous sample, pH adjusted 2 units above pK_a (bases) or below pK_a (acids)
WASHING STEP	1 x CV of 5 % Methanol⁽²⁾ in water
ELUTION STEP	1 x CV of Methanol

Normal Phases

Extraction of compounds from non-polar solvents

Extraction of compounds from non-polar solvents	
CONDITIONING STEP	1 x CV of Isopropanol
EQUILIBRATION STEP	1 x CV of Hexane (or other low polar solvent)
LOADING STEP	Sample diluted in Hexane
WASHING STEP	1 x CV of 5 % Isopropanol in Hexane
ELUTION STEP	1 x CV of 50 - 95 % Isopropanol in Hexane

Notes:

⁽¹⁾ Abbreviation used: CV = Column Volume

⁽²⁾ For polymeric sorbents used in reversed-phase, you can add up to 60 % Methanol in water during the washing step, if your application requires it.

Typical Applications - Ion Exchange Phases

The table below will help you select the right media according to the pK_a of your analyte.

SPE Cartridges & Well Plates Portfolio (<i>Ion Exchange Phases</i>)		
Mode	SiliaPrep Phases	Applications
Ion Exchange Phases: Ionic sorbents	SiliaPrep SAX & SAX-2 (<i>TMA Chloride & Acetate</i>) nec	For weakly acidic molecules (pK_a 3 - 5)
	SiliaPrep Carbonate	For scavenging of TFA, extraction of acids (<i>boronic acids & acidic phenols</i>)
	SiliaPrep Amine (<i>WAX</i>)	For strongly acidic molecules ($pK_a < 3$), structural isomers, saccharides
	SiliaPrep SCX & SCX-2 (<i>Tosic & Propylsulfonic Acids</i>)	For weakly basic molecules (pK_a 7 - 9), catch & release of amines
	SiliaPrep WCX (<i>Carboxylic Acid</i>)	For strongly basic compounds ($pK_a > 9$)
Polymeric Ion Exchange Phases	SiliaPrepX SAX & WAX	For acidic compounds & metabolites, highly stable in organic solvents
	SiliaPrepX SCX & WCX	For basic compounds, highly stable in organic solvents

Experimental Procedures - Ion Exchange Phases

Strong Anion Exchangers (SAX)

Extraction of weak acids (pK_a 3 - 5 ⁽³⁾)	
CONDITIONNING STEP	1 x CV of Methanol
LOADING STEP	Aqueous sample, pH adjusted at 7 - 8
WASHING STEP	1 x CV of Methanol (<i>elution of basic & neutral compounds</i>)
ELUTION STEP	1 x CV of 2 - 5 % HCO₂H in Methanol (<i>elution of weak acidic compounds</i>)

Strong Cation Exchangers (SCX)

Extraction of weak bases (pK_a 7 - 9)	
CONDITIONNING STEP	1 x CV of Methanol
EQUILIBRATION STEP	1 x CV of water
LOADING STEP	Aqueous sample, pH adjusted at 3 - 4
WASHING STEP 1	1 x CV of water
WASHING STEP 2	1 x CV of Methanol (<i>elution of acidic & neutral compounds</i>)
ELUTION STEP	1 x CV of 2 - 5 % NH₄OH⁽⁴⁾ in Methanol (<i>elution of weak basic compounds</i>)

Weak Anion Exchangers (WAX)

Extraction of strong acids ($pK_a < 3$)	
CONDITIONNING STEP	1 x CV of Methanol
EQUILIBRATION STEP	1 x CV of water
LOADING STEP	Aqueous sample, pH adjusted at 4 - 5
WASHING STEP 1	1 x CV of water
WASHING STEP 2	1 x CV of Methanol (<i>elution of basic & neutral compounds</i>)
ELUTION STEP	1 x CV of 2 - 5 % NH₄OH⁽⁴⁾ in Methanol (<i>elution of strong acidic compounds</i>)

Weak Cation Exchangers (WCX)

Extraction of strong bases ($pK_a > 9$)	
CONDITIONNING STEP	1 x CV of Methanol
EQUILIBRATION STEP	1 x CV of water
LOADING STEP	Aqueous sample, pH adjusted at 8
WASHING STEP 1	1 x CV of water
WASHING STEP 2	1 x CV of Methanol (<i>elution of acidic & neutral compounds</i>)
ELUTION STEP	1 x CV of 2 - 5 % HCO₂H in Methanol (<i>elution of strong basic compounds</i>)

Notes:

⁽³⁾ For extraction of Phenol (pK_a 10), we recommend using a polymeric phase (SiliaPrepX SAX) and load the sample with a pH adjusted to 12.

⁽⁴⁾ For silica-based sorbents, NH₄OH can be too aggressive. You can use NH₃ (7M) in Methanol to avoid degrading the phase.

Typical Applications - Specialty Phases & Metal Scavengers

The table below presents our specialty phases, to remove specific compounds from your samples.

SPE Cartridges & Well Plates Portfolio (<i>Specialty Phases & Metal Scavengers</i>)		
Mode	SiliaPrep Phases	Applications
Specialty Phases	SiliaPrep PCB	For extraction of PCB's from waste oil (<i>hexane extract</i>)
	SiliaPrep CleanDRUG	For drugs of abuse applications
	SiliaPrep CleanENVI	For PAH's, PCB's, herbicides and pesticides from waste waters
	SiliaPrep PAH	For PAH's from waste waters
Metal Scavengers	SiliaPrep Cysteine, Diamine, DMT, DOTA, TAAcOH, TAAcONa, Thiol, Thiourea, Imidazole, Triamine	For lowering the residual metal concentration of various metal complexes (<i>Pd, Pt, Rh, Ru, Ni, Sn, etc</i>) to single digit ppm

Experimental Procedures - Specialty Phases & Metal Scavengers

The procedures below are only convenient starting points for method development. Further optimization may be required to tailor the method to your application needs.

Specialty Phases

Extraction of PCBs, drugs and PAHs

PCBs from waste oil with SiliaPrep PCB	
CONDITIONING STEP	1 x CV ⁽¹⁾ of Hexane
LOADING STEP	Diluted sample (<i>with Hexane</i>)
ELUTION STEP	1 x CV of Hexane

Drugs of abuse with SiliaPrep CleanDRUG	
CONDITIONING STEP	1 x CV of Methanol
EQUILIBRATION STEP	1 x CV of water (<i>buffered at pH 6</i>)
LOADING STEP	Aqueous sample (<i>buffered at pH 6</i>)
WASHING STEP	1 x CV of water then 1 x CV of Methanol
ELUTION STEP	1 x CV of Isopropanol:NH₄OH (90:10)

Environmental samples with SiliaPrep CleanENVI & PAH	
CONDITIONING STEP	1 x CV of Dichloromethane then 1 x CV of Methanol
EQUILIBRATION STEP	1 x CV of water
LOADING STEP	Aqueous sample
WASHING STEP	1 x CV of 5 - 50 % Methanol in water
ELUTION STEP	1 x CV of Dichloromethane

Notes:

⁽¹⁾ Abbreviation used: CV = Column Volume

⁽²⁾ Non retentive SPE (*Catch & Release*): analyte won't retain on the sorbent and will elute directly during loading and rinsing steps. Scavenged compounds will remain in the SPE cartridge.

Metal Scavengers

Catch of the metal & release of your analyte

Catch and release of the analyte	
EQUILIBRATION STEP	1 x CV of sample solvent
LOADING STEP ⁽²⁾	Diluted sample
RINSING STEP	1 x CV of sample solvent

Product Selection Guide by Manufacturer

The table below will help you find equivalences to products of well-known SPE manufacturers.

Product Selection Guide by Manufacturer					
SiliaCycle	Waters	Phenomenex	Agilent	Biotage	Macherey-Nagel
SiliaPrep C18 Plus	Sep-Pak® tC18	Strata® C18-E	Bond Elut® C18	Isolute® C18 (EC)	Chromabond® C18 ec
SiliaPrep C18 nec		Strata® C18-U	Bond Elut® C18 OH	Isolute® C18	Chromabond® C18
SiliaPrep C18 WPD	Sep-Pak® C18	Strata® C18-T	Bond Elut® C18 EWP	Isolute® MFC18	Chromabond® C18 ec f
SiliaPrep C8	Sep-Pak® C8	Strata® C8	Bond Elut® C8	Isolute® C8 (EC)	
SiliaPrep C8 nec				Isolute® C8	Chromabond® C8
SiliaPrep Phenyl (PH)		Strata® Phenyl	Bond Elut® PH	Isolute® PH	Chromabond® C ₆ H ₅
SiliaPrep PFP					
SiliaPrep Cyano (CN)	Sep-Pak® Cyanopropyl	Strata® CN	Bond Elut® Cyano (CN)	Isolute® CN	Chromabond® CN
SiliaPrep Diol nec	Sep-Pak® Diol		Bond Elut® Diol (2OH)	Isolute® DIOL	Chromabond® OH (Diol)
SiliaPrep Silica		Strata® Silica (Si-1)	Bond Elut® SI	Isolute® SI	Chromabond® SiOH
SiliaPrep Silica WPD	Sep-Pak® Silica				
SiliaPrep Florisil LP & Florisil PR	Sep-Pak® Florisil®	Strata® FL-PR (Florisil®)	Bond Elut® Florisil	Isolute® FL	Chromabond® Florisil®
SiliaPrep Alumina (Acidic, Neutral, Basic)	Sep-Pak® Alumina (A, N, B)	Strata® Alumina-N (AL-N)	Bond Elut® Alumina (-A, -N, -B)	Isolute® ALUMINA (AL-A, AL-N & AL-B)	Chromabond® AloX (A, N, B)
SiliaPrep SAX nec (TMA Chloride)	Sep-Pak® Accell™ Plus QMA	Strata® SAX	Bond Elut® SAX	Isolute® SAX	Chromabond® SB
SiliaPrep SAX-2 nec (TMA Acetate)				Isolute® PE-AX	
SiliaPrep Carbonate	Accell Plus QMA Carbonate Plus Light			Isolute® Si-Carbonate (Si-TMA-CO ₃)	
SiliaPrep Amine (WAX)	Sep-Pak® Amino	Strata® NH ₂	Bond Elut® NH ₂	Isolute® NH ₂	Chromabond® NH ₂
SiliaPrep Tosic Acid (SCX)		Strata® SCX	Bond Elut® SCX	Isolute® SCX-3	Chromabond® SA
SiliaPrep SCX-2 (Propylsulfonic Acid)			Bond Elut® PRS	Isolute® SCX-2	Chromabond® PSA
SiliaPrep WCX (Carboxylic Acid)	Sep-Pak® Accell™ Plus CM	Strata® WCX	Bond Elut® CBA	Isolute® CBA	Chromabond® PCA
SiliaPrep PCB			Bond Elut® PCB		Chromabond® SA/SiOH
SiliaPrep CleanDRUG		Strata® Screen-C	Bond Elut® Certify	Isolute® HCX	Chromabond® Drug
SiliaPrep CleanENVI & SiliaPrep PAH		Strata® PAH	EnvirElut®	Isolute® PAH	Chromabond® C18 PAH
SiliaPrepX HLB	Oasis® HLB	Strata®-X	Bond Elut® NEXUS		Chromabond® HLB
SiliaPrepX DVB			Bond Elut® ENV	Isolute® 101	Chromabond® HR-X
SiliaPrepX SAX	Oasis® MAX	Strata®-X-A	Bond Elut® Plexa PAX	Evolute® Express AX	Chromabond® HR-XA
SiliaPrepX WAX	Oasis® WAX	Strata®-X-AW		Evolute® Express WAX	Chromabond® HR-XAW
SiliaPrepX SCX	Oasis® MCX	Strata®-X-C	Bond Elut® Plexa PCX	Evolute® Express CX	Chromabond® HR-XC
SiliaPrepX WCX	Oasis® WCX	Strata®-X-CW	Bond Elut® NEXUS WCX	Evolute® Express WCX	Chromabond® HR-XCW

All trademarks and registered trademarks are the property of their respective owners. SiliaCycle takes no responsibility for any error or omission relating to this information.

Ordering Information - Reversed-Phases

We can provide a complete range of SPE cartridge volumes and bed weights, and 96-Well plates.

SiliaPrep Reversed-Phases SPE & Well Plates Formats					
Formats	Qty/Box	SiliaPrep C18 Plus	SiliaPrep C18 WPD	SiliaPrep C18 nec	SiliaPrep C8
SiliaPrep SPE Cartridges					
1 mL / 50 mg	100	SPE-R00830B-01B	SPE-R33229G-01B	SPE-R35530B-01B	SPE-R31030B-01B
1 mL / 100 mg	100	SPE-R00830B-01C	SPE-R33229G-01C	SPE-R35530B-01C	SPE-R31030B-01C
3 mL / 200 mg	50	SPE-R00830B-03G	SPE-R33229G-03G	SPE-R35530B-03G	SPE-R31030B-03G
3 mL / 500 mg	50	SPE-R00830B-03P	SPE-R33229G-03P	SPE-R35530B-03P	SPE-R31030B-03P
6 mL / 500 mg	50	SPE-R00830B-06P	SPE-R33229G-06P	SPE-R35530B-06P	SPE-R31030B-06P
6 mL / 1 g	50	SPE-R00830B-06S	SPE-R33229G-06S	SPE-R35530B-06S	SPE-R31030B-06S
6 mL / 2 g	50	SPE-R00830B-06U	SPE-R33229G-06U	SPE-R35530B-06U	SPE-R31030B-06U
12 mL / 2 g	20	SPE-R00830B-12U	SPE-R33229G-12U	SPE-R35530B-12U	SPE-R31030B-12U
25 mL / 5 g*	20	SPE-R00830B-20X	SPE-R33229G-20X	SPE-R35530B-20X	SPE-R31030B-20X
70 mL / 10 g*	16	FLH-R00830B-70Y	FLH-R33229G-70Y	FLH-R35530B-70Y	FLH-R31030B-70Y
70 mL / 15 g*	16	FLH-R00830B-70i	FLH-R33229G-70i	FLH-R35530B-70i	FLH-R31030B-70i
70 mL / 20 g*	16	FLH-R00830B-70Z	FLH-R33229G-70Z	FLH-R35530B-70Z	FLH-R31030B-70Z
150 mL / 25 g*	10	FLH-R00830B-95K	FLH-R33229G-95K	FLH-R35530B-95K	FLH-R31030B-95K
150 mL / 50 g*	10	FLH-R00830B-95M	FLH-R33229G-95M	FLH-R35530B-95M	FLH-R31030B-95M
150 mL / 70 g*	10	FLH-R00830B-95N	FLH-R33229G-95N	FLH-R35530B-95N	FLH-R31030B-95N
276 mL / 100 g*	12	FLH-R00830B-276F	FLH-R33229G-276F	FLH-R35530B-276F	FLH-R31030B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL / 200 mg	50	SPC-R00830B-10G	SPC-R33229G-10G	SPC-R35530B-10G	SPC-R31030B-10G
10 mL / 500 mg	50	SPC-R00830B-10P	SPC-R33229G-10P	SPC-R35530B-10P	SPC-R31030B-10P
Mini-SiliaPrep SPE Cartridges					
500 mg	50	SPS-R00830B-P	SPS-R33229G-P	SPS-R35530B-P	SPS-R31030B-P
1 g	50	SPS-R00830B-S	SPS-R33229G-S	SPS-R35530B-S	SPS-R31030B-S
SiliaPrep 96-Well Plates					
2 mL / 50 mg	1	96W-R00830B-B	96W-R33229G-B	96W-R35530B-B	96W-R31030B-B
2 mL / 100 mg	1	96W-R00830B-C	96W-R33229G-C	96W-R35530B-C	96W-R31030B-C

* Commercialized under SiliaSep OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

SiliaPrep Reversed-Phases SPE & Well Plates Formats				
Formats	Qty/Box	SiliaPrep C8 nec	SiliaPrep Phenyl (PH)	SiliaPrep PFP
SiliaPrep SPE Cartridges				
1 mL / 50 mg	100	SPE-R31130B-01B	SPE-R34030B-01B	SPE-R67530B-01B
1 mL / 100 mg	100	SPE-R31130B-01C	SPE-R34030B-01C	SPE-R67530B-01C
3 mL / 200 mg	50	SPE-R31130B-03G	SPE-R34030B-03G	SPE-R67530B-03G
3 mL / 500 mg	50	SPE-R31130B-03P	SPE-R34030B-03P	SPE-R67530B-03P
6 mL / 500 mg	50	SPE-R31130B-06P	SPE-R34030B-06P	SPE-R67530B-06P
6 mL / 1 g	50	SPE-R31130B-06S	SPE-R34030B-06S	SPE-R67530B-06S
6 mL / 2 g	50	SPE-R31130B-06U	SPE-R34030B-06U	SPE-R67530B-06U
12 mL / 2 g	20	SPE-R31130B-12U	SPE-R34030B-12U	SPE-R67530B-12U
25 mL / 5 g*	20	SPE-R31130B-20X	SPE-R34030B-20X	SPE-R67530B-20X
70 mL / 10 g*	16	FLH-R31130B-70Y	FLH-R34030B-70Y	FLH-R67530B-70Y
70 mL / 15 g*	16	FLH-R31130B-70i	FLH-R34030B-70i	FLH-R67530B-70i
70 mL / 20 g*	16	FLH-R31130B-70Z	FLH-R34030B-70Z	FLH-R67530B-70Z
150 mL / 25 g*	10	FLH-R31130B-95K	FLH-R34030B-95K	FLH-R67530B-95K
150 mL / 50 g*	10	FLH-R31130B-95M	FLH-R34030B-95M	FLH-R67530B-95M
150 mL / 70 g*	10	FLH-R31130B-95N	FLH-R34030B-95N	FLH-R67530B-95N
276 mL / 100 g*	12	FLH-R31130B-276F	FLH-R34030B-276F	FLH-R67530B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL / 200 mg	50	SPC-R31130B-10G	SPC-R34030B-10G	SPC-R67530B-10G
10 mL / 500 mg	50	SPC-R31130B-10P	SPC-R34030B-10P	SPC-R67530B-10P
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-R31130B-P	SPS-R34030B-P	SPS-R67530B-P
1 g	50	SPS-R31130B-S	SPS-R34030B-S	SPS-R67530B-S
SiliaPrep 96-Well Plates				
2 mL / 50 mg	1	96W-R31130B-B	96W-R34030B-B	96W-R67530B-B
2 mL / 100 mg	1	96W-R31130B-C	96W-R34030B-C	96W-R67530B-C

* Commercialized under SiliaSep OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

Ordering Information - Normal Phases

We can provide a complete range of SPE cartridge volumes and bed weights, and 96-Well plates.

SiliaPrep Normal Phases SPE & Well Plates Formats						
Formats	Qty/Box	SiliaPrep Cyano (CM)	SiliaPrep Diol nec	SiliaPrep Florisil	SiliaPrep Florisil LP	SiliaPrep Florisil PR
SiliaPrep SPE Cartridges						
1 mL / 50 mg	100	SPE-R38030B-01B	SPE-R35030B-01B	SPE-AUT-0014-01B	SPE-AUT-0014LP-01B	SPE-AUT-0015-01B
1 mL / 100 mg	100	SPE-R38030B-01C	SPE-R35030B-01C	SPE-AUT-0014-01C	SPE-AUT-0014LP-01C	SPE-AUT-0015-01C
3 mL / 200 mg	50	SPE-R38030B-03G	SPE-R35030B-03G	SPE-AUT-0014-03G	SPE-AUT-0014LP-03G	SPE-AUT-0015-03G
3 mL / 500 mg	50	SPE-R38030B-03P	SPE-R35030B-03P	SPE-AUT-0014-03P	SPE-AUT-0014LP-03P	SPE-AUT-0015-03P
6 mL / 500 mg	50	SPE-R38030B-06P	SPE-R35030B-06P	SPE-AUT-0014-06P	SPE-AUT-0014LP-06P	SPE-AUT-0015-06P
6 mL / 1 g	50	SPE-R38030B-06S	SPE-R35030B-06S	SPE-AUT-0014-06S	SPE-AUT-0014LP-06S	SPE-AUT-0015-06S
6 mL / 2 g	50	SPE-R38030B-06U	SPE-R35030B-06U	SPE-AUT-0014-06U	SPE-AUT-0014LP-06U	SPE-AUT-0015-06U
12 mL / 2 g	20	SPE-R38030B-12U	SPE-R35030B-12U	SPE-AUT-0014-12U	SPE-AUT-0014LP-12U	SPE-AUT-0015-12U
25 mL / 5 g*	20	SPE-R38030B-20X	SPE-R35030B-20X	SPE-AUT-0014-20X	SPE-AUT-0014LP-20X	SPE-AUT-0015-20X
70 mL / 10 g*	16	FLH-R38030B-70Y	FLH-R35030B-70Y	FLH-AUT-0014-70Y	FLH-AUT-0014LP-70Y	FLH-AUT-0015-70Y
70 mL / 15 g*	16	FLH-R38030B-70i	FLH-R35030B-70i	FLH-AUT-0014-70i	FLH-AUT-0014LP-70i	FLH-AUT-0015-70i
70 mL / 20 g*	16	FLH-R38030B-70Z	FLH-R35030B-70Z	FLH-AUT-0014-70Z	FLH-AUT-0014LP-70Z	FLH-AUT-0015-70Z
150 mL / 25 g*	10	FLH-R38030B-95K	FLH-R35030B-95K	FLH-AUT-0014-95K	FLH-AUT-0014LP-95K	FLH-AUT-0015-95K
150 mL / 50 g*	10	FLH-R38030B-95M	FLH-R35030B-95M	FLH-AUT-0014-95M	FLH-AUT-0014LP-95M	FLH-AUT-0015-95M
150 mL / 70 g*	10	FLH-R38030B-95N	FLH-R35030B-95N	FLH-AUT-0014-95N	FLH-AUT-0014LP-95N	FLH-AUT-0015-95N
276 mL / 100 g*	12	FLH-R38030B-276F	FLH-R35030B-276F	FLH-AUT-0014-276F	FLH-AUT-0014LP-276F	FLH-AUT-0015-276F
SiliaPrep Large Reservoir Volume SPE Cartridges						
10 mL / 200 mg	50	SPC-R38030B-10G	SPC-R35030B-10G	SPC-AUT-0014-10G	SPC-AUT-0014LP-10G	SPC-AUT-0015-10G
10 mL / 500 mg	50	SPC-R38030B-10P	SPC-R35030B-10P	SPC-AUT-0014-10P	SPC-AUT-0014LP-10P	SPC-AUT-0015-10P
Mini-SiliaPrep SPE Cartridges						
500 mg	50	SPS-R38030B-P	SPS-R35030B-P	SPS-AUT-0014-P	SPS-AUT-0014LP-P	SPS-AUT-0015-P
1 g	50	SPS-R38030B-S	SPS-R35030B-S	SPS-AUT-0014-S	SPS-AUT-0014LP-S	SPS-AUT-0015-S
SiliaPrep 96-Well Plates						
2 mL / 50 mg	1	96W-R38030B-B	96W-R35030B-B	96W-AUT-0014-B	96W-AUT-0014LP-B	96W-AUT-0015-B
2 mL / 100 mg	1	96W-R38030B-C	96W-R35030B-C	96W-AUT-0014-C	96W-AUT-0014LP-C	96W-AUT-0015-C

* Commercialized under SiliaSep OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

SiliaPrep Normal Phases SPE & Well Plates Formats						
Formats	Qty/Box	SiliaPrep Silica	SiliaPrep Silica WPD	SiliaPrep Acidic Alumina	SiliaPrep Neutral Alumina	SiliaPrep Basic Alumina
SiliaPrep SPE Cartridges						
1 mL / 50 mg	100	SPE-R10030B-01B	SPE-R10029G-01B	SPE-AUT-0053-01B	SPE-AUT-0054-01B	SPE-AUT-0055-01B
1 mL / 100 mg	100	SPE-R10030B-01C	SPE-R10029G-01C	SPE-AUT-0053-01C	SPE-AUT-0054-01C	SPE-AUT-0055-01C
3 mL / 200 mg	50	SPE-R10030B-03G	SPE-R10029G-03G	SPE-AUT-0053-03G	SPE-AUT-0054-03G	SPE-AUT-0055-03G
3 mL / 500 mg	50	SPE-R10030B-03P	SPE-R10029G-03P	SPE-AUT-0053-03P	SPE-AUT-0054-03P	SPE-AUT-0055-03P
6 mL / 500 mg	50	SPE-R10030B-06P	SPE-R10029G-06P	SPE-AUT-0053-06P	SPE-AUT-0054-06P	SPE-AUT-0055-06P
6 mL / 1 g	50	SPE-R10030B-06S	SPE-R10029G-06S	SPE-AUT-0053-06S	SPE-AUT-0054-06S	SPE-AUT-0055-06S
6 mL / 2 g	50	SPE-R10030B-06U	SPE-R10029G-06U	SPE-AUT-0053-06U	SPE-AUT-0054-06U	SPE-AUT-0055-06U
12 mL / 2 g	20	FLH-R10030B-15U*	SPE-R10029G-12U	SPE-AUT-0053-12U	SPE-AUT-0054-12U	SPE-AUT-0055-12U
25 mL / 5 g*	20	FLH-R10030B-25X	SPE-R10029G-20X	SPE-AUT-0053-20X	SPE-AUT-0054-20X	SPE-AUT-0055-20X
70 mL / 10 g*	16	FLH-R10030B-70Y	FLH-R10029G-70Y	FLH-AUT-0053-70Y	FLH-AUT-0054-70Y	FLH-AUT-0055-70Y
70 mL / 15 g*	16	FLH-R10030B-70i	FLH-R10029G-70i	FLH-AUT-0053-70i	FLH-AUT-0054-70i	FLH-AUT-0055-70i
70 mL / 20 g*	16	FLH-R10030B-70Z	FLH-R10029G-70Z	FLH-AUT-0053-70Z	FLH-AUT-0054-70Z	FLH-AUT-0055-70Z
150 mL / 25 g*	10	FLH-R10030B-95K	FLH-R10029G-95K	FLH-AUT-0053-95K	FLH-AUT-0054-95K	FLH-AUT-0055-95K
150 mL / 50 g*	10	FLH-R10030B-95M	FLH-R10029G-95M	FLH-AUT-0053-95M	FLH-AUT-0054-95M	FLH-AUT-0055-95M
150 mL / 70 g*	10	FLH-R10030B-95N	FLH-R10029G-95N	FLH-AUT-0053-95N	FLH-AUT-0054-95N	FLH-AUT-0055-95N
276 mL / 100 g*	12	FLH-R10030B-276F	FLH-R10029G-276F	FLH-AUT-0053-276F	FLH-AUT-0054-276F	FLH-AUT-0055-276F
SiliaPrep Large Reservoir Volume SPE Cartridges						
10 mL / 200 mg	50	SPC-R10030B-10G	SPC-R10029G-10G	SPC-AUT-0053-10G	SPC-AUT-0054-10G	SPC-AUT-0055-10G
10 mL / 500 mg	50	SPC-R10030B-10P	SPC-R10029G-10P	SPC-AUT-0053-10P	SPC-AUT-0054-10P	SPC-AUT-0055-10P
Mini-SiliaPrep SPE Cartridges						
500 mg	50	SPS-R10030B-P	SPS-R10029G-P	SPS-AUT-0053-P	SPS-AUT-0054-P	SPS-AUT-0055-P
1 g	50	SPS-R10030B-S	SPS-R10029G-S	SPS-AUT-0053-S	SPS-AUT-0054-S	SPS-AUT-0055-S
SiliaPrep 96-Well Plates						
2 mL / 50 mg	1	96W-R10030B-B	96W-R10029G-B	-	-	-
2 mL / 100 mg	1	96W-R10030B-C	96W-R10029G-C	-	-	-

* Commercialized under SiliaSep OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

Ordering Information - Ion Exchange Phases

We can provide a complete range of SPE cartridge volumes and bed weights, and 96-Well plates.

SiliaPrep Ion Exchange Phases SPE & Well Plates Formats					
Formats	Qty/Box	SiliaPrep SAX <i>nec</i>	SiliaPrep SAX-2 <i>nec</i>	SiliaPrep Carbonate	SiliaPrep Amine (WAX)
SiliaPrep SPE Cartridges					
1 mL / 50 mg	100	SPE-R66530B-01B	SPE-R66430B-01B	SPE-R66030B-01B	SPE-R52030B-01B
1 mL / 100 mg	100	SPE-R66530B-01C	SPE-R66430B-01C	SPE-R66030B-01C	SPE-R52030B-01C
3 mL / 200 mg	50	SPE-R66530B-03G	SPE-R66430B-03G	SPE-R66030B-03G	SPE-R52030B-03G
3 mL / 500 mg	50	SPE-R66530B-03P	SPE-R66430B-03P	SPE-R66030B-03P	SPE-R52030B-03P
6 mL / 500 mg	50	SPE-R66530B-06P	SPE-R66430B-06P	SPE-R66030B-06P	SPE-R52030B-06P
6 mL / 1 g	50	SPE-R66530B-06S	SPE-R66430B-06S	SPE-R66030B-06S	SPE-R52030B-06S
6 mL / 2 g	50	SPE-R66530B-06U	SPE-R66430B-06U	SPE-R66030B-06U	SPE-R52030B-06U
12 mL / 2 g	20	SPE-R66530B-12U	SPE-R66430B-12U	SPE-R66030B-12U	SPE-R52030B-12U
25 mL / 5 g*	20	SPE-R66530B-20X	SPE-R66430B-20X	SPE-R66030B-20X	SPE-R52030B-20X
70 mL / 10 g*	16	FLH-R66530B-70Y	FLH-R66430B-70Y	FLH-R66030B-70Y	FLH-R52030B-70Y
70 mL / 15 g*	16	FLH-R66530B-70i	FLH-R66430B-70i	FLH-R66030B-70i	FLH-R52030B-70i
70 mL / 20 g*	16	FLH-R66530B-70Z	FLH-R66430B-70Z	FLH-R66030B-70Z	FLH-R52030B-70Z
150 mL / 25 g*	10	FLH-R66530B-95K	FLH-R66430B-95K	FLH-R66030B-95K	FLH-R52030B-95K
150 mL / 50 g*	10	FLH-R66530B-95M	FLH-R66430B-95M	FLH-R66030B-95M	FLH-R52030B-95M
150 mL / 70 g*	10	FLH-R66530B-95N	FLH-R66430B-95N	FLH-R66030B-95N	FLH-R52030B-95N
276 mL / 100 g*	12	FLH-R66530B-276F	FLH-R66430B-276F	FLH-R66030B-276F	FLH-R52030B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL / 200 mg	50	SPC-R66530B-10G	SPC-R66430B-10G	SPC-R66030B-10G	SPC-R52030B-10G
10 mL / 500 mg	50	SPC-R66530B-10P	SPC-R66430B-10P	SPC-R66030B-10P	SPC-R52030B-10P
Mini-SiliaPrep SPE Cartridges					
500 mg	50	SPS-R66530B-P	SPS-R66430B-P	SPS-R66030B-P	SPS-R52030B-P
1 g	50	SPS-R66530B-S	SPS-R66430B-S	SPS-R66030B-S	SPS-R52030B-S
SiliaPrep 96-Well Plates					
2 mL / 50 mg	1	96W-R66530B-B	96W-R66430B-B	96W-R66030B-B	96W-R52030B-B
2 mL / 100 mg	1	96W-R66530B-C	96W-R66430B-C	96W-R66030B-C	96W-R52030B-C

* Commercialized under SiliaSep OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

SiliaPrep Ion Exchange Phases SPE & Well Plates Formats				
Formats	Qty/Box	SiliaPrep Tonic Acid (SCX)	SiliaPrep SCX-2	SiliaPrep WCX
SiliaPrep SPE Cartridges				
1 mL / 50 mg	100	SPE-R60530B-01B	SPE-R51230B-01B	SPE-R70030B-01B
1 mL / 100 mg	100	SPE-R60530B-01C	SPE-R51230B-01C	SPE-R70030B-01C
3 mL / 200 mg	50	SPE-R60530B-03G	SPE-R51230B-03G	SPE-R70030B-03G
3 mL / 500 mg	50	SPE-R60530B-03P	SPE-R51230B-03P	SPE-R70030B-03P
6 mL / 500 mg	50	SPE-R60530B-06P	SPE-R51230B-06P	SPE-R70030B-06P
6 mL / 1 g	50	SPE-R60530B-06S	SPE-R51230B-06S	SPE-R70030B-06S
6 mL / 2 g	50	SPE-R60530B-06U	SPE-R51230B-06U	SPE-R70030B-06U
12 mL / 2 g	20	SPE-R60530B-12U	SPE-R51230B-12U	SPE-R70030B-12U
25 mL / 5 g*	20	SPE-R60530B-20X	SPE-R51230B-20X	SPE-R70030B-20X
70 mL / 10 g*	16	FLH-R60530B-70Y	FLH-R51230B-70Y	FLH-R70030B-70Y
70 mL / 15 g*	16	FLH-R60530B-70i	FLH-R51230B-70i	FLH-R70030B-70i
70 mL / 20 g*	16	FLH-R60530B-70Z	FLH-R51230B-70Z	FLH-R70030B-70Z
150 mL / 25 g*	10	FLH-R60530B-95K	FLH-R51230B-95K	FLH-R70030B-95K
150 mL / 50 g*	10	FLH-R60530B-95M	FLH-R51230B-95M	FLH-R70030B-95M
150 mL / 70 g*	10	FLH-R60530B-95N	FLH-R51230B-95N	FLH-R70030B-95N
276 mL / 100 g*	12	FLH-R60530B-276F	FLH-R51230B-276F	FLH-R70030B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges				
10 mL / 200 mg	50	SPC-R60530B-10G	SPC-R51230B-10G	SPC-R70030B-10G
10 mL / 500 mg	50	SPC-R60530B-10P	SPC-R51230B-10P	SPC-R70030B-10P
Mini-SiliaPrep SPE Cartridges				
500 mg	50	SPS-R60530B-P	SPS-R51230B-P	SPS-R70030B-P
1 g	50	SPS-R60530B-S	SPS-R51230B-S	SPS-R70030B-S
SiliaPrep 96-Well Plates				
2 mL / 50 mg	1	96W-R60530B-B	96W-R51230B-B	96W-R70030B-B
2 mL / 100 mg	1	96W-R60530B-C	96W-R51230B-C	96W-R70030B-C

* Commercialized under SiliaSep OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

Ordering Information - Polymeric Phases

We can provide a complete range of SPE cartridge volumes and bed weights, and 96-Well plates.

SiliaPrepX Polymeric Phases SPE & Well Plates Formats				
Formats	Qty/Box	SiliaPrepX DVB	SiliaPrepX HLB	SiliaPrepX SCX
SiliaPrepX SPE Cartridges				
1 mL / 30 mg	100	SPE-P0001-01AA	SPE-P0002-01AA	SPE-P0005-01AA
3 mL / 30 mg	50	SPE-P0001-03AA	SPE-P0002-03AA	SPE-P0005-03AA
3 mL / 60 mg	50	SPE-P0001-03BB	SPE-P0002-03BB	SPE-P0005-03BB
6 mL / 100 mg	30	SPE-P0001-06C	SPE-P0002-06C	SPE-P0005-06C
6 mL / 200 mg	30	SPE-P0001-06G	SPE-P0002-06G	SPE-P0005-06G
6 mL / 500 mg	30	SPE-P0001-06P	SPE-P0002-06P	SPE-P0005-06P
12 mL / 500 mg	20	SPE-P0001-12P	SPE-P0002-12P	SPE-P0005-12P
12 mL / 1 g	20	SPE-P0001-12S	SPE-P0002-12S	SPE-P0005-12S
25 mL / 1 g*	20	SPE-P0001-20S	SPE-P0002-20S	SPE-P0005-20S
25 mL / 2 g*	20	SPE-P0001-20U	SPE-P0002-20U	SPE-P0005-20U
70 mL / 5 g*	16	FLH-P0001-70X	FLH-P0002-70X	FLH-P0005-70X
SiliaPrepX 96-Well Plates				
2 mL / 10 mg	1	96W-P0001-1A	96W-P0002-1A	96W-P0005-1A
2 mL / 30 mg	1	96W-P0001-AA	96W-P0002-AA	96W-P0005-AA
2 mL / 60 mg	1	96W-P0001-BB	96W-P0002-BB	96W-P0005-BB

SiliaPrepX Polymeric Phases SPE & Well Plates Formats				
Formats	Qty/Box	SiliaPrepX SAX	SiliaPrepX WCX	SiliaPrepX WAX
SiliaPrepX SPE Cartridges				
1 mL / 30 mg	100	SPE-P0010-01AA	SPE-P0015-01AA	SPE-P0020-01AA
3 mL / 30 mg	50	SPE-P0010-03AA	SPE-P0015-03AA	SPE-P0020-03AA
3 mL / 60 mg	50	SPE-P0010-03BB	SPE-P0015-03BB	SPE-P0020-03BB
6 mL / 100 mg	30	SPE-P0010-06C	SPE-P0015-06C	SPE-P0020-06C
6 mL / 200 mg	30	SPE-P0010-06G	SPE-P0015-06G	SPE-P0020-06G
6 mL / 500 mg	30	SPE-P0010-06P	SPE-P0015-06P	SPE-P0020-06P
12 mL / 500 mg	20	SPE-P0010-12P	SPE-P0015-12P	SPE-P0020-12P
12 mL / 1 g	20	SPE-P0010-12S	SPE-P0015-12S	SPE-P0020-12S
25 mL / 1 g*	20	SPE-P0010-20S	SPE-P0015-20S	SPE-P0020-20S
25 mL / 2 g*	20	SPE-P0010-20U	SPE-P0015-20U	SPE-P0020-20U
70 mL / 5 g*	16	FLH-P0010-70X	FLH-P0015-70X	FLH-P0020-70X
SiliaPrepX 96-Well Plates				
2 mL / 10 mg	1	96W-P0010-1A	96W-P0015-1A	96W-P0020-1A
2 mL / 30 mg	1	96W-P0010-AA	96W-P0015-AA	96W-P0020-AA
2 mL / 60 mg	1	96W-P0010-BB	96W-P0015-BB	96W-P0020-BB

* Commercialized under SiliaSepX OT branding.

Note:

- Custom formats available on request.
- Add "-J" at the end of the part number to order a box of 200 cartridges.

www.greyhoundchrom.com Email: sales@greyhoundchrom.com Tel: +44 (0) 151 649 4000

Ordering Information - Specialty Phases

We can provide a complete range of SPE cartridge volumes and bed weights.

SiliaPrep Specialty Phases SPE Formats					
Formats	Qty/Box	SiliaPrep PCB	SiliaPrep CleanDRUG	SiliaPrep CleanENVI	SiliaPrep PAH
SiliaPrep SPE Cartridges					
1 mL / 50 mg	100	SP2-R00650030B-01B	SPEC-R651230B-01B	SPEC-R31930B-01B	SP2-R0610030B-01B
1 mL / 100 mg	100	SP2-R00650030B-01C	SPEC-R651230B-01C	SPEC-R31930B-01C	SP2-R0610030B-01C
3 mL / 200 mg	50	SP2-R00650030B-03G	SPEC-R651230B-03G	SPEC-R31930B-03G	SP2-R0610030B-03G
3 mL / 500 mg	50	SP2-R00650030B-03P	SPEC-R651230B-03P	SPEC-R31930B-03P	SP2-R0610030B-03P
6 mL / 500 mg	50	SP2-R00650030B-06P	SPEC-R651230B-06P	SPEC-R31930B-06P	SP2-R0610030B-06P
6 mL / 1 g	50	SP2-R00650030B-06S	SPEC-R651230B-06S	SPEC-R31930B-06S	SP2-R0610030B-06S
6 mL / 2 g	50	SP2-R00650030B-06U	SPEC-R651230B-06U	SPEC-R31930B-06U	SP2-R0610030B-06U
12 mL / 2 g	20	SP2-R00650030B-12U	SPEC-R651230B-12U	SPEC-R31930B-12U	SP2-R0610030B-12U
25 mL / 5 g*	20	SP2-R00650030B-20X	SPEC-R651230B-20X	SPEC-R31930B-20X	SP2-R0610030B-20X
70 mL / 10 g*	16	FLH-R00650030B-70Y	FLH-R651230B-70Y	FLH-R31930B-70Y	FLH-R0610030B-70Y
70 mL / 15 g*	16	FLH-R00650030B-70i	FLH-R651230B-70i	FLH-R31930B-70i	FLH-R0610030B-70i
70 mL / 20 g*	16	FLH-R00650030B-70Z	FLH-R651230B-70Z	FLH-R31930B-70Z	FLH-R0610030B-70Z
150 mL / 25 g*	10	FLH-R00650030B-95K	FLH-R651230B-95K	FLH-R31930B-95K	FLH-R0610030B-95K
150 mL / 50 g*	10	FLH-R00650030B-95M	FLH-R651230B-95M	FLH-R31930B-95M	FLH-R0610030B-95M
150 mL / 70 g*	10	FLH-R00650030B-95N	FLH-R651230B-95N	FLH-R31930B-95N	FLH-R0610030B-95N
276 mL / 100 g*	12	FLH-R00650030B-276F	FLH-R651230B-276F	FLH-R31930B-276F	FLH-R0610030B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL / 200 mg	50	SPC-R00650030B-10G	SPC-R651230B-10G	SPC-R31930B-10G	SPC-R0610030B-10G
10 mL / 500 mg	50	SPC-R00650030B-10P	SPC-R651230B-10P	SPC-R31930B-10P	SPC-R0610030B-10P
Mini-SiliaPrep SPE Cartridges					
500 mg	50	SPS-R00650030B-P	SPS-R651230B-P	SPS-R31930B-P	SPS-R0610030B-P
1 g	50	SPS-R00650030B-S	SPS-R651230B-S	SPS-R31930B-S	SPS-R0610030B-S

* Commercialized under SiliaSep OT branding.

Note: Custom formats available on request.

Ordering Information - Metal Scavenger Phases

To find out which SiliaPrep Metal Scavenger will better suit your need, we recommend performing some screening using the SiliaPrep Metal Scavenger Kit (PN: SPE-K30730B-03P) to quickly determine which scavenger presents the highest efficiency, and optimize the results.

We can provide a complete range of SPE cartridge volumes and bed weights.

SiliaPrep Metal Scavenger Phases SPE Formats					
Formats	Qty/Box	Cysteine	DMT	TAAcOH	TAAcONa
SiliaPrep SPE Cartridges					
1 mL / 50 mg	100	SPE-R80530B-01B	SPE-R79030B-01B	SPE-R69030B-01B	SPE-R69230B-01B
1 mL / 100 mg	100	SPE-R80530B-01C	SPE-R79030B-01C	SPE-R69030B-01C	SPE-R69230B-01C
3 mL / 200 mg	50	SPE-R80530B-03G	SPE-R79030B-03G	SPE-R69030B-03G	SPE-R69230B-03G
3 mL / 500 mg	50	SPE-R80530B-03P	SPE-R79030B-03P	SPE-R69030B-03P	SPE-R69230B-03P
6 mL / 500 mg	50	SPE-R80530B-06P	SPE-R79030B-06P	SPE-R69030B-06P	SPE-R69230B-06P
6 mL / 1 g	50	SPE-R80530B-06S	SPE-R79030B-06S	SPE-R69030B-06S	SPE-R69230B-06S
6 mL / 2 g	50	SPE-R80530B-06U	SPE-R79030B-06U	SPE-R69030B-06U	SPE-R69230B-06U
12 mL / 2 g	20	SPE-R80530B-12U	SPE-R79030B-12U	SPE-R69030B-12U	SPE-R69230B-12U
25 mL / 5 g*	20	SPE-R80530B-20X	SPE-R79030B-20X	SPE-R69030B-20X	SPE-R69230B-20X
70 mL / 10 g*	16	FLH-R80530B-70Y	FLH-R79030B-70Y	FLH-R69030B-70Y	FLH-R69230B-70Y
70 mL / 15 g*	16	FLH-R80530B-70i	FLH-R79030B-70i	FLH-R69030B-70i	FLH-R69230B-70i
70 mL / 20 g*	16	FLH-R80530B-70Z	FLH-R79030B-70Z	FLH-R69030B-70Z	FLH-R69230B-70Z
150 mL / 25 g*	10	FLH-R80530B-95K	FLH-R79030B-95K	FLH-R69030B-95K	FLH-R69230B-95K
150 mL / 50 g*	10	FLH-R80530B-95M	FLH-R79030B-95M	FLH-R69030B-95M	FLH-R69230B-95M
150 mL / 70 g*	10	FLH-R80530B-95N	FLH-R79030B-95N	FLH-R69030B-95N	FLH-R69230B-95N
276 mL / 100 g*	12	FLH-R80530B-276F	FLH-R79030B-276F	FLH-R69030B-276F	FLH-R69230B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL / 200 mg	50	SPC-R80530B-10G	SPC-R79030B-10G	SPC-R69030B-10G	SPC-R69230B-10G
10 mL / 500 mg	50	SPC-R85030B-10P	SPC-R79030B-10P	SPC-R69030B-10P	SPC-R69230B-10P
Mini-SiliaPrep SPE Cartridges					
500 mg	50	SPS-R80530B-P	SPS-R79030B-P	SPS-R69030B-P	SPS-R69230B-P
1 g	50	SPS-R85030B-S	SPS-R79030B-S	SPS-R69030B-S	SPS-R69230B-S

* Commercialized under SiliaSep OT branding.

Note: Custom formats available on request.






SiliaPrep Metal Scavenger Phases SPE Formats

Formats	Qty/Box	Thiol	Thiourea	Imidazole	Triamine
SiliaPrep SPE Cartridges					
1 mL / 50 mg	100	SPE-R51030B-01B	SPE-R69530B-01B	SPE-R79230B-01B	SPE-R48030B-01B
1 mL / 100 mg	100	SPE-R51030B-01C	SPE-R69530B-01C	SPE-R79230B-01C	SPE-R48030B-01C
3 mL / 200 mg	50	SPE-R51030B-03G	SPE-R69530B-03G	SPE-R79230B-03G	SPE-R48030B-03G
3 mL / 500 mg	50	SPE-R51030B-03P	SPE-R69530B-03P	SPE-R79230B-03P	SPE-R48030B-03P
6 mL / 500 mg	50	SPE-R51030B-06P	SPE-R69530B-06P	SPE-R79230B-06P	SPE-R48030B-06P
6 mL / 1 g	50	SPE-R51030B-06S	SPE-R69530B-06S	SPE-R79230B-06S	SPE-R48030B-06S
6 mL / 2 g	50	SPE-R51030B-06U	SPE-R69530B-06U	SPE-R79230B-06U	SPE-R48030B-06U
12 mL / 2 g	20	SPE-R51030B-12U	SPE-R69530B-12U	SPE-R79230B-12U	SPE-R48030B-12U
25 mL / 5 g*	20	SPE-R51030B-20X	SPE-R69530B-20X	SPE-R79230B-20X	SPE-R48030B-20X
70 mL / 10 g*	16	FLH-R51030B-70Y	FLH-R69530B-70Y	FLH-R79230B-70Y	FLH-R48030B-70Y
70 mL / 15 g*	16	FLH-R51030B-70i	FLH-R69530B-70i	FLH-R79230B-70i	FLH-R48030B-70i
70 mL / 20 g*	16	FLH-R51030B-70Z	FLH-R69530B-70Z	FLH-R79230B-70Z	FLH-R48030B-70Z
150 mL / 25 g*	10	FLH-R51030B-95K	FLH-R69530B-95K	FLH-R79230B-95K	FLH-R48030B-95K
150 mL / 50 g*	10	FLH-R51030B-95M	FLH-R69530B-95M	FLH-R79230B-95M	FLH-R48030B-95M
150 mL / 70 g*	10	FLH-R51030B-95N	FLH-R69530B-95N	FLH-R79230B-95N	FLH-R48030B-95N
276 mL / 100 g*	12	FLH-R51030B-276F	FLH-R69530B-276F	FLH-R79230B-276F	FLH-R48030B-276F
SiliaPrep Large Reservoir Volume SPE Cartridges					
10 mL / 200 mg	50	SPC-R51030B-10G	SPC-R69530B-10G	SPC-R79230B-10G	SPC-R48030B-10G
10 mL / 500 mg	50	SPC-R51030B-10P	SPC-R69530B-10P	SPC-R79230B-10P	SPC-R48030B-10P
Mini-SiliaPrep SPE Cartridges					
500 mg	50	SPS-R51030B-P	SPS-R69530B-P	SPS-R79230B-P	SPS-R48030B-P
1 g	50	SPS-R51030B-S	SPS-R69530B-S	SPS-R79230B-S	SPS-R48030B-S

* Commercialized under SiliaSep OT branding.

Note: Custom formats available on request.



SiliaPrep™ and SiliaPrepX™ Applications

Forensics

Opioids	19
Antidepressants.....	20
Pharmaceutical drugs.....	21
Anti-inflammatory drugs.....	22
Performance-enhancing drugs	23
Cannabinoids.....	24
Stimulants.....	25
Varia	26

Food

Orange & apple juices	27
Milk & honey	28
Salmon & beef.....	29
Potatoes	30

Environment

Pesticides in water.....	31
Bisphenol A in water	34
Drugs in water	35
Other analytes in water.....	36
Other matrices.....	38





Extraction of Methadone from Human Urine and Serum

CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G
SAMPLE PRETREATMENT	200 µL of Phosphoric Acid 2 % was added to 1 mL of urine / serum sample
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	6 mL of Hydrochloric Acid 0.1N then 6 mL of Methanol, dry the cartridge
ELUTION STEP	2 x 3 mL of 10 % Ammonia in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 1 µg/mL
	Methadone in urine 90 %
	Methadone in serum 95 %

Extraction of Methadone and EDDP from Human Urine

CARTRIDGE	SiliaPrepX HLB 1 mL / 30 mg Part Number: SPE-P0002-01AA
SAMPLE PRETREATMENT	40 µL of internal standard (<i>Methadone-d9</i> at 200 ng/mL in Methanol) was added to 200 µL of urine sample and 200 µL of Ammonium Hydroxide 4 %
CONDITIONNING STEP	1 mL of Methanol
EQUILIBRATION STEP	1 mL of Ammonium Hydroxide 2 %
LOADING STEP	Urine sample was slowly aspirated through the cartridge
WASHING STEP	1 mL of Methanol / Ammonium Hydroxide 2 % (50:50) then 1 mL of Methanol / Ammonium Hydroxide 2 % (80:20)
ELUTION STEP	1 mL of Methanol / water (80:20)
FURTHER TREATMENT	Quantification by LDTD/MS/MS (<i>collaboration with Phytronix</i>)
RECOVERY	at 1,000 ng/mL
	Methadone 91 %
	EDDP 85 %

Extraction of Fentanyl and Norfentanyl from Urine

CARTRIDGE	SiliaPrep CleanDRUG 1 mL / 100 mg Part Number: SPEC-R651230B-01C
SAMPLE PRETREATMENT	200 µL of urine was added to 600 µL of Sodium Acetate in water and 40 µL of internal standard (200 ng/mL in Methanol)
CONDITIONNING STEP	1 mL of Methanol
EQUILIBRATION STEP	1 mL of water and 1 mL of Sodium Acetate in water (100 mM, pH 6.0)
LOADING STEP	Urine sample was slowly aspirated through the cartridge
WASHING STEP	1 mL of water then 1 mL of Methanol
ELUTION STEP	1 mL of Ethyl Acetate / Isopropanol / Ammonium Hydroxide (78:20:2)
FURTHER TREATMENT	Evaporation, reconstitution and quantification by LDTD/MS/MS (<i>collaboration with Phytronix</i>)
RECOVERY	at 500 ng/mL
	Fentanyl 96 %
	Norfentanyl 98 %

Extraction of Codeine from Human Urine and Serum

CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G
SAMPLE PRETREATMENT	200 µL of Phosphoric Acid 2 % was added to 1 mL of urine / serum sample
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	6 mL of Hydrochloric Acid 0.1N then 6 mL of Methanol, dry the cartridge
ELUTION STEP	2 x 3 mL of 5 % Ammonia in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS
RECOVERY	at 1 µg/mL
	Codeine in urine 70 %
	Codeine in serum 92 %



Extraction of Tricyclic Antidepressants from Serum

CARTRIDGE	SiliaPrepX WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB
SAMPLE PRETREATMENT	250 µL of serum were diluted with 1 mL of 10 % Formic Acid in water
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	1 mL of 5 % Formic Acid in water then 1 mL Methanol, dry the cartridge
ELUTION STEP	3 mL of 5 % Formic Acid in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS
RECOVERY	at 1 µg/mL
	Doxepine 79 %
	Imipramine 79 %
	Amitriptyline 91 %
	Trimipramine 98 %

Extraction of Pharmaceutical Drugs from Serum

CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G
SAMPLE PRETREATMENT	200 µL of Phosphoric Acid 2 % was added to 1 mL of serum sample
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	6 mL of Chlorhydric Acid 0.1N, dry the cartridge
ELUTION STEP	2 x 3 mL of Methanol (<i>acidic and neutrals analytes</i>) and 2 x 3 mL of 10 % Ammonia in Methanol (<i>basic analytes</i>)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS
RECOVERY	at 1 µg/mL
	Indomethacin 33 %
	Tolmetin 73 %
	Hexobarbital 80 %
	Naproxen 85 %
	Suprofen 108 %
	Phenobarbital 108 %
	Trimipramine 92 %
	Amitriptyline 94 %
	Imipramine 95 %
	Doxepin 101 %

Ropinirole & Amitriptyline Detection in Human Plasma

CARTRIDGE	SiliaPrep CleanDRUG 3 mL / 500 mg Part Number: SPEC-R651230B-03P
SAMPLE PRETREATMENT	Mix 0.1 mL of plasma with 0.1 mL of Methanol and water (50:50) and 2 mL of 1 % Acetic Acid
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of water
LOADING STEP	Plasma sample was slowly aspirated through the cartridge
WASHING STEP	3 mL of water then 3 mL of Methanol
ELUTION STEP	3 mL of 5 % Ammonium Hydroxide in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 10 ng/mL
	Ropinirole 94 %
	Amitriptyline 90 %

Extraction of Pharmaceutical Drugs from Serum

CARTRIDGE	SiliaPrepX SAX 6 mL / 200 mg Part Number: SPE-P0010-06G
SAMPLE PRETREATMENT	pH of serum sample was adjusted to basic value with Sodium Hydroxide 1N
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	6 mL of water, dry the cartridge
ELUTION STEP	2 x 3 mL of Methanol (<i>basic analytes</i>) and 2 x 3 mL of Formic Acid 10 % in Methanol (<i>acidic analytes</i>)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS
RECOVERY	at 1 µg/mL
	Nortriptyline 69 %
	Doxepine 72 %
	Trimipramine 73 %
	Protriptyline 75 %
	Amitriptyline 78 %
	Imipramine 80 %
	Tolmetin 85 %
	Naproxen 86 %
	Suprofen 96 %



Extraction of Tricyclic Antidepressants from Serum

CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G	
CONDITIONING STEP	5 mL of Methanol	
EQUILIBRATION STEP	5 mL of water	
LOADING STEP	1 mL of serum sample (<i>pH value adjusted with 25 µL of Phosphoric Acid</i>) was slowly aspirated through the cartridge	
WASHING STEP	5 mL of water, dry the cartridge	
ELUTION STEP	2 x 3 mL of Methanol	
FURTHER TREATMENT	Quantification by LC/MS	
RECOVERY	Protriptyline	80 %
	Nortriptyline	75 %
	Doxepine	91 %
	Imipramine	88 %
	Amitriptyline	88 %
	Trimipramine	88 %

Extraction of Barbiturates from Serum

CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G	
CONDITIONING STEP	6 mL of Methanol	
EQUILIBRATION STEP	6 mL of water	
LOADING STEP	1 mL of serum sample was slowly aspirated through the cartridge	
WASHING STEP	6 mL of water, dry the cartridge	
ELUTION STEP	6 x 1 mL of Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS	
RECOVERY	at 100 ng/ml	
	Phenobarbital	99 %
	Pentobarbital	69 %
	Hexobarbital	86 %

Extraction of Antibacterial Drugs from Serum

CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G	
SAMPLE PRETREATMENT	Mix 0.1 mL of plasma with 2 mL of 1 % Acetic Acid	
CONDITIONING STEP	6 mL of Methanol	
EQUILIBRATION STEP	6 mL of water	
LOADING STEP	1 mL of serum sample was slowly aspirated through the cartridge	
WASHING STEP	6 mL of water, dry the cartridge	
ELUTION STEP	2 x 3 mL of Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / Water and quantification by LC/MS	
RECOVERY	Cinoxacin	100 %
	Penicillin G	94 %
	Penicillin V	90 %
	Cloxacillin	88 %

Sibutramine Detection in Human Plasma

CARTRIDGE	SiliaPrep CleanDRUG 3 mL / 500 mg Part Number: SPEC-R651230B-03P	
SAMPLE PRETREATMENT	Mix 0.1 mL of plasma with 2 mL of 1 % Acetic Acid	
CONDITIONING STEP	3 mL of Methanol	
EQUILIBRATION STEP	3 mL of water	
LOADING STEP	Plasma sample was slowly aspirated through the cartridge	
WASHING STEP	3 mL of water then 3 mL of Methanol	
ELUTION STEP	3 mL of 5 % Ammonium Hydroxide in Methanol	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS	
RECOVERY	at 5 ng/mL: 82 %	



Extraction of Acidic Pharmaceuticals from Serum

CARTRIDGE	SiliaPrepX SAX 6 mL / 200 mg Part Number: SPE-P0010-06G			
SAMPLE PRETREATMENT	pH of serum sample was adjusted to basic value with Sodium Hydroxide 1N			
CONDITIONING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of water			
LOADING STEP	Treated sample was slowly aspirated through the cartridge			
WASHING STEP	6 mL of water, then 6 mL of Sodium Hydroxide 0.1N and 6 mL of Methanol, dry the cartridge			
ELUTION STEP	6 mL of 1 % Formic Acid in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS			
RECOVERY	at 1 µg/mL			
	Carprofen	69 %	Diclofenac	95 %
	Ibuprofen	88 %	Fenoprofen	98 %
	Ketoprofen	90 %	Fenoprop	104 %
	Meclofenamic Acid	92 %	Flurbiprofen	106 %



Extraction of Anti-inflammatory Drugs From Serum

CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G			
CONDITIONING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of water			
LOADING STEP	1 mL of serum sample (<i>pH value adjusted with 25 µL of Phosphoric Acid</i>) was slowly aspirated through the cartridge			
WASHING STEP	6 mL of 5 % Metanol in water, dry the cartridge			
ELUTION STEP	2 x 3 mL of Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS			
RECOVERY	Suprofen	89 %	Naproxen	87 %
	Tolmetin	89 %	Flurbiprofen	87 %
	Sulindac	84 %	Indomethazin	85 %
	Piroxicam	86 %	Acetyl Salicylic Acid	72 %



Determination of Clenbuterol in Human Plasma	
CARTRIDGE	SiliaPrep CleanDRUG 1 mL / 100 mg Part Number: SPEC-R651230B-01C
SAMPLE PRETREATMENT	50 µL of internal standard (<i>Clenbuterol-d9 at 20 ng/mL in Methanol</i>) was added to 500 µL of plasma and 500 µL of Sodium Acetate (100 mM, pH 6)
CONDITIONNING STEP	1 mL of Methanol
EQUILIBRATION STEP	1 mL of water and 1 mL of Sodium Acetate (100 mM, pH 6)
LOADING STEP	Plasma sample was slowly aspirated through the cartridge
WASHING STEP	1 mL of water, then 1 mL of Acetic Acid 1M and 2 x 1 mL of Methanol
ELUTION STEP	1 mL of Ethyl Acetate/Isopropanol/Ammonium Hydroxide (78:20:2)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LDTD/MS/MS (<i>collaboration with Phytronix</i>)
RECOVERY	at 100 pg/mL: 94 %

Extraction of Atenolol from Human Urine	
CARTRIDGE	SiliaPrepX WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB
SAMPLE PRETREATMENT	9 mL of urine was mixed with 1 mL of a solution of Atenolol in Methanol / Water (10:90)
CONDITIONNING STEP	2 mL of Methanol
EQUILIBRATION STEP	2 mL of water
LOADING STEP	1 mL of sample solution was slowly aspirated through the cartridge
WASHING STEP	2 mL of Monopotassium Phosphate 25 mM (pH 5) then 2 mL of Methanol, dry the cartridge
ELUTION STEP	2 mL of 2 % Formic Acid in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 10 µg/mL: 90 %

Determination of Testosterone in Human Urine	
CARTRIDGE	Mini-SiliaPrep C18 WPD 500 mg Part Number: SPS-R33229G-P
CONDITIONNING STEP	5 mL of Methanol
EQUILIBRATION STEP	5 mL of water
LOADING STEP	2mL of urine sample was slowly aspirated through the cartridge
WASHING STEP	5 mL of water then 5 mL of Hexane
ELUTION STEP	5 mL of Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, derivatization using Girard-P and quantification by LC/MS/MS
RECOVERY	at 250 ng/mL: 95 %

Extraction of Steroids from Serum									
CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G								
CONDITIONNING STEP	5 mL of Methanol								
EQUILIBRATION STEP	5 mL of water								
LOADING STEP	1 mL of serum sample was slowly aspirated through the cartridge								
WASHING STEP	5 mL of 5 % Metanol in water, dry the cartridge								
ELUTION STEP	2 x 3 mL of Methanol								
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS								
RECOVERY	<table border="0"> <tr> <td>Methyl-6a-hydroxy-11β-progesterone</td> <td>89 %</td> </tr> <tr> <td>Methyl-6a-hydroxy-17a-progesterone</td> <td>86 %</td> </tr> <tr> <td>Methyl-6a-hydroxy-17a-progesterone acetate</td> <td>84 %</td> </tr> <tr> <td>Hydrocortisone-21-acetate</td> <td>31 %</td> </tr> </table>	Methyl-6a-hydroxy-11β-progesterone	89 %	Methyl-6a-hydroxy-17a-progesterone	86 %	Methyl-6a-hydroxy-17a-progesterone acetate	84 %	Hydrocortisone-21-acetate	31 %
Methyl-6a-hydroxy-11β-progesterone	89 %								
Methyl-6a-hydroxy-17a-progesterone	86 %								
Methyl-6a-hydroxy-17a-progesterone acetate	84 %								
Hydrocortisone-21-acetate	31 %								



Isolation of Synthetic Cannabinoid Metabolites from Urine

CARTRIDGE	SiliaPrep CleanDRUG 1 mL / 30 mg Part Number: SPEC-R651230B-03G					
SAMPLE PRETREATMENT	1 mL of synthetic urine was spiked with the metabolites and deuterated internal standard, then diluted with 2 mL of a Phosphate buffer solution (pH 6)					
CONDITIONING STEP	3 mL of Methanol					
EQUILIBRATION STEP	3 mL of water and 1 mL of Phosphate buffer					
LOADING STEP	Urine sample was slowly aspirated through the cartridge					
WASHING STEP	3 mL of water then 3 mL of Phosphate buffer / Acetonitrile (80:20)					
ELUTION STEP	6 mL of Ethyl Acetate / Methanol (90:10)					
FURTHER TREATMENT	Evaporation under Nitrogen, derivatization using BSTFA and TMCS, and quantification by GC/MS					
RECOVERY	at 1,000 ng/mL	SiliaPrepX Clean DRUG	Bond Elut® Certify II	HyperSep™ Verify AX	Clean Screen® CSTHC	
		JWH-018	102 %	109 %	112 %	97 %
		JWH-122	96 %	72 %	111 %	80 %
		JWH-250	101 %	71 %	118 %	89 %
CONCLUSION	Our SiliaPrep CleanDRUG performs as well as competitive products to extract cannabinoid metabolites from urine.					

Source: Thesis "An Evaluation of Commercially Available Solid Phase Extraction Cartridges for the Isolation of Synthetic Cannabinoid Metabolites from urine", by Amanda Marie Forni, B.S., Ohio University, 2011



Detection of Δ⁹-Tetrahydrocannabinol in Human Plasma

CARTRIDGE	SiliaPrep CleanENVI 3 mL / 500 mg Part Number: SPEC-R31930B-03P				
SAMPLE PRETREATMENT	250 µL of plasma was added to 1 mL Phosphate buffer (0.1M, pH 6.0)				
CONDITIONING STEP	3 mL of Methanol, then 3 mL of Hydrochloric Acid 1M and 3 mL of water				
EQUILIBRATION STEP	5 mL of water				
LOADING STEP	Plasma sample was slowly aspirated through the cartridge				
WASHING STEP	2 mL of water, then 1 mL of Acetic Acid 1M and 2 mL of 20 % Methanol in water				
ELUTION STEP	5 mL of Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, derivatization using Dansyl Chloride, liquid-liquid extraction, centrifugation, evaporation under Nitrogen, reconstitution with Formic Acid / Acetone and quantification by LC/MS				
RECOVERY	at 2 ng/mL				
		THC	80 %		
		THC-COOH	99 %		
		THC-OH	92 %		



Extraction of Phencyclidine (PCP) from Human Urine

CARTRIDGE	SiliaPrepX HLB 1 mL / 30 mg Part Number: SPE-P0002-01AA
SAMPLE PRETREATMENT	40 µL of internal standard (<i>PCP-d5 at 200 ng/mL in Methanol</i>) was added to 200 µL of urine sample and 200 µL of Ammonium Hydroxide 4 %
CONDITIONNING STEP	1 mL of Methanol
EQUILIBRATION STEP	1 mL of Ammonium Hydroxide 2 %
LOADING STEP	Urine sample was slowly aspirated through the cartridge
WASHING STEP	1 mL of Methanol / Ammonium Hydroxide 2 % (50:50) then 1 mL of Methanol / Ammonium Hydroxide 2 % (80:20)
ELUTION STEP	1 mL of Methanol / Hydrochloric Acid 0.02N (80:20)
FURTHER TREATMENT	Quantification by LDTD/MS/MS (collaboration with Phytronix)
RECOVERY	at 25 ng/mL: 99 %

Drugs of Abuse Determination in Human Urine

CARTRIDGE	SiliaPrep CleanDRUG 3 mL / 200 mg Part Number: SPEC-R651230B-03G
SAMPLE PRETREATMENT	0.5 mL of urine sample was mixed with 2.5 mL Sulfuric Acid 0.1M
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Sulfuric Acid 0.1M
LOADING STEP	2 mL of urine sample was slowly aspirated through the cartridge
WASHING STEP	3 mL of Phosphate buffer (pH 7), then 3 mL of Sulfuric Acid 0.1M and 3 mL of Methanol
ELUTION STEP	2 x 3 mL of Ammonium Hydroxide (5 % in Methanol)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Methanol and quantification by LC/MS
RECOVERY	at 25 ng/mL
	MDMA 92 %
	MDEA 89 %
	Amphetamine 82 %

Amphetamine Quantification in Human Urine

CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB				
SAMPLE PRETREATMENT	100 µL of TFA was added to 10 mL of urine				
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of water				
LOADING STEP	1 mL of urine sample was slowly aspirated through the cartridge				
WASHING STEP	3 mL of (5:95) Methanol / water with 2 % Ammonium Hydroxide; then 3 mL of (20:80) Methanol / water with 2 % Ammonium Hydroxide and 1 mL of (80:20) Methanol / water				
ELUTION STEP	3 mL of Methanol then 3 mL of 2 % Formic Acid in Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water (70:30) and quantification by LC/MS				
RECOVERY	at 100 ng/mL	SiliaPrepX HLB	Bond Elut® Plexa	Oasis® HLB	Strata™-X
	Amphetamine	91 %	88 %	75 %	87 %
	MDA	86 %	86 %	91 %	98 %
	MDEA	95 %	97 %	90 %	101 %
	MDMA	92 %	94 %	91 %	101 %
	Methamphetamine	92 %	95 %	86 %	101 %
	Phentermine	99 %	93 %	90 %	97 %
CONCLUSION	SiliaPrepX HLB allows to extract amphetamines from urine with recoveries as good as competitive products.				



Extraction of Camphorsulfonic Acid from Serum

CARTRIDGE	SiliaPrepX WAX 3 mL / 60 mg Part Number: SPE-P0020-03BB
SAMPLE PRETREATMENT	5 mL of serum was mixed with 5 µL of a solution of Camphorsulfonic Acid (0.5 mg/mL) and 5 mL of Phosphoric Acid 4 %
CONDITIONNING STEP	2 mL of Methanol
EQUILIBRATION STEP	2 mL of water
LOADING STEP	2 mL of sample solution was slowly aspirated through the cartridge
WASHING STEP	2 mL of 2 % Formic Acid in water then 2 mL of Methanol, dry the cartridge
ELUTION STEP	2 mL of 5 % Ammonia in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 0.25 µg/mL: 99 %

Extraction of Alkaloids from Serum

CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	1 mL of serum sample was slowly aspirated through the cartridge
WASHING STEP	6 mL of Methanol, dry the cartridge
ELUTION STEP	2 x 3 mL of Acetone
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 2 µg/mL
	Atropine 99 %
	Papaverine 97 %
	Noscapine 95 %
	Strychnine 94 %
	Quinine 60 %

Extraction of Caffeine, Cotinine & Nicotine from Human Urine

CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB				
SAMPLE PRETREATMENT	500 µL of urine was mixed with 1.5 mL of Sodium Hydroxide 0.1M				
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of water				
LOADING STEP	1 mL of urine sample was slowly aspirated through the cartridge				
WASHING STEP	3 mL of water and dry the cartridge				
ELUTION STEP	3 mL of Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS				
RECOVERY	at 100 ng/mL	SiliaPrepX HLB	Bond Elut® Plexa	Oasis® HLB	Strata™-X
	Caffeine	97 %	99 %	96 %	97 %
	Cotinine	99 %	100 %	98 %	99 %
	Nicotine	89 %	86 %	90 %	89 %
CONCLUSION	SiliaPrepX HLB is as efficient as competitive products to extract caffeine, cotinine and nicotine from urine.				



Extraction of Fungicides in Apple Juice					
CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G				
SAMPLE PRETREATMENT	0.5 mL of Sodium Hydroxide 0.1N was added to 5 mL of apple juice				
CONDITIONNING STEP	6 mL of Methanol				
EQUILIBRATION STEP	6 mL of Ammonia 2 %				
LOADING STEP	Treated sample was slowly aspirated through the cartridge				
WASHING STEP	3 mL of Ammonia 2 %, 3 mL of 30 % Methanol in Ammonia 5 %, 3 mL of Hydrochloric Acid 0.1N and 3 mL of Methanol, dry the cartridge				
ELUTION STEP	6 mL of 30 % Methanol in Ammonia 5 %				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Methanol and quantification by LC/MS				
RECOVERY	at 1 µg/mL				
	<table border="1"> <tr> <td>Carbendazime</td> <td>89 %</td> </tr> <tr> <td>Thiabendazole</td> <td>92 %</td> </tr> </table>	Carbendazime	89 %	Thiabendazole	92 %
Carbendazime	89 %				
Thiabendazole	92 %				

Extraction of Patulin from Apple Juice	
CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB
SAMPLE PRETREATMENT	100 µL of internal standard (250 µg/mL of Patulin-13C (3) in water) and 75 µL of Pectinase Aspergillus Aculeatus were added to 9 mL of apple juice, centrifugate at 3,000 rpm
CONDITIONNING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of water
LOADING STEP	2 mL of sample supernatant was slowly aspirated through the cartridge
WASHING STEP	3 mL of 1 % Sodium Bicarbonate and 1 mL of 0.1 % Acetic Acid, dry the cartridge
ELUTION STEP	2 x 1.5 mL of Ethyl Acetate
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 150 ng/kg: 85 %

Determination of Carbendazim in Orange Juice									
CARTRIDGE	SiliaPrepX SCX 3 mL / 60 mg Part Number: SPE-P0005-03BB								
SAMPLE PRETREATMENT	Centrifugate 5 mL of orange juice 5 min at 3,000 rpm. Sample 1 mL of the supernatant. Add 2 mL of Acetic Acid 10 % and vortex 1 min								
CONDITIONNING STEP	3 mL of Methanol								
EQUILIBRATION STEP	3 mL of Acetic Acid 10 %								
LOADING STEP	3 mL of the treated sample was slowly aspirated through the cartridge								
WASHING STEP	2 mL of Acetic Acid 10 % then 2 mL of Methanol								
ELUTION STEP	3 mL of 5 % Ammonium Hydroxide in Methanol								
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Methanol and quantification by LC/MS								
RECOVERY	at 100 ng/mL								
	<table border="1"> <tr> <td>SiliaPrepX SCX</td> <td>93 %</td> </tr> <tr> <td>Bond Elut® Plexa PCX</td> <td>92 %</td> </tr> <tr> <td>Oasis® MCX</td> <td>92 %</td> </tr> <tr> <td>Strata™-X-C</td> <td>91 %</td> </tr> </table>	SiliaPrepX SCX	93 %	Bond Elut® Plexa PCX	92 %	Oasis® MCX	92 %	Strata™-X-C	91 %
SiliaPrepX SCX	93 %								
Bond Elut® Plexa PCX	92 %								
Oasis® MCX	92 %								
Strata™-X-C	91 %								
CONCLUSION	SiliaPrepX SCX performs as well as Waters, Phenomenex & Agilent products to extract carbendazim from orange juice.								



FOOD APPLICATIONS

Enrichment of Streptomycin in Honey	
CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G
SAMPLE PRETREATMENT	Add 2 g of honey to 8 mL of phosphate buffer (pH 2), filter, dilute to 16 mL (with the same phosphate buffer) and adjust pH value to 7.5
CONDITIONING STEP	5 mL of Methanol
EQUILIBRATION STEP	3 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	5 mL of water, dry the cartridge
ELUTION STEP	5 mL of 3 % Formic Acid in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Acetonitrile and identification by LC/UV
RECOVERY	at 10 µg/kg: 30 %

Extraction of Melamine from Milk	
CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G
SAMPLE PRETREATMENT	1 mL of Hydrochloric Acid 1N was added to 10 mL of milk sample, then mixed with 10 mL of Methylene Chloride. After 15 min centrifugation, remove aqueous layer and extract again organic layer 2 times with 5 mL of Hydrochloric Acid 0.1N. Combine the 3 aqueous fractions.
CONDITIONING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	Combined aqueous fractions were slowly aspirated through the cartridge
WASHING STEP	6 mL of Hydrochloric Acid 0.1N then 6 mL of Methanol, dry the cartridge
ELUTION STEP	2 x 6 mL of 30 % Methanol in Ammonia 5 %
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Methanol and quantification by LC/MS
RECOVERY	at 1 µg/mL: 99 %

Sulfonamides, Tetracyclines & Pyrimethamine Determination in Milk						
CARTRIDGES	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB	OR	SiliaPrepX DVB 3 mL / 60 mg Part Number: SPE-P0001-03BB			
SAMPLE PRETREATMENT	Vortex 2 min 600 µL of bovine milk with 250 µL of 20 % Trichloroacetic Acid in water. Add 2.5 mL of McIlvain buffer (vortex 3 min). Adjust pH of the solution at 5.5 with 1M Sodium Hydroxide. Centrifugate at 3,000 rpm for 5 min.					
CONDITIONING STEP	3 mL of Methanol					
EQUILIBRATION STEP	3 mL of water					
LOADING STEP	1 mL of the treated sample was slowly aspirated through the cartridge					
WASHING STEP	2 x 3 mL of 10 % Methanol in Ammonium Acetate buffer (pH 5.5), dry the cartridge					
ELUTION STEP	3 mL of Methanol					
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS					
RECOVERY	at 1,000 pg/mL	SiliaPrepX HLB	SiliaPrepX DVB	Bond Elut® Plexa	Oasis® HLB	Strata®-X
	Sulfathiazol	84 %	83 %	85 %	83 %	86 %
	Sulfadiazine	90 %	89 %	88 %	87 %	85 %
	Sulfamethoxypyridazine	87 %	89 %	85 %	83 %	87 %
	Sulfamethazole	88 %	84 %	87 %	89 %	82 %
	Sulfamethazine	83 %	84 %	86 %	86 %	84 %
	Pyrimethamine	90 %	90 %	91 %	89 %	86 %
	Tetracycline	96 %	96 %	95 %	84 %	88 %
	Oxytetracycline	96 %	96 %	93 %	80 %	87 %
CONCLUSION	SiliaPrepX HLB and DVB are both equivalent to competitive products to extract sulfonamides, tetracyclines and pyrimethamine from milk.					



Extraction of Marbofloxacin & Sarafloxacin from Salmon

CARTRIDGE	SiliaPrepX SCX 3 mL / 60 mg Part Number: SPE-P0005-03BB				
SAMPLE PRETREATMENT	Add 2 g of salmon and 15 mL of 3 % H ₃ PO ₄ aqueous solution in a 50 mL tube. Shake the tube in a horizontal position for 15 min. Add 5 mL of hexane and vortex for 2 min. Centrifugate at 3,000 rpm for 5 min. Recuperate the aqueous phase from the gelled organic phase by filtration.				
CONDITIONING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of Hydrochloric Acid 1M and 3 mL of water				
LOADING STEP	3 mL of the filtered sample was slowly aspirated through the cartridge				
WASHING STEP	2 mL of Hydrochloric Acid 2M then 1 mL of Methanol				
ELUTION STEP	3 mL of 10 % Ammonium Hydroxide in Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Methanol and quantification by LC/MS				
RECOVERY	at 50 ppb <table border="1"> <tr> <td>Marbofloxacin</td> <td>97 %</td> </tr> <tr> <td>Sarafloxacin</td> <td>87 %</td> </tr> </table>	Marbofloxacin	97 %	Sarafloxacin	87 %
Marbofloxacin	97 %				
Sarafloxacin	87 %				



Extraction of Clenbuterol and Ractopamine from Beef

CARTRIDGE	SiliaPrepX WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB				
SAMPLE PRETREATMENT	100 µL of internal standard (250 µg/mL of Ractopamine d-6 and 250 µg/mL of Clenbuterol-d9 in Methanol) were added to 1g of chopped beef. Add 5 mL of 0.2N Sodium Acetate (pH 5.2) and 50 µL of Beta-Glucuronidase / Arylsulfatase. Add 2.5 mL of 0.1M Perchloric Acid, 2 mL of Phosphoric Acid 4 % in Acetonitrile and 5 mL of 0.5M Glycine (pH 10.5). Adjust to pH 10,50. Add 10 mL of Acetonitrile, 4g of MgSO ₄ and 1g of NaCl. Evaporation and reconstitution with 0.1M Perchloric Acid.				
CONDITIONING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of water				
LOADING STEP	2 mL of treated sample was slowly spirated through the cartridge				
WASHING STEP	1.5 mL of Phosphate buffer 25 mM (pH 7), then 3 mL of water and 1 mL of Methanol				
ELUTION STEP	3 mL of Formic Acid 2 % in Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS				
RECOVERY	at 70 ppb <table border="1"> <tr> <td>Clenbuterol</td> <td>92 %</td> </tr> <tr> <td>Ractopamine</td> <td>91 %</td> </tr> </table>	Clenbuterol	92 %	Ractopamine	91 %
Clenbuterol	92 %				
Ractopamine	91 %				



FOOD APPLICATIONS



Extraction of Glycoalkaloids from Potatoes

CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G
SAMPLE PRETREATMENT	Extract 3 g of potato powder with 20 mL of water / Acetic Acid / Sodium Metabisulfite (95:5:0.5). Centrifugate for 10 min and filter.
CONDITIONNING STEP	5 mL of Acetonitrile
EQUILIBRATION STEP	5 ml of water / Acetic Acid / Sodium Metabisulfite (95:5:0.5)
LOADING STEP	10 mL of treated sample was slowly aspirated through the cartridge
WASHING STEP	4 mL of 0.5 % Ammonium Hydroxide, then 4 mL of water and 4 mL of Acetonitrile / water (15:85). Dry the cartridge.
ELUTION STEP	5 mL of Acetonitrile / Potassium Dihydrogen Phosphate 10mM (60:40), pH 7.6
FURTHER TREATMENT	Qualitative analysis by TLC



Acrylamide Determination in Fried Potato Chips

CARTRIDGES	<u>Step 1:</u> SiliaPrepX HLB 6 mL / 200 mg Part Number: SPE-P0002-06G <u>Step 2:</u> SiliaPrepX SCX 3 mL / 60 mg Part Number: SPE-P0005-03BB
SAMPLE PRETREATMENT	<u>Extraction 1:</u> vortex for 1 min 1g of potato chips and 8 mL of Sodium Chloride aqueous 4M. Incubate 30 min at 60°C (<i>vortex 10 sec every 10 min</i>). Centrifugate for 10 min at 4,500 rpm and collect the supernatant. <u>Extraction 2:</u> repeat previous 3 steps with same potato chips. Add 1 mL of solution Cirraz 1 (15 g of $K_3[Fe(CN)_6]$ in 100 mL water) and 1 mL of solution Cirraz 2 (30 g of $Zn(O_2CCH_3)_2$ in 100 mL water).
CONDITIONNING STEP (1) [SILIAPREP X HLB]	3 mL of Methanol
EQUILIBRATION STEP (1)	3 mL of water
LOADING STEP (1)	1.5 mL of the treated sample was slowly aspirated through the cartridge
WASHING STEP (1)	1.5 mL of water
ELUTION STEP (1)	3 mL of Methanol
CONDITIONNING STEP (2) [SILIAPREP X SCX]	3 mL of Methanol
LOADING STEP (2)	The treated sample eluted from SiliaPrepX HLB was slowly aspirated through the cartridge (<i>collect this fraction</i>)
WASHING STEP (2)	1 mL of Methanol (<i>mix this fraction with the one previously collected</i>)
FURTHER TREATMENT	Evaporation to dryness, reconstitution with water / Methanol and quantification by LC/MS
RECOVERY	at 100 µg/kg: 88 %



Determination of Pesticides in Drinking Water

CARTRIDGE	SiliaPrepX HLB 6 mL / 200 mg Part Number: SPE-P0002-06G					
CONDITIONING STEP	6 mL of Methanol					
EQUILIBRATION STEP	6 mL of water (HPLC grade)					
LOADING STEP	100 mL of drinking water was slowly aspirated through the cartridge					
WASHING STEP	6 mL of water (HPLC grade)					
ELUTION STEP	2 x 6 mL of Methanol					
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol and quantification by LC/MS					
RECOVERY	at 1,000 pg/mL	Atrazine	Benalaxyl	Carbendazim	Chloroxuron	Imazalil
	SiliaPrepX HLB	75 %	76 %	103 %	91 %	78 %
	Oasis® HLB	66 %	48 %	103 %	99 %	78 %
		Methalaxyl	Myclobutanil	Propoxur	Simazine	Thiambazole
	SiliaPrepX HLB	87 %	91 %	70 %	98 %	91 %
	Oasis® HLB	61 %	101 %	42 %	79 %	80 %
CONCLUSION	SiliaPrepX HLB compares favorably with Oasis® HLB for the extraction of 8 pesticides out of 10.					

Pesticides Determination in Drinking Water

CARTRIDGE	SiliaPrep CleanENVI 6 mL / 1 g Part Number: SPEC-R31930B-06S	
CONDITIONING STEP	10 mL of Methanol	
EQUILIBRATION STEP	10 mL of water (HPLC grade)	
LOADING STEP	10 mL of drinking water was slowly aspirated through the cartridge	
WASHING STEP	2 x 5 mL of water (HPLC grade)	
ELUTION STEP	2 x 3 mL of Acetone	
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Methanol and quantification by LC/MS	
RECOVERY	at 50 ng/mL	
	Atrazine	84 %
	Simazine	95 %
	Alachlor	68 %

Pesticides Determination in Water

CARTRIDGE	SiliaPrepX LRV SAX 10 mL / 60 mg Part Number: SPC-P0010-10BB
CONDITIONING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of water (HPLC grade)
LOADING STEP	100 mL of sample water was slowly aspirated through the cartridge
WASHING STEP	3 mL of water (HPLC grade)
ELUTION STEP	3 mL of Methanol then 3 mL of Methanol with Formic Acid 2 %
FURTHER TREATMENT	Quantification by LC/MS/MS
RECOVERY	at 1,000 pg/mL: > 80 % for 23 pesticides



ENVIRONMENT APPLICATIONS

Glyphosate & AMPA Determination in Water

CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB				
SAMPLE PRETREATMENT	Derivatization with FMOC-Cl: to 5 mL of sample water was added 325 µL of Sodium Borate 50mM, 200 µL of EDTA 0.1M, 4.5 mL of Acetonitrile and 0.6 mL of FMOC-Cl 50 mg/mL. Evaporate. Aqueous supernatant was mixed with 2 mL Ethyl Acetate. Adjust pH of the aqueous layer to 3 by adding 100 µL of Formic Acid 5 %.				
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of water and 3 mL of Formic Acid 0.1 %				
LOADING STEP	Derivatized sample was slowly aspirated through the cartridge				
WASHING STEP	1 mL of Formic Acid 0.1 % then 2 x 500 µL of water, dry the cartridge				
ELUTION STEP	3 mL of Methanol				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Acetonitrile and quantification by LC/MS/MS				
RECOVERY	at 5 ng/mL				
	<table border="1"> <tr> <td>Glyphosate</td> <td>120 %</td> </tr> <tr> <td>AMPA</td> <td>106 %</td> </tr> </table>	Glyphosate	120 %	AMPA	106 %
Glyphosate	120 %				
AMPA	106 %				

Diquat & Paraquat Determination in Water

CARTRIDGE	SiliaPrepX WCX 3 mL / 60 mg Part Number: SPE-P0015-03BB				
CONDITIONNING STEP	3 mL of Methanol				
EQUILIBRATION STEP	3 mL of water				
LOADING STEP	100 mL of sample water was slowly aspirated through the cartridge				
WASHING STEP	3 mL of water then 3 mL of Methanol				
ELUTION STEP	2 x 3 mL of Acetonitrile / Isopropanol / Formic Acid (40:40:20)				
FURTHER TREATMENT	Quantification by LC/MS/MS				
RECOVERY	at 10 ppb				
	<table border="1"> <tr> <td>Diquat</td> <td>90 %</td> </tr> <tr> <td>Paraquat</td> <td>90 %</td> </tr> </table>	Diquat	90 %	Paraquat	90 %
Diquat	90 %				
Paraquat	90 %				

Determination of Pesticides in Water (by GC/ECD)

CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB			
CONDITIONNING STEP	3 mL of 30 % Acetone in Toluene then 3 mL of Methanol			
EQUILIBRATION STEP	3 mL of distilled water			
LOADING STEP	100 mL of sample water was slowly aspirated through the cartridge			
WASHING STEP	3 mL of distilled water, dry the cartridge			
ELUTION STEP	500 µL of Acetone, then 2 mL of 30 % Acetone in Toluene and 2.5mL of 30 % Acetone in Toluene			
FURTHER TREATMENT	Quantification by GC/ECD			
RECOVERY	Trifluralin	90 %	Endrin	95 %
	Lindane	88 %	4,4'-DDT	75 %
	Aldrin	78 %	Diclofop-methyl	90 %
	Heptachlor Epoxide	88 %	Methoxychlor	86 %
	Dieldrin	90 %	Chlordane	79 %



Extraction of Desphenyl Chloridazon from Water	
CARTRIDGE	SiliaPrepX SAX 3 mL / 60 mg Part Number: SPE-P0010-03BB
SAMPLE PRETREATMENT	100 µL of Ammonium Hydroxide 26 % was added to 1 mL of water sample
CONDITIONNING STEP	1 mL of Methanol
EQUILIBRATION STEP	1 mL of Ammonium Hydroxide 5 %
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	1 mL of Ammonium Hydroxide 5 % then 1 mL of Methanol
ELUTION STEP	2 x 1 mL of 5 % Formic Acid in Ethyl Acetate
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile / water and quantification by LC/MS
RECOVERY	at 10 µg/mL: 83 %

Quantification of Acidic Herbicides	
CARTRIDGE	SiliaPrepX SAX 6 mL / 200 mg Part Number: SPE-P0010-06G
SAMPLE PRETREATMENT	pH of sample was adjusted to basic value with Sodium Hydroxide 1N
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of water
LOADING STEP	Treated sample was slowly aspirated through the cartridge
WASHING STEP	3 mL of Sodium Acetate then 3 mL of Methanol, dry the cartridge
ELUTION STEP	2 x 3 mL of Formic Acid 10 % in Methanol
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS
RECOVERY	at 1 µg/mL
	Bentazon 79 %
	Dicamba 87 %
	2,4-Dichlorophenoxy Acetic Acid 82 %
	4-Chloro-2-methylphenoxy Acetic Acid 76 %

Isothiazolinone Biocides in an Aqueous Sample	
CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G
SAMPLE PRETREATMENT	5 mL of isothiazolinones standard solution (1 µg/mL) are diluted in 50 mL water and 500 µL Formic Acid. The solution is filled up to 100 mL.
CONDITIONNING STEP	6 mL of Methanol
EQUILIBRATION STEP	6 mL of 0.1 % Formic Acid in water
LOADING STEP	5 mL of sample was slowly aspirated through the cartridge
WASHING STEP	6 mL 0.1 % Formic Acid in water, dry the cartridge
ELUTION STEP	3 mL of Methanol then 6 mL of Acetonitrile
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS
RECOVERY	at 50 ng/L
	Methylisothiazolinone 93 %
	Chloromethylisothiazolinone 96 %
	Benzisothiazolinone 85 %
	Butylbenzisothiazolinone 88 %
	Octylisothiazolinone 90 %
	Dichlorooctylisothiazolinone 83 %



Extraction of Bisphenol A, Triclosan & Ethynyl Estradiol from Water

CARTRIDGE	SiliaPrepX HLB 3 mL / 60 mg Part Number: SPE-P0002-03BB						
SAMPLE PRETREATMENT	To 25 mL of sample water was added 250 µL of internal standard (1 ppb of 17α-Ethynyl Estradiol d-6, 1 ppb of Bisphenol A d-16 and 0.4 ppb of Triclosan d-3 in Methanol)						
CONDITIONING STEP	3 mL of Methanol						
EQUILIBRATION STEP	3 mL of water and 1 mL of Acetic Acid 100mM						
LOADING STEP	Treated sample was slowly aspirated through the cartridge						
WASHING STEP	3 mL of water, 1 mL of Acetic Acid 100 mM and 2 mL of 20 % Methanol in water, dry the cartridge						
ELUTION STEP	2 x 3 mL of Dichloromethane / Acetone (50:50)						
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Sodium Carbonate in water, derivatization with Dansyl Chloride and quantification by LC/MS/MS						
RECOVERY	<table><tr><td>17α-Ethynyl Estradiol</td><td>93 %</td></tr><tr><td>Bisphenol A</td><td>115 %</td></tr><tr><td>Triclosan</td><td>98 %</td></tr></table>	17α-Ethynyl Estradiol	93 %	Bisphenol A	115 %	Triclosan	98 %
17α-Ethynyl Estradiol	93 %						
Bisphenol A	115 %						
Triclosan	98 %						



Analysis of Bisphenol A in Bottled Water

CARTRIDGE	SiliaPrep C18 Plus 6 mL / 200 mg (glass) Part Number: SPE-R00830B-06G
SAMPLE PRETREATMENT	100 µL of internal standard (Bisphenol A-d16 in methanol, 1 µg/mL) was added to 50 mL of bottled water
CONDITIONING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of water (HPLC grade) and 1 mL of Acetic Acid 100 mM
LOADING STEP	The whole sample was aspirated through the cartridge using SiliCycle MiniBlock equipment (2 drops / second)
WASHING STEP	5 mL of water (HPLC grade), dry the cartridge
ELUTION STEP	3 mL of Methanol
FURTHER TREATMENT	Evaporation to dryness, derivatization using Dansyl Chloride, liquid-liquid extraction, evaporation, reconstitution with Methanol and quantification by LC/MS/MS
RECOVERY	at 3,000 pg/mL: 97 %



Pharmaceutical Drugs Determination in Water

CARTRIDGE	Silia <i>PrepX</i> HLB (200 mg) + SAX (60 mg) / 10 mL Part Number: SPC-P0210-10i			
SAMPLE PRETREATMENT	100 mL of sample water was mixed with 5 mL of Sodium Acetate 10 %. pH was adjusted to 9.5 with a buffer solution (NH_4Cl 0.5M and NH_4OH 0.5M in water).			
CONDITIONING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of water and 6 mL of buffer pH 9.5			
LOADING STEP	Treated sample was slowly aspirated through the cartridge			
WASHING STEP	3 mL of buffer pH 9.5 and 3 mL of water, dry the cartridge			
ELUTION STEP	2 mL of Methanol and 2 mL of Formic Acid 2 % in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with water / Acetonitrile and quantification by LC/MS/MS			
RECOVERY	at 100 ppt			
	Trimethoprim	105 %	Caffeine C13	96 %
	Sulphamethoxazole	100 %	Acetaminophen	93 %
	Naproxen	100 %	Norfloracin	70 %
	Ibuprofen	85 %	Maprotiline	79 %
	Carbamazepine	102 %		

Determination of Tricyclic Antidepressants in Water

CARTRIDGE	Silia <i>PrepX</i> DVB 3 mL / 60 mg Part Number: SPE-P0001-03BB		
CONDITIONING STEP	1 mL of Methanol		
EQUILIBRATION STEP	1 mL of water		
LOADING STEP	1 mL of sample was slowly aspirated through the cartridge		
WASHING STEP	1 mL of water		
ELUTION STEP	1 mL of Acetonitrile		
FURTHER TREATMENT	Quantification by LC/MS		
RECOVERY	at 1 µg/mL		
	Protriptyline	93 %	
	Nortriptyline	90 %	



ENVIRONMENT APPLICATIONS

Determination of Explosives in Well Water																																					
CARTRIDGE	SiliaPrepX DVB 6 mL / 200 mg Part Number: SPE-P0001-06G																																				
CONDITIONING STEP	6 mL of Methanol, 6 mL of Acetonitrile																																				
EQUILIBRATION STEP	10 mL of water																																				
LOADING STEP	1 L of well water (with 5 g of Sodium Chloride) was slowly aspirated through the cartridge																																				
WASHING STEP	10 mL of water, DO NOT dry the cartridge																																				
ELUTION STEP	6 of mL Methanol / Acetonitrile (50:50)																																				
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS																																				
RECOVERY	at 1 µg/L																																				
	<table border="1"> <tbody> <tr> <td>Hexanitrodiphenylamine</td> <td>96 %</td> <td>4-Amino-2,6-dinitrotoluene</td> <td>95 %</td> </tr> <tr> <td>Diphenylamine</td> <td>100 %</td> <td>2-Amino-4,6-dinitrotoluene</td> <td>94 %</td> </tr> <tr> <td>Pentaerythritol Tetranitrate</td> <td>108 %</td> <td>2,4,6-Trinitrotoluene</td> <td>92 %</td> </tr> <tr> <td>3-Nitrotoluene</td> <td>78 %</td> <td>Nitroglycerine</td> <td>88 %</td> </tr> <tr> <td>4-Nitrotoluene</td> <td>81 %</td> <td>1,3-Dinitrobenzene</td> <td>86 %</td> </tr> <tr> <td>2-Nitrotoluene</td> <td>67 %</td> <td>1,3,5-Trinitrobenzene</td> <td>96 %</td> </tr> <tr> <td>2,6-Dinitrotoluene</td> <td>94 %</td> <td>Ethylene Glycol Dinitrate</td> <td>95 %</td> </tr> <tr> <td>2,4-dinitrotoluene</td> <td>85 %</td> <td>Picric Acid</td> <td>92 %</td> </tr> <tr> <td>Octogen</td> <td>94 %</td> <td>Diethylene Glycol Dinitrate</td> <td>74 %</td> </tr> </tbody> </table>	Hexanitrodiphenylamine	96 %	4-Amino-2,6-dinitrotoluene	95 %	Diphenylamine	100 %	2-Amino-4,6-dinitrotoluene	94 %	Pentaerythritol Tetranitrate	108 %	2,4,6-Trinitrotoluene	92 %	3-Nitrotoluene	78 %	Nitroglycerine	88 %	4-Nitrotoluene	81 %	1,3-Dinitrobenzene	86 %	2-Nitrotoluene	67 %	1,3,5-Trinitrobenzene	96 %	2,6-Dinitrotoluene	94 %	Ethylene Glycol Dinitrate	95 %	2,4-dinitrotoluene	85 %	Picric Acid	92 %	Octogen	94 %	Diethylene Glycol Dinitrate	74 %
Hexanitrodiphenylamine	96 %	4-Amino-2,6-dinitrotoluene	95 %																																		
Diphenylamine	100 %	2-Amino-4,6-dinitrotoluene	94 %																																		
Pentaerythritol Tetranitrate	108 %	2,4,6-Trinitrotoluene	92 %																																		
3-Nitrotoluene	78 %	Nitroglycerine	88 %																																		
4-Nitrotoluene	81 %	1,3-Dinitrobenzene	86 %																																		
2-Nitrotoluene	67 %	1,3,5-Trinitrobenzene	96 %																																		
2,6-Dinitrotoluene	94 %	Ethylene Glycol Dinitrate	95 %																																		
2,4-dinitrotoluene	85 %	Picric Acid	92 %																																		
Octogen	94 %	Diethylene Glycol Dinitrate	74 %																																		

Determination of Surfactants in Water													
CARTRIDGE	SiliaPrepX WAX 3 mL / 60 mg Part Number: SPE-P0020-03BB												
CONDITIONING STEP	2 mL of 5 % Ammonia in Methanol then 2 mL of Methanol												
EQUILIBRATION STEP	2 mL of water												
LOADING STEP	500 mL of water sample was slowly aspirated through the cartridge												
WASHING STEP	2 mL of water, then 2 mL of Acetone / Acetonitrile / Formic Acid (50:50:1) and 2 mL of Methanol												
ELUTION STEP	2 mL of 5 % Ammonia in Methanol												
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS												
RECOVERY	at 20 µg/L												
	<table border="1"> <tbody> <tr> <td>Perfluorooctane Sulfonate Potassium Salt</td> <td>81 %</td> </tr> <tr> <td>Perfluoropentanoic Acid</td> <td>94 %</td> </tr> <tr> <td>Perfluorohexanoic Acid</td> <td>94 %</td> </tr> <tr> <td>Perfluorooctanoic Acid</td> <td>95 %</td> </tr> <tr> <td>Perfluoropropionic Acid</td> <td>103 %</td> </tr> <tr> <td>Perfluorododecanoic Acid</td> <td>82 %</td> </tr> </tbody> </table>	Perfluorooctane Sulfonate Potassium Salt	81 %	Perfluoropentanoic Acid	94 %	Perfluorohexanoic Acid	94 %	Perfluorooctanoic Acid	95 %	Perfluoropropionic Acid	103 %	Perfluorododecanoic Acid	82 %
Perfluorooctane Sulfonate Potassium Salt	81 %												
Perfluoropentanoic Acid	94 %												
Perfluorohexanoic Acid	94 %												
Perfluorooctanoic Acid	95 %												
Perfluoropropionic Acid	103 %												
Perfluorododecanoic Acid	82 %												

Quantification of Phenolic Acids							
CARTRIDGE	SiliaPrepX SAX 6 mL / 200 mg Part Number: SPE-P0010-06G						
SAMPLE PRETREATMENT	pH of sample was adjusted to basic value with Sodium Hydroxide 1N						
CONDITIONING STEP	6 mL of Methanol						
EQUILIBRATION STEP	6 mL of water						
LOADING STEP	Treated sample was slowly aspirated through the cartridge						
WASHING STEP	3 mL of water, then 3 mL of Sodium Hydroxide 0.1N and 3 mL of Methanol, dry the cartridge						
ELUTION STEP	2 x 3 mL of Formic Acid 5 % in Methanol						
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS						
RECOVERY	at 1 µg/mL						
	<table border="1"> <tbody> <tr> <td>Syringic Acid</td> <td>70 %</td> </tr> <tr> <td>Vanillic Acid</td> <td>86 %</td> </tr> <tr> <td>p-Hydroxybenzoic Acid</td> <td>97 %</td> </tr> </tbody> </table>	Syringic Acid	70 %	Vanillic Acid	86 %	p-Hydroxybenzoic Acid	97 %
Syringic Acid	70 %						
Vanillic Acid	86 %						
p-Hydroxybenzoic Acid	97 %						



Extraction of Amines from an Aqueous Sample

CARTRIDGE	SiliaPrepX SCX 6 mL / 200 mg Part Number: SPE-P0005-06G			
SAMPLE PRETREATMENT	200 µL of Phosphoric Acid 2 % was added to 1 mL of aqueous sample			
CONDITIONING STEP	6 mL of Methanol			
EQUILIBRATION STEP	6 mL of water			
LOADING STEP	Treated sample was slowly aspirated through the cartridge			
WASHING STEP	6 mL of Hydrochloric Acid 0.1N then 6 mL of Methanol, dry the cartridge			
ELUTION STEP	2 x 3 mL of 10 % Ammonia in Methanol			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Methanol / water and quantification by LC/MS			
RECOVERY	at 100 ppm			
	2-Naphthylamine	65 %	4,4'-Methylene-bis-(2-chloro-aniline)	75 %
	Benzidine	104 %	4,4'-Oxydianiline	104 %
	5-Nitro-o-toluidine	80 %	4,4'-Methylenedianiline	109 %
	Xenylamine	89 %	4,4'-Thiodianiline	100 %
	o-Aminoazotoluene	89 %	4,4'-Methylendi-o-toluidine	110 %
	4-Aminoazobenzene	99 %	3,3-Dichlorobenzidine	110 %

Extraction of PAHs from Drinking Water

CARTRIDGE	SiliaPrep PAH 6 mL / 1.5 g Part Number: SP2-R0610030B-06T			
CONDITIONING STEP	5 mL of 2-Propanol			
EQUILIBRATION STEP	5 mL of water / 2-Propanol (92:8)			
LOADING STEP	500 mL of drinking water was slowly aspirated through the cartridge			
WASHING STEP	3 mL of Dichloromethane HPLC grade, soak the sorbent for 10 minutes before eluting. Repeat a second time.			
ELUTION STEP	2 mL of Dichloromethane HPLC grade, soak the sorbent for 10 minutes before eluting. Combine the 3 eluates.			
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetone / water and qualification by HPLC (<i>Fluorescence</i>)			
RECOVERY		Benzo[b]fluoranthene	Benzo[k]fluoranthene	Benzo[a]pyrene
	SiliaPrep PAH	118 %	99 %	94 %
	BAKERBOND PAH Aqua	117 %	102 %	100 %
		Benzo[ghi]perylene	Indeno[1,2,3-cd]pyrene	
	SiliaPrep PAH	117 %	126 %	
	BAKERBOND PAH Aqua	115 %	114 %	
CONCLUSION	SiliaPrep PAH performs as well as BAKERBOND PAH Aqua for the extraction of PAHs from water.			



Analysis of Pesticides in Oats, after a Fatty Acids Cleanup

CARTRIDGE	SiliaPrep Diamine 6 mL / 500 mg Part Number: SPE-R49030B-06P
SAMPLE PRETREATMENT	10 g of oat was added to 100 mL of water and 200 mL of Acetone. 35 g of NaCl and 100 mL of 50 % Ethylacetate in Cyclohexane were added for liquid-liquid extraction. The organic layer (200 mL) was dried with NaSO ₄ , filtered, evaporated and reconstituted with 10 mL of 50 % Ethylacetate in Cyclohexane.
CONDITIONING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Acetone and 3 mL of 50 % Ethyl Acetate in Cyclohexane
LOADING STEP	1 mL of treated sample was slowly aspirated through the cartridge (<i>collect the eluted solvent</i>)
WASHING STEP	6 mL of Hydrochloric Acid 0.1N then 6 mL of Methanol, dry the cartridge
ELUTION STEP	15 mL of 50 % Ethyl Acetate in Cyclohexane (<i>mix with the fraction previously collected</i>)
FURTHER TREATMENT	Evaporation under Nitrogen, reconstitution with Acetonitrile, derivatization using HMDS and TFA, and quantification by GC/MS
RECOVERY	> 80 % for 84 pesticides < 1 % for fatty acids

Source: P. Steinbach, W. Schwack, "Comparison of different solid-phase-extraction cartridges for a fatty acid cleanup of the ethyl acetate / cyclohexane based multi-pesticide residue method EN 12393", *J. Chromatogr. A*, **2014**, 1323, 28 - 38



Triacylglycerols Profiling of Marine Microalgae

CARTRIDGE	SiliaPrep Silica 3 mL / 500 mg Part Number: SPE-R10030B-03P
SAMPLE PRETREATMENT	Algal extracts were extracted with Hexane, washed with water and evaporated
CONDITIONING STEP	3 mL of Hexane
EQUILIBRATION STEP	3 mL of distilled water
LOADING STEP	50 mg of lipid sample in 300 µL of Hexane was slowly aspirated through the cartridge
ELUTION STEP	Elution 1 (<i>for triacylglycerols</i>): Hexane / Diethyl Ether / Acetic Acid (80:20:1) Elution 2 (<i>for polar lipids and chlorophyll</i>): Acetone
FURTHER TREATMENT	Evaporation, reconstitution with Hexane and quantification by LC/MS/MS

Source: M. Danielewicz, L. Anderson, A. Franz, "Triacylglycerol profiling of marine microalgae by mass spectrometry", *Journal of Lipid Research*, **2011**, 52, 2101 - 2108



Extraction of Allantoin from a Cosmetic Product

CARTRIDGE	SiliaPrepX SAX 3 mL / 60 mg Part Number: SPE-P0010-03BB
SAMPLE PRETREATMENT	1g of cosmetic was diluted in 100 mL of water, pH was adjusted to 10 with Ammonium Hydroxide 5 %
CONDITIONING STEP	3 mL of Methanol
EQUILIBRATION STEP	3 mL of Ammonium Hydroxide 5 %
LOADING STEP	1 mL of treated sample was slowly aspirated through the cartridge
WASHING STEP	3 mL of Ammonium Hydroxide 5 % then 3 mL of Methanol
ELUTION STEP	2 x 1 mL of Hydrochloric Acid 0.6 %
FURTHER TREATMENT	Add Acetonitrile / Ammonium Chloride 30mM, and qualification by HPLC

Share your complete SiliaPrep / SiliaPrepX application with SiliCycle... ... and get any standard* SPE box for free!

Send us:

- the detailed method (*cartridge used, conditioning, equilibration, loading, washing & elution steps*)
- analytes & matrix
- pre- and post-treatments
- recoveries
- HPLC or GC analysis

* box of 200 cartridges are not available for this promotion

Contact us: sampleprep@silicycle.com





SiliaPrep™ Accessories

Simplify your Solid-Phase Extractions

- Vacuum Manifolds
- Empty Tubes
- Adapters & Vacuum Adapters
- 96-Well Collection Plates
- Phase Separator Cartridges

Maximize your Productivity with SiliaPrep Accessories

SiliCycle offers various accessories for SPE Cartridges and Well Plates to simplify method development and expedite high throughput analysis.

SiliaPrep SPE Vacuum Manifolds

Run multiple samples simultaneously, with a controlled flow rate for higher reproducibility, with SiliaPrep SPE Vacuum Manifolds. These manifolds are available in 12 and 24-Positions configurations and allow consistent extraction. No possibility of cross-contamination from one sample to another.

The design consists in a clear glass chamber equipped with replaceable individual stopcocks (*also known as control valves*) and solvent guide needles. The adjustable rack allows the use of a wide variety of collection vessels including 13 and 16 mm test tubes, autosampler vials and volumetric flasks.

Simply apply a vacuum source to elute sample through a cartridge directly to the collection vessel of choice.

Complete sets include:

- Glass chamber, vacuum gauge & bleed valve
- Cover, gasket, male and female luer fittings
- Individual stopcocks and needles
- Collection rack with posts, shelves and retaining clips.



SiliaPrep SPE Vacuum Manifolds (Complete Sets)

Product Number	Description
AUT-0128-12	12-Positions SiliaPrep SPE Vacuum Manifold
AUT-0129-24	24-Positions SiliaPrep SPE Vacuum Manifold

SiliaPrep Vacuum Manifold Accessories

Various replacement parts are available for the two SiliaPrep Vacuum Manifolds offered by SiliCycle.

SiliaPrep Vacuum Manifold Accessories			
Description	12-Positions Vacuum Manifold	24-Positions Vacuum Manifold	
SiliaPrep Vacuum Manifold Complete Set	AUT-0128-12 (1/box)	AUT-0129-24 (1/box)	
Spare Parts Ordering Information			
INCLUDED	Glass chamber [Dimensions: Length x Width x Height]	AUT-0182-2 (1/box) [7" x 5.25" x 7"]	AUT-0185 (1/box) [12" x 5.25" x 7"]
	Vacuum gauge, valve & glass chamber kit	AUT-0187 (1/box)	AUT-0189 (1/box)
	Top cover gasket	AUT-0174 (2/box)	AUT-0193 (2/box)
	Polypropylene stopcocks	AUT-0146 (12/box)	AUT-0147 (24/box)
	Polypropylene needles	AUT-0154 (12/box)	AUT-0155 (24/box)
	Collection rack kit (posts, shelves and retaining clips included)	AUT-0202 (1/box)	AUT-0204 (1/box)
	Plate for 13 mm test tubes	AUT-0205 (1/box)	AUT-0207 (1/box)
	Plate for 16 mm test tubes	AUT-0208 (1/box)	AUT-0210 (1/box)
	Plate for autosampler vials	AUT-0213 (1/box)	-
Plate for volumetric flasks	AUT-0214 (1/box)	-	

Note: Stainless Steel needles and Teflon® needles are available upon request.

SiliaPrep Waste Containers

Disposable solvent resistant polypropylene containers are available for the 12-Positions manifold. These waste containers greatly simplify sample preparation, solvent disposal and clean-up. Depending on the nature of the solvent used, the waste container can be reused many times prior to discarding.

Note: Waste containers not available for the 24-Positions vacuum manifold.



AUT-0176 (10/box)

SiliaPrep Drying Manifold Covers

SiliaPrep Drying Manifold Covers can be used to concentrate samples with a flow of air or gas (nitrogen).

SiliaPrep Drying Manifold Covers	
Product Number	Description
AUT-0215-12	12-Positions SiliaPrep Drying Manifold Cover (1/box)
AUT-0215-24	24-Positions SiliaPrep Drying Manifold Cover (1/box)



AUT-0215-12

SiliaPrep Adapters

Enable cartridge stacking and easy SPE cartridge connection with syringe or gas lines (for positive pressure).

SiliaPrep Adapters	
Product Number	Description
AUT-0172	SiliaPrep Adapter for 1, 3, 6 & 12 mL SPE (10/box)
AUT-0173	SiliaPrep Adapter for 25 & 70 mL SPE (10/box)



AUT-0172



AUT-0173

SiliaPrep Vacuum Adapters

Fast, user-friendly, and economical adapters for SPE cartridges. Only a vacuum source is needed.

SiliaPrep Vacuum Adapter - Flasks		
Joint	PN	Description
24/40	AUT-0043	24/40 - SiliaPrep Vacuum Adapter (1/box)
19/22	AUT-0044	19/22 - SiliaPrep Vacuum Adapter (1/box)
14/22	AUT-0045	14/22 - SiliaPrep Vacuum Adapter (1/box)

SiliaPrep Vacuum Adapter - Screw Thread Vials		
Thread	PN	Description
22/400	AUT-0046	22/400 Vial - SiliaPrep Vacuum Adapter Without Vial Connector (1/box)
22/400	AUT-0047	22/400 Vial - SiliaPrep Vacuum Adapter With Vial Connector (1/box)



AUT-0043



AUT-0044



AUT-0045



AUT-0046



AUT-0047

SiliaPrep Empty Tubes

You can use our SiliaPrep Empty Tubes to pack your own SPE cartridges with bulk sorbents of your choice.

SiliaPrep Empty Tubes	
Formats	Description
SIM-0007-001	Empty 1 mL SPE tube with 2 frits (100/box)
SIM-0008-003	Empty 3 mL SPE tube with 2 frits (100/box)
SIM-0002-006	Empty 6 mL SPE tube with 2 frits (100/box)
SIM-0003-012	Empty 12 mL SPE tube with 2 frits (100/box)
SIM-0004-020	Empty 25 mL SPE tube with 2 frits (100/box)
SIM-0006-060	Empty 60 mL SPE tube with 2 frits (100/box)
SIM-0009-150	Empty 150 mL SPE tube with 2 frits (20/box)


SiliaPrep 96-Well Collection Plates

SiliCycle offers SiliaPrep 96-Well Collection Plates, made from polypropylene with extremely low extractable levels.



These collection plates are available with square deep shape in both 1.0 mL and 2.0 mL well volume (*V-shaped bottom*), and with round bottom in 1 mL only. Cap mats are available for all of these collection plates.

SiliaPrep 96-Well Collection Plates	
Product Number	Description
96W-0009	SiliaPrep 96-Well Collection Plate Square Bottom, 2 mL (50/box)
96W-0010	SiliaPrep 96-Well Collection Plate Square Bottom, 1 mL (50/box)
96W-0011	SiliaPrep 96-Well Collection Plate Round Bottom, 1 mL (50/box)


96-Well Collection Plates Square Shape




96-Well Collection Plates Round Shape



1.0 mL Well Volume



2.0 mL Well Volume



1.0 mL Well Volume

SiliaPrep Disposable Reservoir Trays for 96-Well Plates




SiliCycle offers SiliaPrep Disposable Reservoir Trays to collect waste solvents used during activation, loading and washing steps. These disposable trays are made of polycarbonate and are compatible with all manifolds used with well plates.

SiliaPrep Disposable Reservoir Trays	
Product Number	Description
96W-0012	SiliaPrep Disposable Reservoir Trays (25/box)




SiliaPrep 96-Well Plate Cap Mats

SiliCycle offers SiliaPrep 96-Well Plate Cap Mats compatible with most 96-Well Plates available on the market. These cap mats are made of premium-quality silicone, with a PTFE coating for ultra low bleed. Slit and 384-Well Plate cap mats are available under request.



SiliaPrep 96-Well Plate Square
Silicone / PTFE Cap Mats



SiliaPrep 96-Well Plate Round
Silicone / PTFE Cap Mats

SiliaPrep 96-Well Plate Cap Mats			
Well Shape	Quantity	Product Number	Description
Square	5/box	96M-0001S	SiliaPrep 96-Well Plate Square Silicone / PTFE Cap Mats (to be used with 96W-0009 & 96W-0010 collection plates)
	25/box	96M-0001S-25	
	50/box	96M-0001S-50	
	100/box	96M-0001S-100	
Round	5/box	96M-0001R	SiliaPrep 96-Well Plate Round Silicone / PTFE Cap Mats (to be used with 96W-0011 collection plates)
	25/box	96M-0001R-25	
	50/box	96M-0001R-50	
	100/box	96M-0001R-100	

Note: Contact us if you are looking for a cap mat not listed above.

SiliaPrep Phase Separator Cartridges


SiliCycle offers SiliaPrep Phase Separator Cartridges to separate the aqueous phase from heavier chlorinated solvents, under gravity. These ready-to-use cartridges are fitted with a proprietary hydrophobic frit and are a great alternative to liquid-liquid extraction, the most popular technique to do this separation. However, this last method is time consuming, requires the use of a glass funnel (*which needs to be washed between each separate extraction*) and is not suitable for multiple extractions. SiliaPrep Phase Separator Cartridges solve these drawbacks and offer many advantages:

- Ease of use
- Efficient and cost saving
- Compatible with automated systems

SiliaPrep Phase Separator Cartridges	
Product Number	Description
PS-012	SiliaPrep Phase Separator Cartridges, 12 mL (100/box)
PS-060	SiliaPrep Phase Separator Cartridges, 60 mL (50/box)
PS-150	SiliaPrep Phase Separator Cartridges, 150 mL (25/box)

Typical Experimental Procedure

1. Select the appropriate size of SiliaPrep Phase Separator Cartridge to hold your entire sample volume (*both aqueous and chlorinated phases*).
2. Connect the SiliaPrep Phase Separator Cartridge on a vacuum manifold. Ensure the collection vessel volume is large enough to entirely recover the organic layer.
Note: Do not connect the manifold to a vacuum source
3. Transfer the biphasic sample on top of the SiliaPrep Phase Separator Cartridge.
4. After a few seconds (*under gravity*), the water immiscible chlorinated solvent will start to pass through the frit.
5. The proprietary frit used in the SiliaPrep Phase Separator Cartridge allows the aqueous layer to be left on the column for at least 48 hours without passing through the frit.



Important Advice

- **Process under gravity only - Do not apply vacuum or positive pressure**

The SiliaPrep Phase Separator Cartridges are designed to be used under gravity only. The use of a vacuum or positive pressure source can lead to a loss of separation efficiency.

- **Biphasic system required**

The sample needs to contain water and a water immiscible solvent (*with greater density than water, to form the lower layer*). Most common solvents are dichloromethane, chloroform and other chlorinated solvents.

Try to minimize the presence of water miscible solvent (*i.e. methanol, ethanol or acetone*), which can cause problems in obtaining a truly biphasic system. The phase separator may not work effectively if the two phases are merging.

- **More efficient compound partition**

To obtain a more efficient compound partition between aqueous and organic layers, a liquid-liquid extraction can be done prior to use the phase separator column.



SiliaPrep™ Tips

Micro-SPE Cartridges

- Simple, fast analyte retention & elution with minimal loss
- Sorbents directly embedded into inner cartridge wall
- High binding capacity
- No back-pressure


SiliaPrep Tips for Micro Sample Preparation

SiliaPrep Tips Micro-SPE Cartridges are designed for micro-purification and micro-extraction of femtomole (*fmol*) to picomole (*pmol*) quantities of analytes prior to the analysis by chromatographic techniques and / or mass spectrometry.

The constant improvement in these techniques of analysis has allowed scientists to decrease the limit of quantification in several applications. This lower limit has pushed SPE manufacturers to design new SPE cartridges accepting smaller volumes of analyte.

These tips are specially designed to achieve extraction and purification of small molecules, peptides, phosphopeptides and proteins. They are packed with our SiliaBond functionalized silica gels and specialty phases to cover the broadest spectrum of applications requiring small volume of analytes.

The phases are embedded directly in the inner surface of the tip to provide consistent flow rates. Finally, no glue is used during packing procedures in order to prevent any contamination of the analyte.



SiliaPrep Tips Micro-SPE Cartridges Sizes

SiliaPrep Tips Micro-SPE Cartridges are available in 3 different cartridge formats, based on the binding capacity of each embedded sorbent.

SiliaPrep Tips Micro-SPE Cartridges Specifications				
Tip Volume (μL)	Sample Volume (μL)	Binding Capacity (μg)	Sorbent Weight (μg)	Product Number
1 - 10	0.5 - 10	1	30	-T1
10 - 200	2 - 25	2.5	75	-T2
10 - 200	5 - 50	15	400	-T3

SiliaPrep Tips Micro-SPE Cartridges are sold in box of 96.

SiliaPrep Tips General Experimental Procedure

The following lines present the general experimental procedure for the purification and enrichment of small molecules, peptides and proteins using SiliaPrep Tips Micro-SPE Cartridges.

1. Conditioning Step:

Attach the SiliaPrep Tips to a micropipette.
Aspirate / expel the elution solution 5 times and the binding solution 3 times.

2. Loading Step:

Aspirate / expel the sample 20 to 50 times to allow compounds to adsorb onto the sorbent.


3. Washing Step:

Aspirate / expel the binding solution 10 times and discard the expelled solution each time.

4. Elution Step:

Aspirate / expel the elution solution 10 times and collect the expelled solution in a suitable clean tube. Repeat with a fresh portion of elution solution if you want to be sure to collect all of the adsorbed compounds.

Note: repeat 3 - 5 times for carbon black sorbent



SiliaPrep Tips Micro-SPE Application

Micro-Extraction of Dextromethorphan from Plasma	
CARTRIDGE	SiliaPrepX Tips C18 10 µL/30 µg Part Number: SPET-C18-T1
SAMPLE PRETREATMENT	8 µL of plasma sample was mixed with 2 µL of internal standard (Dextromethorphan-d3 at 10 ng/mL in Methanol)
CONDITIONING STEP	8 µL of Methanol (10 aspirate / expel)
EQUILIBRATION STEP	8 µL of water (10 aspirate / expel)
LOADING STEP	Plasma sample (30 aspirate / expel)
WASHING STEP	8 µL of water (10 aspirate / expel) then 8 µL of 25 % Methanol in water (10 aspirate / expel)
ELUTION STEP	8 µL of Acetonitrile (30 aspirate / expel)
FURTHER TREATMENT	Quantification by LDTD/MS/MS (collaboration with Phytronix)
RECOVERY	at 10 ng/mL
	Dextromethorphan 86 %
	Dextromethorphan d-3 80 %



SiliaPrep Tips Sorbent Selection Guide

SiliaPrep Tips Sorbent Selection Guide				
Molecule	Application	Analyte	Sorbent	
Small Molecules	Desalting	All	C18; Carbon Black	
	Protein removal	All	C18; HILIC	
	Metal scavenging	All	Cysteine; DMT; Imidazole; PSA; TAAcOH; TAAcONa; Thiol; Thiourea; Triamine	
	Enrichment	Hydrophobic		C18; HLB; DVB; Carbon Black; HILIC
		Hydrophilic		Silica; Cyano; Carbon Black; HILIC
		Neutral		C18; HLB; DVB; Carbon Black; HILIC; Cyano
		Cationic		SCX; WCX; Polymeric SCX & WCX
Anionic		SAX; NH ₂ ; Polymeric SAX & WAX		
Fluorinated Compounds		Fluoro		
Peptides	Desalting	All	C4; C8; C18; Carbon Black; HILIC	
	SDS removal	All	SDS Removal	
	Enrichment	Glycopeptide		Carbon Black; HILIC; TiO ₂
		Phosphopeptide		TiO ₂ ; ZrO ₂ ; TiO ₂ /ZrO ₂ ; SAX; NH ₂ ; Polymeric SAX & WAX
Other peptide		SAX; NH ₂ ; SCX; WCX; Polymeric SAX, WAX, SCX & WCX		
Proteins	SDS removal	All	SDS Removal	
	Tryptic digestion	All	Trypsin	
	Enrichment	Phosphoprotein		TiO ₂ ; ZrO ₂ ; TiO ₂ /ZrO ₂ ; SAX; NH ₂ ; Polymeric SAX & WAX
		Other protein		SAX; NH ₂ ; SCX; WCX; Polymeric SAX, WAX, SCX & WCX
Oligo-saccharides	Desalting	All	Carbon Black	
	Enrichment	Sulfated glycan		SAX; XSAX
		Sialo-glycan		SAX; XSAX
Other oligosaccharide		Carbon Black; HILIC; TiO ₂		

SiliaPrep XL Tips for Bigger Volumes

For bigger volumes, we also offer SiliaPrep XL Tips Micro-SPE Cartridges, in 3 different formats. Please note these tips are top loading instead of by aspiration.


SiliaPrep XL Tips Micro-SPE Cartridges Specifications				
Tip Volume (μL)	Sample Volume (μL)	Binding Capacity (μg)	Sorbent Weight (mg)	Product Number
1 - 10	1 - 10	400	4	-T1
10 - 200	2 - 25	1,000	10	-T2
100 - 1,000	20 - 1,000	5,000	50	-T3

SiliaPrep XL Tips Micro-SPE Cartridges T1 and T2 are sold in box of 96, T3 in box of 20.

SiliaPrep Tips Sorbent Descriptions & Ordering Information

SiliaPrep Tips Sorbent Descriptions and Ordering Information				
SiliaPrep Tips Sorbent	Description	Product Number		
		10 μ L / 30 μ g	200 μ L / 75 μ g	200 μ L / 400 μ g
C18	Highest hydrophobic character sorbent. Mainly used for small peptides and small molecules purification & enrichment or protein / peptide desalting analysis.	SPET-C18-T1	SPET-C18-T2	SPET-C18-T3
C8	Mid-level hydrophobic sorbent. Mainly used for sample treatment of proteins and peptides requiring a lower hydrophobic capacity, and protein / peptide desalting analysis.	SPET-C8-T1	SPET-C8-T2	SPET-C8-T3
C4	Lowest hydrophobic character sorbent. Mainly used for protein purification & enrichment or protein / peptide desalting analysis.	SPET-C4-T1	SPET-C4-T2	SPET-C4-T3
HLB	Polymeric sorbent with an hydrophilic-lipophilic balance. Mainly used for hydrophobic and neutral molecules enrichment.	SPET-HLB-T1	SPET-HLB-T2	SPET-HLB-T3
DVB	Highly hydrophobic polymeric sorbent. Mainly used for hydrophobic and neutral molecules enrichment.	SPET-DVB-T1	SPET-DVB-T2	SPET-DVB-T3
Carbon Black	Hydrophilic and hydrophobic character. Mainly used for purification of oligosaccharides and other macromolecules containing sugar functions, or protein / peptide desalting.	SPET-CB-T1	SPET-CB-T2	SPET-CB-T3
HILIC	Moderately polar sorbent. Mainly used for proteins removal, peptides desalting, small molecules enrichment and detergent removal (<i>broad spectrum of detergents</i>).	SPET-HI-T1	SPET-HI-T2	SPET-HI-T3
Cyano (CM)	Both polar and hydrophobic character. Mainly used for hydrophilic and neutral molecules enrichment.	SPET-CN-T1	SPET-CN-T2	SPET-CN-T3
Silica	Most polar sorbent. Mainly used for hydrophilic molecules enrichment.	SPET-SI-T1	SPET-SI-T2	SPET-SI-T3
SAX	Strong anion exchanger sorbent. Mainly used for weak acids enrichment.	SPET-SAX-T1	SPET-SAX-T2	SPET-SAX-T3
NH₂ (WAX)	Weak anion exchanger sorbent. Mainly used for strong acids enrichment (<i>phosphopeptides and phosphoproteins</i>).	SPET-NH2-T1	SPET-NH2-T2	SPET-NH2-T3
SCX	Strong cation exchanger sorbent. Mainly used for weak bases enrichment.	SPET-SCX-T1	SPET-SCX-T2	SPET-SCX-T3
WCX	Weak cation exchanger sorbent. Mainly used for strong bases enrichment.	SPET-WCX-T1	SPET-WCX-T2	SPET-WCX-T3
SAX Polymeric	Polymeric sorbent functionalized by a strong anion exchanger. Mainly used for weak acids enrichment.	SPET-XSAX-T1	SPET-XSAX-T2	SPET-XSAX-T3
WAX Polymeric	Polymeric sorbent functionalized by a weak anion exchanger. Mainly used for strong acids enrichment (<i>phosphopeptides and phosphoproteins</i>).	SPET-XWAX-T1	SPET-XWAX-T2	SPET-XWAX-T3
SCX Polymeric	Polymeric sorbent functionalized by a strong cation exchanger. Mainly used for weak bases enrichment.	SPET-XSCX-T1	SPET-XSCX-T2	SPET-XSCX-T3
WCX Polymeric	Polymeric sorbent functionalized by a weak cation exchanger. Mainly used for strong bases enrichment.	SPET-XWCX-T1	SPET-XWCX-T2	SPET-XWCX-T3
TiO₂	High selectivity for multiple phosphorylated peptides. Mainly used for phosphopeptide enrichment and phospholipid removal.	SPET-TI-T1	SPET-TI-T2	SPET-TI-T3
ZrO₂	High selectivity for mono-phosphorylated peptides. Mainly used for phosphopeptide enrichment and phospholipid removal.	SPET-ZR-T1	SPET-ZR-T2	SPET-ZR-T3
TiO₂ / ZrO₂	Excellent alternative for the enrichment of a broad spectrum of phosphopeptides (<i>literature suggests only 30 % overlap in phosphopeptides isolated by TiO₂ versus ZrO₂</i>) and phospholipid removal.	SPET-TIZR-T1	SPET-TIZR-T2	SPET-TIZR-T3
SDS Removal	Used to remove SDS from peptides and proteins.	SPET-SDS-T1	SPET-SDS-T2	SPET-SDS-T3
Trypsin	Used to cleave proteins and peptides at the C-terminal side, with minimal protease contaminants.	SPET-TRYP-T1	SPET-TRYP-T2	SPET-TRYP-T3
Fluoro	Fluorinated sorbent. Mainly used for fluorine containing molecules enrichment.	SPET-FL-T1	SPET-FL-T2	SPET-FL-T3
Metal Scavengers	Mainly used to lower the residual metal concentration of various metal complexes (<i>Pd, Pt, Rh, Ru, Ni, Sn, etc</i>). Choice of 9 metal scavenging sorbents: Cysteine, DMT, Imidazole, Diamine, TAAcOH, TAAcONa, Thiol, Thiourea and Triamine.	SPET-CYS-T1 SPET-DMT-T1 SPET-IMIDAZ-T1 SPET-PSA-T1 SPET-TAAcOH-T1 SPET-TAAcONa-T1 SPET-THIOL-T1 SPET-THIOUREA-T1 SPET-TRINH2-T1	SPET-CYS-T2 SPET-DMT-T2 SPET-IMIDAZ-T2 SPET-PSA-T2 SPET-TAAcOH-T2 SPET-TAAcONa-T2 SPET-THIOL-T2 SPET-THIOUREA-T2 SPET-TRINH2-T2	SPET-CYS-T3 SPET-DMT-T3 SPET-IMIDAZ-T3 SPET-PSA-T3 SPET-TAAcOH-T3 SPET-TAAcONa-T3 SPET-THIOL-T3 SPET-THIOUREA-T3 SPET-TRINH2-T3

Note: Add "XL" after "SPET" for ordering SiliaPrep XL Tips. For example SPETXL-C18-T1.



SiliaQuick™ & SiliaFast™

Sample Preparation & Pesticide Analysis

Two comprehensive solutions available from SiliCycle to simplify your sample prep and analysis: SiliaQuick™ QuEChERS and SiliaFast™ FaPEX®.

If you are frustrated with time and expenses of your sample prep & cleanup procedures, we have simple, economical, highly performant new alternatives to share with you!

SiliaQuick™ QuEChERS

The QuEChERS technique was developed in 2003 by USDA (*United States Department of Agriculture*) scientists to simplify and accelerate the analysis of pesticides in various fruit and vegetable samples. The name QuEChERS is formed by an acronym of the properties that are observed with this technique: **Quick, Easy, Cheap, Effective, Rugged and Safe**.

Initially popularized for the detection and analysis of traces of pesticides in a high throughput environment, scientists have expanded the use of this method to the analysis of a vast array of herbicides, fungicides, antibiotics, drugs and any other compounds present in a myriad of food, beverage, animal and human matrices.

The QuEChERS technique can be summarized as a three-step methodology, starting with a **Liquid Extraction**, followed by a dispersive **Solid-Phase Extraction** clean-up and completed by a **LC or GC Analysis**.

In comparison to traditional sample preparation analysis – a combination of Liquid-Liquid Extraction & Solid-Phase Extraction – the QuEChERS methodology is about 6 times faster, uses 6 - 9 times less solvent, is a safer, greener, much less costly technique, and requires no additional and cumbersome apparatus (*funnels, rotary evaporators, etc.*).



SiliaFast™ FaPEX®

AVAILABLE SOON!







One of the fastest extraction / clean-up approaches for pesticide residue analysis

FaPEX stands for "**F**ast **P**esticide **E**xtraction" and may be considered as "*QuEChERS made even easier*".

This 1-step extraction method preceding LC/MS/MS or GC/MS/MS analysis will ensure you:

- Extraction of thousands of pesticides simultaneously
- Reduction by at least 60 % of labor cost
- Up to 120 X faster than existing methods
- Less operating equipment, less organic solvents and waste than QuEChERS
- Impressive versatility
- High reliability

SiliaFast FaPEX Portfolio		
Cartridge	Name	Matrices
	SiliaFast™ FaPEX-gen	General matrices and all forms
	SiliaFast™ FaPEX-chl	Matrices based on Chlorophyll
	SiliaFast™ FaPEX-cer	Cereal, rice and grains
	SiliaFast™ FaPEX-dry	Tea and dried herbs

Want to learn more?

Contact us: sampleprep@silicycle.com

SiliaQuick QuEChERS Portfolio

Step 1: Liquid Extraction

SiliaQuick QuEChERS Liquid Extraction Step		
Original Method	Buffered Methods	
	AOAC 2007.01 Method	EN 15662 Method
10 g Sample	15 g Sample	10 g Sample
4 g MgSO ₄ ; 1.5 g NaCl	6 g MgSO ₄ ; 1.5 g NaOAc	4 g MgSO ₄ ; 1 g NaCl ; 1 g SCTD ; 0.5 g SCDS
PN: QE-0001-100P (packets only) PN: QE-0001-100K (packets & tubes)	PN: QE-0002-100P (packets only) PN: QE-0002-100K (packets & tubes)	PN: QE-0003-100P (packets only) PN: QE-0003-100K (packets & tubes)

Step 2: dSPE (dispersive Solid Phase Extraction)

SiliaQuick QuEChERS dSPE Step					
Cap Color for 2 mL tubes	Matrix	2 mL tubes for small extract volumes		15 mL tubes for large extract volumes	
		AOAC 2007.01	EN 15662	AOAC 2007.01	EN 15662
Clear	General matrices <ul style="list-style-type: none"> Apples Bananas Broccoli ... 	150 mg MgSO ₄ 50 mg PSA PN: QD-1000-2T	150 mg MgSO ₄ 25 mg PSA PN: QD-1001-2T	1200 mg MgSO ₄ 400 mg PSA PN: QD-2000-15T	900 mg MgSO ₄ 150 mg PSA PN: QD-2001-15T
Pink	Pigmented matrices <ul style="list-style-type: none"> Lettuces Peppers Strawberries ... 	150 mg MgSO ₄ 50 mg PSA 50 mg GCB PN: QD-1002-2T	150 mg MgSO ₄ 25 mg PSA 2.5 mg GCB PN: QD-1003-2T	1200 mg MgSO ₄ 400 mg PSA 400 mg GCB PN: QD-2002-15T	900 mg MgSO ₄ 150 mg PSA 15 mg GCB PN: QD-2003-15T
Green	Highly pigmented matrices <ul style="list-style-type: none"> Urine Avocados Coffee ... 	150 mg MgSO ₄ 50 mg PSA 50 mg GCB 50 mg C18 PN: QD-1004-2T	150 mg MgSO ₄ 25 mg PSA 7.5 mg GCB PN: QD-1005-2T	1200 mg MgSO ₄ 400 mg PSA 400 mg GCB 400 mg C18 PN: QD-2004-15T	900 mg MgSO ₄ 150 mg PSA 45 mg GCB PN: QD-2005-15T
Blue	Fatty and waxed matrices <ul style="list-style-type: none"> Milk Shrimps Blood Liver ... 	150 mg MgSO ₄ 50 mg PSA 50 mg C18 PN: QD-1006-2T	150 mg MgSO ₄ 25 mg PSA 25 mg C18 PN: QD-1007-2T	1200 mg MgSO ₄ 400 mg PSA 400 mg C18 PN: QD-2006-15T	900 mg MgSO ₄ 150 mg PSA 150 mg C18 PN: QD-2007-15T

Bulk Sorbents Available for Your Own Recipe Creation

Bulk Sorbents for QuEChERS			
Product	Product Number	Available Quantities	
SiliaQuick™ Anhydrous Magnesium Sulfate (MgSO ₄)	AUT-0310	<ul style="list-style-type: none"> 5 g 10 g 25 g 50 g 100 g 250 g 500 g 1 kg 5 kg 10 kg 25 kg and more 	
SiliaQuick™ C18	AUT-1313		
SiliaQuick™ Primary Secondary Amine (PSA)	Endcapped		AUT-0312
	Non-endcapped		AUT-1312
SiliaQuick™ Amine	AUT-0412		
SiliaQuick™ Graphitized Carbon Black (GCB)	AUT-0311		



Contact Us

SiliCycle, your worldwide partner

SiliCycle Headquarters & Factory



SiliCycle America Team

america@silicycle.com

SiliCycle European Office & Storage Facility


europa@silicycle.com

SiliCycle Shanghai

contact@silicyclechina.com

SiliCycle Indian Office

india@silicycle.com



Technical Support

At SiliCycle, we are committed to providing the best technical support possible. Our worldwide Technical Support Group is comprised of a team of highly qualified M. Sc., Ph. D. Chemists and Engineers, Technical Support Professionals and Service Coordinators who are prepared to troubleshoot, answer questions and provide solutions for your service and applications needs.


In order to better respond to your technical inquiries, feel free to contact us in three different ways:

E-mail: support@silicycle.com

Phone: • International **+1 418.874.0054**
• USA and Canada **+1 877.745.4292 (Toll-Free)**

Founded in 1995, SiliCycle® is specialized in the development, manufacturing and commercialization of high value silica-based and specialty products for chromatography, analytical and organic chemistry.

E-PAK® : FLOW CARTRIDGES FOR METAL REMOVAL NEW




- Eliminates the use of insoluble particulates in reactors
- High adsorption capacity and flow rate
- Various sizes available for easy scale-up from lab to industrial scale

METAL & ORGANIC SCAVENGING



- Removal of:
- Metals
 - Electrophiles & Nucleophiles
 - Potential Genotoxic Impurities (PGI)
 - Other organic residues

CATALYSIS & SYNTHESIS



- Couplings (*Suzuki, Stille, Heck, ...*)
- Debenzylations & Hydrogenations
- Oxidations
- And Many More Reactions

ACIDS, BASES & REAGENTS



- Acido-basic Reactions
- Amide Couplings
- Reductive Aminations
- Other Reactions

SAMPLE PREPARATION



- SPE & Well Plates
- Micro-SPE Tips
- QuEChERS & FaPEX NEW
- SPE Hardware & Manifold

HIGH PRESSURE CHROMATOGRAPHY



- Bulk Sorbents
- HPLC Columns
- SFC Columns
- Guard Cartridges & Accessories

EQUIPMENTS




- Personal Parallel Synthesizer: MiniBlock® & MiniBlock® XT
- Vacuum Manifold

LOW PRESSURE CHROMATOGRAPHY



- Bulk Silica Gels (*Irregular & Spherical*)
- Bonded Phases
- TLC Plates
- Pre-packed Flash Cartridges

EXTRACTION & PURIFICATION NEW



- Extraction & Purification Services
- Essential Oils & Hydrosols
- Purified Natural Extracts
- Extra-Pure Omega-3

R&D SERVICES



- Scavenging Screening
- Method Development & Optimization
- Impurities Identification
- Custom Column Packing



Distributed By:

Greyhound Chromatography and Allied Chemicals
6 Kelvin Prk, Birkenhead, Merseyside,
CH41 1LT, UK

Orders: By Phone +44 9)0 151 649 4000
By Email: sales@greyhoundchrom.com
Online: www.greyhoundchrom.com