MACHEREY-NAGEL

This combined document includes a brochure for all of the following areas for Macherey Nagel:

- MACHEREY NAGEL Alugram Xtra TLC Aluminium Sheets
- MACHEREY NAGEL SIL HD and Nano-SIL HD
- MACHEREY NAGEL TLC Mikro-Sets/TLC micro-sets
- MACHEREY NAGEL TLC Thin Layer Chromatography
- MACHEREY NAGEL TLC Catalogue

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Chromatography

ALUGRAM[®] Xtra

State of the Art TLC Aluminium Sheets

*NEW
now with
nano silica and
concentrating zone*

200 μ m

outstanding wettability
easy and reliable
excellent separation efficiency

MACHEREY-NAGEL

www.mn-net.com

Since 1911

ALUGRAM® Xtra SIL G · aluminium sheets unmodified standard silica layers on aluminium for TLC

- outstanding wettability for precise colorization results, even with 100 % aqueous eluents
- excellent separation efficiency and reproducibility from lot to lot
- easy and reliable cutting due to an optimized binder system, no flaking of silica

Silica 60, specific surface (BET) ~ 500 m²/g, mean pore size 60 Å, specific pore volume 0.75 mL/g, **particle size 5–17 µm**


Indicator: manganese activated zinc silicate with green fluorescence for short-wave UV (254 nm); special inorganic fluorescent pigment with blue fluorescence for long-wave UV (366 nm)

Binder: highly polymeric product, which is stable in almost all organic solvents and resistant towards aggressive visualization reagents; binder system for ALUGRAM® Xtra is also completely stable in purely aqueous eluents.

Ordering information

Plate size [cm]	2.5 x 7.5	4 x 8	5 x 7.5	5 x 10	5 x 20	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	200	50	20	50	50	20	25		
ALUGRAM® Xtra SIL G · aluminium sheets with standard silica									
SIL G			818230.20	818261	818232		818233	0.20 mm	–
SIL G/UV ₂₅₄	818329	818331	818330.20	818360	818332	818362	818333	0.20 mm	UV ₂₅₄

Separation of nutmeg ingredients


Sample solution: shake 1.0 g freshly powdered drug for 3 min with 4 mL methanol and filter; apply 10 µL
 Layer: ALUGRAM® Xtra SIL G UV₂₅₄
 Eluent: toluene – ethyl acetate (95:5, v/v)
 Migration distance: 15 cm
 Detection: 254 nm: underivatized
 daylight and 366 nm: spray with 5 % ethanolic sulphuric acid, 1 % vanillic acid and heat to 105 °C

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Separation of saffron ingredients

Sample solution: stir 10 mg drug with 50 µL water in a small glass reaction tube. After 3 minutes add 1 mL methanol and store the solution for 20 minutes in the dark. Afterwards filter through a CHROMAFIL® Xtra GF-100/25 filter; apply 10 µL

Layer: ALUGRAM® Xtra SIL G UV₂₅₄
 Eluent: ethyl acetate – 2-propanol – water (65:25:10, v/v/v)
 Migration distance: 10 cm
 Detection: the sheet is dried with a hair dryer and analyzed under daylight, UV 254 nm and 366 nm.

The chromatograms show as main compound Naphtol yellow S. Other colored zones may appear.

MN Appl. No. 403600

NEW

ALUGRAM® Xtra SILGUR · aluminium sheets standard silica layers with concentrating zone for TLC

- 🔸 **concentrating zone:** valuable aid for manual application and time saving
- 🔸 **excellent separation efficiency**
- 🔸 **easy cutting and outstanding wettability**

Silica 60, specific surface (BET) ~ 500 m²/g, mean pore size 60 Å, specific pore volume 0.75 mL/g, **particle size 5–17 µm**

Kieselguhr zone for rapid sample application: because kieselguhr is completely inert towards a large number of compounds, the samples always form a narrow band at the interface of the two adsorbents, irrespective of shape, size or position of the spots in the concentrating zone (see figure). Separation then takes place in the silica layer.

Ordering information


Plate size [cm]	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	20	25		

NEW

ALUGRAM® Xtra · aluminium sheets standard silica layers with concentrating zone

SILGUR	818412	818413	0.20 mm	–
SILGUR UV ₂₅₄	818422	818423	0.20 mm	UV ₂₅₄

Save your precious time!

A valuable aid for manual application especially of large volumes of very dilute samples is the concentrating zone ①, which consists of a chromatography inactive adsorbent (kieselguhr). The substances to be separated are concentrated to a narrow band in the concentrating zone ③. The separation starts at the beginning of the chromatographically active adsorbent silica ④.

Concentrating zone as “rapid application zone” - quantitative evaluation of chromatograms is possible, even if samples applied irregular ②. TLC layers with concentrating zone facilitate the handling and also save time in analysis.

NEW

ALUGRAM® Xtra Nano-SIL G · aluminium sheets nano silica layers for HPTLC

- sharper separations in shorter development time and shorter migration distances
- smaller samples and an increased detection sensitivity
- easy cutting and outstanding wettability

Nano silica 60, specific surface (BET) ~ 500 m²/g, mean pore size 60 Å, specific pore volume 0.75 mL/g, particle size 2–10 µm

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NEW

ALUGRAM® Xtra Nano-SILGUR · aluminium sheets nano silica layers with concentrating zone for HPTLC

- sharper separations in shorter development time and shorter migration distances
- concentrating zone: valuable aid for manual application and time saving
- easy cutting and outstanding wettability

Nano silica 60, specific surface (BET) ~ 500 m²/g, mean pore size 60 Å, specific pore volume 0.75 mL/g, particle size 2–10 µm

Kieselguhr zone for rapid sample application: because kieselguhr is completely inert towards a large number of compounds, the samples always form a narrow band at the interface of the two adsorbents, irrespective of shape, size or position of the spots in the concentrating zone. Separation then takes place in the silica layer.

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NEW


ALUGRAM® Xtra · Nano-SILGUR aluminium sheets nano silica layers with concentrating zone for HPTLC

Nano-SILGUR	818432	0.20 mm	–
Nano-SILGUR UV ₂₅₄	818442	0.20 mm	UV ₂₅₄

MACHEREY-NAGEL

SIL HD and Nano-SIL HD

Chromatography

For precise colorization in TLC

- High luminosity
- Brilliant staining properties
- Excellent separation efficiency

MACHEREY-NAGEL

www.mn-net.com

SIL HD and Nano-SIL HD

Improved TLC glass plates from MACHEREY-NAGEL

- Outstanding dyeability and abrasion resistance due to an optimized binder system
- Good wettability for precise colorization results, even with 100 % aqueous detection reagents
- Excellent separation efficiency due to an optimized particle size distribution
- High suitability for trace analyses resulting from a UV indicator with increased brilliance and a low-noise background of the layer
- Available as glass plates with or without fluorescent indicator

Technical data

SIL HD

- Silica 60, mean pore size 60 Å
- Specific pore volume 0.75 mL/g
- Particle size 5–17 µm (TLC)

Nano-SIL HD

- Silica 60, mean pore size 60 Å
- Specific pore volume 0.75 mL/g
- Particle size 2–10 µm (HPTLC)

Separation of steroids

MN Appl. No. 403810

Plates: SIL HD UV₂₅₄, 20 x 20 cm (REF 809223)
 Competitor M silica gel 60 F₂₅₄, glass backed, 20 x 20 cm

Sample: 1 µL of 0.2 % in acetone


Eluent: Chloroform – methanol (97:3 v/v)

Migration: 6 cm in 15 min

Detection: Treatment with potassium permanganate staining solution and subsequent drying at 120 °C for 10 min

SIL HD UV₂₅₄

Competitor M

Compound	R _f SIL HD UV ₂₅₄	R _f Competitor M
Hydrocortisone	0.07	0.07
Cortexolone	0.21	0.22
Methyltestosterone	0.44	0.41

Improved contrast

The excellent staining property of SIL HD provides a better visibility.

Separation of vitamins

MN Appl. No. 403800

Plates: SIL HD UV₂₅₄, 20 x 20 cm (REF 809223)
 Competitor M silica gel 60 F₂₅₄, glass backed, 20 x 20 cm

Sample: 1 µL of 0.1 % in chloroform – methanol (60:40 v/v)


Eluent: Chloroform – methanol (60:40 v/v)

Migration: SIL HD UV₂₅₄: 7 cm in 25 min
 Competitor M: 7 cm in 20 min

Detection: TLC/HPTLC scanner UV 254 nm

SIL HD UV₂₅₄

Competitor M

Compound	R _f SIL HD UV ₂₅₄	R _f Competitor M
L-ascorbic acid	0.02	0.06
Niacin	0.19	0.16
Nicotinamide	0.34	0.32

Separation of aflatoxins

MN Appl. No. 403740


Plates: SIL HD, 20 x 20 cm (REF 809213)
Competitor M silica gel 60, glass backed, 20 x 20 cm

Sample: 1 μ L
AFB₁, AFG₁: 2 μ g/mL in acetonitrile
AFB₂, AFG₂: 0.5 μ g/mL in acetonitrile

Eluent: Chloroform – acetone (90:10 v/v)

Migration: 7 cm in 20 min

Detection: TLC/HPTLC scanner UV 366 nm

Compound	R _f SIL HD	R _f Competitor M
AFG ₂	0.24	0.25
AFG ₁	0.28	0.30
AFB ₂	0.33	0.34
AFB ₁	0.38	0.41

Separation of aflatoxins (HPTLC)

MN Appl. No. 403750


Plates: Nano-SIL HD, 10 x 20 cm (REF 811213)
Competitor M silica gel 60, glass backed, 10 x 20 cm

Sample: 1 μ L
AFB₁, AFG₁: 2 μ g/mL in acetonitrile
AFB₂, AFG₂: 0.5 μ g/mL in acetonitrile

Eluent: Chloroform – acetone (90:10 v/v)

Migration: 5 cm in 15 min

Detection: TLC/HPTLC scanner UV 366 nm

Compound	R _f Nano-SIL HD	R _f Competitor M
AFG ₂	0.22	0.21
AFG ₁	0.25	0.24
AFB ₂	0.30	0.26
AFB ₁	0.34	0.31

Good to know

For a fast and cost efficient separation of aflatoxins we recommend Nano-SIL HD. Best results can be achieved with UV light of wavelength 366 nm.

Separation of pesticides

MN Appl. No. 403780

Plates: SIL HD UV₂₅₄, 20 x 20 cm (REF 809223)
Competitor M silica gel 60 F₂₅₄, glass backed, 20 x 20 cm


Sample: 1 µL of 0.1 % in acetone

Eluent: *n*-hexane – acetone (80:20 v/v)


Migration: SIL HD UV₂₅₄: 7 cm in 10 min
Competitor M: 7 cm in 8 min

Detection: TLC/HPTLC scanner UV 254 nm

SIL HD UV₂₅₄

Competitor M

Compound	R _f SIL HD UV ₂₅₄	R _f Competitor M
Hexazinone	0.05	0.04
Metoxuron	0.07	0.07
Monuron	0.15	0.13
Azinphos-methyl	0.24	0.21
Prometryn	0.38	0.37
Pyridate	0.50	0.47
Trifluralin	0.60	0.57

Separation of insecticides

MN Appl. No. 403770

Plates: SIL HD UV₂₅₄, 20 x 20 cm (REF 809223)
Competitor M silica gel 60 F₂₅₄, glass backed, 20 x 20 cm

Sample: 1 µL of 0.1 % in dichloromethane

Eluent: *n*-heptane

Migration: 7 cm in 10 min

Detection: TLC/HPTLC scanner UV 254 nm

SIL HD UV₂₅₄

Competitor M

Compound	R _f SIL HD UV ₂₅₄	R _f Competitor M
Endrin	0.02	0.02
DDT	0.24	0.23
Aldrin	0.41	0.41

Separation of sulfonamides

MN Appl. No. 403790

Plates: SIL HD UV₂₅₄, 20 x 20 cm (REF 809223)
Competitor M silica gel 60 F₂₅₄, glass backed, 20 x 20 cm

Sample: 1 µL of 0.1 % in methanol

Eluent: Ethyl acetate – methanol – ammonia sol. (25 %)
(70:15:15 v/v/v)


Migration: 7 cm in 20 min

Detection: TLC/HPTLC scanner UV 254 nm

SIL HD UV₂₅₄

Competitor M

Compound	R _f SIL HD ₂₅₄	R _f Competitor M
Sulfadiazine	0.16	0.16
Sulfamerazine	0.21	0.20
Sulfisoxazole	0.25	0.24
Sulfapyridine	0.40	0.39
Sulfanilamide	0.63	0.59

Separation of flavonoids

MN Appl. No. 403760

Plates: SIL HD UV₂₅₄, 20 x 20 cm (REF 809223)
Competitor M silica gel 60 F₂₅₄, glass backed, 20 x 20 cm

Sample: 1 µL of 0.1 % in methanol


Eluent: Ethyl acetate – butanone – formic acid – water
(50:30:10:10 v/v/v/v)

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
We Meet your Needs

If you have any questions concerning TLC micro-sets or further TLC or chromatography products, please contact us.

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Visit our chromatography website!

- Service and product information about TLC products
- Information about further products for SPE (sample preparation), HPLC and GC
- Free application database

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- E. Hahn-Deinstrop Applied Thin-Layer Chromatography
2nd edition, Wiley-VCH, Weinheim, Germany, 2006
- P. E. Wall Thin-Layer Chromatography:
A Modern Practical Approach
Royal Society of Chemistry, London, UK, 2005
- J. Sherma, B. Fried Handbook of Thin-Layer Chromatography
3rd edition, CRC Press, Boca Raton, USA, 2003
- B. Fried, J. Sherma Thin-Layer Chromatography
4th edition, CRC Press, Boca Raton, USA, 1999
- B. Fried, J. Sherma Practical Thin-Layer Chromatography:
A Multidisciplinary Approach
CRC Press, Boca Raton, USA, 1996
- J. C. Touchstone Practice of Thin Layer Chromatography
Wiley-Interscience, New York, USA, 1992
- N. Grinberg Modern Thin-Layer Chromatography
Marcel Dekker, New York, USA, 1990
- F. Geiss Fundamentals of Thin Layer Chromatography /
Planar Chromatography
Dr. Alfred Hüthig Verlag, Heidelberg, Germany, 1987
- R. E. Kaiser (Editor) Planar Chromatography, Volume 1
Dr. Alfred Hüthig Verlag, Heidelberg, Germany, 1986
- J. G. Kirchner Thin Layer Chromatography
2nd edition, Wiley-Interscience, New York, USA, 1973
- J. C. Touchstone Quantitative Thin Layer Chromatography
Wiley-Interscience, New York, USA, 1973
- E. Stahl Thin Layer Chromatography
1st edition, Springer-Verlag, Berlin, Germany /
Academic Press, New York, USA, 1965
- G. B. Marini-Bettòlo Scientific Reports of the Istituto Superiore di Sanita:
Thin-Layer Chromatography
Elsevier Publishing Company, Amsterdam, Netherlands, 1964

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Contents

German instructions	1
English instructions	
Introduction	40
Principles and basic steps of TLC	41
Sample preparation · Sample application	41
Development of a chromatogram · Visualization of separated substances	42
Evaluation of thin layer chromatograms	43
TLC micro-sets	44
TLC micro-sets for thin layer chromatography	44
TLC micro-set A	45
Description of TLC micro-set A for beginners · Ordering information	45
Experiments with TLC micro-set A	46
Separation of a mixture of fat soluble / lipophilic dyes	46
Separation of a mixture of anthraquinone dyes	47
Separation of a mixture of food dyes	48
Separation of dyes from felt tip pens	49
TLC micro-set F1	50
Description of TLC micro-set F1 · Ordering information	51
Experiments with TLC micro-set F1	52
Separation of amino acids (test mixture and urine)	52
Separation of heavy metal cations	54
TLC micro-set F2	56
Description of TLC micro-set F2 · Ordering information	57
Experiments with TLC micro-set F2	58
Analysis of edible fats	58
Analysis of fats and cholesterol in blood	60
TLC micro-set F3	62
Description of TLC micro-set F3 · Ordering information	63
Experiments with TLC micro-set F3	64
Separation of analgetics (pain relievers)	64
Analysis of drugs using the example of cinchona bark	66
TLC micro-set M	68
Description of TLC micro-set M · Ordering information	69
Further TLC products	70
Ordering information	70
MN layers for TLC	71
Literature	72
Safety instructions	73
Hazard phrases · Precautionary statements	74

Introduction

The chemical analyses of raw materials, the control during production, the final quality and purity control of manufactured products, as well as pharmaceutical, food and environmental analyses have become increasingly important in the last decades.

One important method to solve these problems is thin layer chromatography (TLC). It is an efficient method for analytical as well as for preparative separations and has almost completely replaced historic paper chromatography. In many fields it has been developed into a standard procedure for routine analyses. Besides speed, accuracy and high sensitivity another important advantage is its economy.

From its beginnings in the sixties, MACHEREY-NAGEL has actively supported TLC by developing and manufacturing numerous adsorbents and innovative binder systems. These include mainly silica gel, cellulose, aluminium oxide and polyamide as well as materials modified for special applications.

With these well-proven adsorbents we manufacture machine-coated TLC glass plates as well as polyester or aluminium pre-coated sheets. This eliminates the need for time-consuming and costly hand-coating of TLC plates or sheets, and in addition mechanical production guarantees better uniformity of the layers.

Especially for the education in schools and industry, we have assembled TLC micro-sets, which demonstrate both the speed of the method and the high performance of our POLYGRAM® pre-coated sheets. For this purpose as well as for rapid product control during manufacturing processes in scientific and industrial laboratories pre-coated sheets are particularly suited. Furthermore, we wish to encourage further research and help find solutions to individual problems.

The beginner's set (TLC micro-set A, see page 45–49) features separations with simple developing solvents, samples that are colored thus eliminating the need for visualization and the advantage that all equipment needed is contained in the set.

The advanced sets (TLC micro-sets F1, F2, F3, see page 50–67) are suitable for the more advanced user, because some experience and skill is required. For these sets, some of the samples have to be prepared before separation, and for the identification of substances spray reagents have to be used. Each advanced TLC micro-set is designed for specific groups of separations from physiological and pharmaceutical analysis.

Prerequisite for all advanced experiments is a materials kit (TLC micro-sets M, see page 68–69), which can also be used as basic equipment for individual TLC experiments.

The additional equipment required for the advanced experiments is usually present in any laboratory (e.g., fume hood, safety gloves, hairdryer, drying oven, mortar).

This extension of our programme leads you further into an interesting area of analytical chemistry. All types of TLC micro-sets make use of pre-coated POLYGRAM® sheets, size 4 x 8 cm, which also are used in routine analysis.

Thin layer chromatography (TLC) is a multi-stage distribution process. This process involves solvents or solvent mixtures, sample molecules and a suitable adsorbent. For TLC this adsorbent is applied to a suitable support (a glass plate, polyester or aluminium sheet) in a thin layer – for analytical purposes usually 0.1 to 0.25 mm thick.

For the experiments with our TLC micro-sets the mixture of substances to be separated is applied to a pre-coated polyester sheet and placed into a developing chamber containing the developing solvent. In many cases the separation is completed in a few minutes. Colored substances are visible immediately, colorless substances are visualized by spraying with suitable reagents (formation of colored compounds).

If a UV lamp is available and the adsorbent layer is provided with a fluorescent indicator, substances can often be localized by viewing the chromatogram under UV radiation.

To illustrate how a thin layer chromatographic separation is made, we describe below five basic steps:

1. Sample preparation

A sample for chromatographic separation has to meet several conditions, if good separations are to be achieved. It is not possible to go into exact details here. It should be mentioned, however, that in some cases several steps for concentration or pre-cleaning of a sample are necessary.

The examples in our beginner's set do not require complicated procedures. The dyes or dye mixtures chosen for this set can be applied without further preparation. The last example (dyes from felt tip pens) demonstrates the application of TLC to a product in every-day-life, but in this case, too, the sample can be applied without prior preparation.

The advanced sets require the user to carry out some additional steps such as crushing of a sample, extraction, filtration, and concentration as well as sampling. These individual steps in sample preparation introduce the user to techniques often performed in industrial laboratories. Thorough preparation is an important prerequisite for the success of a thin layer chromatographic separation.

Contrary to the beginner's set, which contains samples of selected compositions, the advanced sets have examples relating to practical tasks. For details, please refer to the individual experiments.


2. Sample application

The aim of a chromatographic separation determines how the sample should be applied to the TLC plate or sheet. For analytical separations with qualitative evaluation samples are spotted with a capillary or micro pipette, for quantitative evaluation the sample is applied as spots or short streaks (mostly 1 cm). A streak covering the whole width of the plate is recommended for preparative separations.

For the separations described in this manual the separation mixtures or reference solutions are applied to a pre-coated sheet as spots by means of glass capillaries. Use each capillary for one sample only to avoid contamination of the following samples.

The capillaries fill themselves quickly when dipped into organic test solutions, for aqueous solutions filling will be much slower; in some cases it may be necessary to use a rubber cap. For emptying the capillaries, place the end of the capillary on the layer vertically and carefully; vertically so that the capillary empties itself, and carefully to avoid damage to the layer. Damaged layers result in unevenly shaped spots. To keep spots as small and compact as possible, it is advisable to apply a solution in several portions with intermediate drying (blow with cold or hot air, cabinet dryer). This is especially important for aqueous solutions.


To facilitate application we recommend a spotting guide designed especially for our pre-coated sheets. The following figures demonstrate the clean and easy application of samples with the above-mentioned spotting guide.

Principles and basic steps of TLC

Depending on whether three or four substances are to be applied, place the pre-coated sheet under the appropriate frame of the spotting guide on a plain anti-glide surface. Thus the starting line and spotting points are fixed at the right distances. Additionally, the spotting guide has a groove on the lower edges to prevent solutions smearing under the guide during application.

For the experiments of the beginner's set only three starting points are needed, whereas some of the advanced experiments require four starting points.

For a later evaluation of thin layer chromatograms, starting points and line are marked with a weak pencil on the TLC layer before sample application.

3. Development of a chromatogram

After application of a sample solution allow the solvent to evaporate completely (about 10 min) or blow dry with cold or hot air. Subsequently place the pre-coated sheet into the developing chamber (glass bottle with screw cap) containing approximately 10 mL of developing solvent. Screw the cap on to the bottle while you hold the glass firmly on the table, so that the surface of the developing solvent is not disturbed.

Please take care that the solvent does not unevenly touch, or splash over the applied spots. The solvent ascends through the layer by capillary action and causes the substances to separate.

When the solvent front has reached the desired height (about 6 to 7 cm) stop the separation by removing the sheet from the developing chamber. Immediately mark the solvent front with a soft pencil.

The atmosphere in the developing chamber should be given some thought. For the examples chosen for our TLC micro-sets saturation of the chamber with developing solvent vapor is not necessary. If, however, you want to achieve reproducible retention, saturation of the chamber with solvent vapor is in many cases advantageous. To obtain saturation, line the developing chamber with an absorptive chromatography paper (MN 260, REF 814030). In this case use 15 mL of solvent instead of the 10 mL mentioned above.

4. Visualization of separated substances

After chromatographic separation all colorless substances have to be made visible. If the adsorbent layer contains a fluorescent indicator, viewing under UV light results in non-specific visualization. Another non-specific identification method is staining in an iodine chamber.

You can obtain more specific color reactions by spraying with suitable spray reagents. Combinations of these methods can also be used. Alternatively, TLC sheets can also be dipped shortly into the reagent solution.

For our beginner's set we have chosen dyes, i.e. colored substances which can be recognized without further treatment.

The substances separated with the advanced sets are colorless as is usually the case in practice. For visualization spray reagents are used (e.g., spray reagent according to Dragendorff-Munier, molybdophosphoric acid, ninhydrin, rubenic acid).

Spray reagent according to Dragendorff-Munier is generally used for detection of nitrogen compounds like alkaloids and tertiary amines. It is composed of basic bismuth nitrate, tartaric acid and potassium iodide and forms a potassium tetraiodo bismutate complex. Nitrogen compounds, which are to be analysed, are protonated by tartaric acid and form orange colored (in some cases yellow or red to brown) ion pairs of the formula $[BiI_4][HNR_3]^+$.

Principles and basic steps of TLC

Molybdato-phosphoric acid, dissolved in ethanol, is applied as spray reagent for the detection of reducing substances (e.g., fats, fatty acids, alcohols). A dark blue coloration of the substance spots, also called »molybdenum blue«, results by heating to 120 °C (heat gun, drying oven) after spraying onto the developed TLC sheet.

For detection of amino acids, amines and amino-sugars an ethanolic ninhydrin solution can be used, which forms reddish spots with the analytes after heating to 110 °C.

Heavy metals can be visualized by the utilization of *rubeanic acid*. Subsequent to the separation rubeanic acid in ethanol is sprayed onto the TLC plates. After a following vaporization with ammonia (vapor of 25% ammonia solution from watchglass or in beaker) the heavy metals will be visible.

More details about spray reagents are described in the experimental part of the micro-sets. But here are some general remarks concerning spraying. Use all spray reagents under a fume hood. For spraying, place the TLC foil on a sheet of filter paper. In most cases it is sufficient to fill the sprayer with 5–10 mL of reagent solution. Spray the sheets from a distance of about 15 cm with the aid of a rubber bulb or compressed air if available. It is always better to spray a layer twice very thinly and evenly (with intermediate drying) than to saturate the layer with excessive spray reagent. In the latter case spots tend to become diffuse. After visualization mark the outlines of color zones with a pencil, because some spots tend to fade after a while.


5. Evaluation of thin layer chromatograms

The evaluation depends on the purpose of a chromatographic separation. For qualitative determinations often localization of a substance is sufficient. This can be achieved by parallel runs with reference substances.

The value often used for evaluation is the R_f value (retention factor) or the hundred fold value hR_f . The R_f value is defined as follows:

$$R_f = \frac{\text{distance starting line} - \text{middle of spot}}{\text{distance starting line} - \text{solvent front}} = \frac{b}{a}$$

i.e. R_f values are between 0 and 1, best between 0.1 and 0.8, corresponding to 10–80 for hR_f .

If reproducible R_f values are to be obtained it is, however, essential that several parameters such as saturation of the development chamber with solvent vapor, constant composition of solvent mixtures, constant temperature etc. are strictly controlled.

Quantitative evaluations are possible by comparison with suitable standards. For this purpose either the areas of spots are measured or substances are measured photometrically directly on the layer. The latter procedure, however, requires complex instrumentation.

Note the safety instructions on labels and in this manual (see page 73–75) before every experimental procedure. The disposal of used chemicals must follow international, national and local environmental protection guidelines.

TLC micro-sets

TLC micro-sets for thin layer chromatography

The TLC micro-sets are assemblies of materials and chemicals for easy separations with thin layer chromatography. They are especially recommended as introduction to TLC for pupils, students and trainees. Each set contains detailed explanations and instructions for the respective separations.

● **Beginner's set (TLC micro-set A)**

This set features separations with simple developing solvents and colored samples.

Thus, a visualization is not necessary.

All equipment needed is contained in the set.

● **Advanced sets (TLC micro-sets F1, F2, F3)**

Advanced sets require some experience and skill from the user.

Some of the samples have to be pretreated before separation.

For identification of substances spray reagents have to be used.

For separations with kits F1 to F3 the materials kit TLC micro-set M is required.

● **Materials kit (TLC micro-set M)**

This kit is prerequisite for the separations with kits F1 to F3.

In addition, it serves as basic equipment for the individual study

of further thin layer chromatographic experiments.

TLC micro-set A

Description of TLC micro-set A for beginners

This kit contains all chemicals, accessories and instructions for the following experiments:

- ✓ Separation of a mixture of fat-soluble/lipophilic dyes (test dye mixture 1): butter yellow, indophenol, sudan blue II, sudan red G
- ✓ Separation of a mixture of anthraquinone dyes (test dye mixture 2): blue 1, blue 3, green, green blue, red, violet 1, violet 2
- ✓ Separation of a mixture of food dyes (test dye mixture 3): brilliant black BN (E151), fast red E, erythrosine (E127), yellow orange S (sunset yellow CFC, E110), naphthol red S, ponceau 4 R (E124), tartrazine (E102)
- ✓ Separation of dyes from felt tip pens

Contents of TLC micro-set A for beginners

- 1 manual (also as pdf download at www.mn-net.com/tlc)
- 3 developing chambers, 1 measuring cylinder 10 mL
- 50 glass capillaries 1 μL , 1 spotting guide
- 50 polyester sheets each of POLYGRAM® SIL G/UV₂₅₄, Alox N/UV₂₅₄ and CEL 300, 4 x 8 cm
- 8 mL each of test dye mixture 1 (4 lipophilic dyes), test dyes sudan blue II (A), and sudan red G (B) (each in toluene)*
- 8 mL each of test dye mixture 2 (7 anthraquinone dyes), test dyes blue 1 (C) and violet 2 (D) (each in toluene – cyclohexane (2:1, v/v))*
- 8 mL each of test dye mixture 3 (7 food dyes), test dyes yellow orange S (E), and brilliant black BN (F) (each in water)
- 100 mL each of toluene, toluene – cyclohexane (2:1, v/v), ethanol, 25 % ammonia solution – 2-propanol (5:3, v/v)* and 2.5 % sodium citrate solution
- 2 felt tip pens

Ordering information

Designation	Pack of	REF
TLC micro-set A for beginners*	1 kit	814000
Replacement parts for TLC micro-set A		
Test dye mixture 1*, solution of 4 lipophilic dyes in toluene (components see above)	8 mL	814001
Test dye mixture 2*, solution of 7 anthraquinone dyes in toluene – cyclohexane (2:1, v/v) (components see above)	8 mL	814002
Test dye mixture 3, aqueous solution of 7 food dyes (components see above)	8 mL	814003
Collection of 4 individual components of test dye mixture 1*	4 x 8 mL	814011
Collection of 7 individual components of test dye mixture 2*	7 x 8 mL	814012
Collection of 7 individual components of test dye mixture 3	7 x 8 mL	814013
Sodium citrate, 2.5 g in 100 mL bottles to fill up with distilled water	2.5 g	814029
TLC polyester sheets POLYGRAM®, 4 x 8 cm: 100 x SIL G/UV ₂₅₄ ; 50 x Alox N/UV ₂₅₄ ; 50 x CEL 300	1 set	814028

* These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see safety instructions (see page 73–75) or MSDS (www.mn-net.com/msds). Our TLC micro-sets are carefully assembled and controlled before leaving our manufactory. In the unlikely case, that by delivery the low boiling solvent in the test dye mixtures should be vaporised, you can refill the vial easily by yourself with the fitting solvent.

- ✓ Test dye mixture 1: 8 mL toluene each
- ✓ Test dye mixture 2: 8 mL toluene – cyclohexane (2:1, v/v) each

Further accessories see page 69.

Experiments with TLC micro-set A

Separation of a mixture of fat soluble / lipophilic dyes

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

Fat soluble / lipophilic dyes are used for coloration of natural and synthetic oils, fats and waxes, benzines, mineral oils as well as for the transparent coloration of plastics. Butter yellow was formerly permitted as food dye.

Chemically these dyes belong to several different classes of substances: azo dyes, azines, indophenols and anthraquinones.

Equipment:

1 x POLYGRAM® Alox N/UV₂₅₄, 1 x developing chamber, 3 x glass capillaries, 1 x measuring cylinder, 1 x spotting guide

Sample solution:

test dye mixture 1
(consisting of butter yellow, sudan blue II, sudan red G, indophenol, already dissolved in toluene)

Reference solutions A and B:

2 individual components (sudan blue II (A) and sudan red G (B), each already dissolved in toluene)

Developing solvent:

toluene (approx. 6 mL)

Separation time:

approx. 7 min

Visualization:


no further visualization required – substances are colored

Procedure:

1. Use one filled (1 µL) glass capillary each and apply as follows using the spotting guide (see page 41: 2. Sample application):
starting point 1: reference solution A (sudan blue II)
starting point 2: test dye mixture 1
starting point 3: reference solution B (sudan red G)
2. Applied dye solutions are air-dried for 5 min.
3. Open the developing chamber and fill in approx. 6 mL toluene using the measuring cylinder.
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).

5. Wait for separation (approx. 7 min), then open the chamber in the fume hood, remove TLC sheet and evaporate the solvent.

Chromatogram:

Legend:

(a) butter yellow, (b) sudan blue II,
(c) sudan red G, (d) indophenol

Remarks:

Separation can also be achieved on a silica gel layer, e.g., POLYGRAM® SIL G/UV₂₅₄. Also use toluene as developing solvent.

Experiments with TLC micro-set A

Separation of a mixture of anthraquinone dyes

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

Anthraquinone dyes are important dyestuffs for the textile and leather industries. Their coloration is very bright and light-fast. They are derived from anthraquinone by substitution of the basic structure or condensation of further ring systems.

Equipment:

1 x POLYGRAM® SIL G/UV₂₅₄, 1 x developing chamber, 3 x glass capillaries, 1 x measuring cylinder, 1 x spotting guide

Sample solution:

test dye mixture 2
(consisting of 7 anthraquinone dyes: blue 1, blue 3, green, green blue, red, violet 1 and violet 2, already dissolved in toluene – cyclohexane (2:1, v/v))

Reference solutions C and D:

2 individual components (blue 1 (C) and violet 2 (D), each already dissolved in toluene – cyclohexane (2:1, v/v))

Developing solvent:

toluene – cyclohexane (2:1, v/v)

Separation time:

approx. 2 x 12 min + 10 min intermediate drying

Visualization:


no further visualization required – substances are colored

Procedure:

1. Use one glass capillary each, and apply as follows $\frac{3}{4}$ of capillary volume using the spotting guide (see page 41: 2. Sample application):
starting point 1: reference solution C (blue 1)
starting point 2: test dye mixture 2
(7 anthraquinone dyes)
starting point 3: reference solution D (violet 2)
2. Applied dye solutions are air-dried for 5 min.
3. Open the developing chamber and fill in approx. 10 mL toluene – cyclohexane using the measuring cylinder.
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 12 min or development height of 65–70 mm).

6. Open the separation chamber, remove TLC sheet and close chamber again.
7. Evaporate the developing solvent from the TLC sheet in the fume hood.
8. Again place the dried TLC sheet into the developing chamber and repeat the separation.

Chromatogram:

Legend:

(a) violet 1, (b) green blue, (c) blue 1, (d) green, (e) red, (f) violet 2, (g) blue 3

Remarks:

In this experiment we used a double development without changing the solvent. This procedure is recommended whenever numerous substances have to be separated, and a single development does not give satisfactory results. For the second development another solvent or solvent system may be chosen. Thorough intermediate drying in a fume hood is, however, essential in all these cases.

Use of the different adsorbent layers POLYGRAM® Alox N/UV₂₅₄ (aluminium oxide) or POLYGRAM® CEL 300 (cellulose), as well as variation of the developing solvent provide interesting information on chromatographic relations. Further experiments can be carried out at the discretion of the user.

Experiments with TLC micro-set A

Separation of a mixture of food dyes

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

A multiplicity of colors is not only found in chromatography, but also in our everyday food. Consumers sometimes even judge the quality of foods by their color. This is why even today colors are added to some foods.

Food additives are subject to government regulations in most countries. Before any dye (whether natural or synthetic) is permitted for food use, it has to be proven that it is non-toxic and non-carcinogenic.

Our test mixture contains seven dyes, from which not all are permitted for use in food. Thus, the importance of the separation and identification of such substances is evident.

Equipment:

1 x POLYGRAM® CEL 300, 1 x developing chamber, 3 x glass capillaries, 1 x measuring cylinder, 1 x spotting guide

Sample solution:

test dye mixture 3
(consisting of brilliant black BN, erythrosine, fast red E, naphthol red S, ponceau 4R, tartrazine, yellow orange S, already dissolved in water)

Reference solution E and F:

2 individual components (yellow orange S (E) and brilliant black BN (F), each already dissolved in water)

Developing solvent:

2.5% sodium citrate solution – 25% ammonia solution – 2-propanol (20:5:3, v/v/v), prepared by mixing 10 mL of 2.5% sodium citrate solution with 4 mL of ready-mixed ammonia – 2-propanol solution

Separation time:

20–25 min (development height of approx. 65 mm)


Visualization:

no further visualization required – substances are colored

Procedure:

1. Use one glass capillary each, and apply as follows $\frac{1}{4}$ of capillary volume using the spotting guide (see page 41: 2. Sample application):
starting point 1: reference solution E (yellow orange S)
starting point 2: test dye mixture 3
starting point 3: reference solution F (brilliant black BN)
2. Applied dye solutions are air-dried for 5 min.
3. Open the developing chamber and fill in the solvent (see above).
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 20–25 min or development height of 65 mm).
6. Open the separation chamber, remove TLC sheet and evaporate the developing solvent from TLC sheet in the fume hood.

Chromatogram:

Legend:

(a) tartrazine, (b) ponceau 4R, (c) yellow orange S, (d) naphthol red S, (e) fast red E, (f) brilliant black BN, (g₁) erythrosine (secondary spots), (g₂) erythrosine (main spot)

Remarks:

The solvent mixture has to be prepared fresh for every separation. When it stands for a while, it will alter its composition and thus its separating properties.

Erythrosine shows three spots on the chromatogram: one stronger spot (g₂) and two weaker secondary spots (g₁). Eventually another weak spot may appear for this compound.

Experiments with TLC micro-set A

Separation of dyes from felt tip pens

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

In the previous examples test mixtures prepared specifically for this application were separated. This experiment – the separation of colors from felt tip pens – shows, that industrial dye composition are more complex mixtures and thus more difficult to separate. Black felt tip pens are especially suitable, because they are often composed of many different colors. A red line does by no means indicate that it contains only one red dye.

The chromatograms can be used for quality control (brilliance of color) as well as for product identification. Apart from the two pens supplied with the set other samples can be separated. It is interesting to compare the chromatographic behaviour of colors from felt tip pens to inks, for example.

Equipment:

1 x POLYGRAM® SIL G/UV₂₅₄, 1 x developing chamber, 3 x glass capillaries, 1 x measuring cylinder, 1 x spotting guide

Samples:

2 commercial felt tip pens

Developing solvent:

ethanol (approx. 10 mL)

Separation time:

approx. 20 min

Visualization:


no further visualization required – substances are colored

Procedure:

1. Apply directly with each of the felt tip pens – weakly impressed dot
starting point 1: black felt tip pen
starting point 3: red felt tip pen
The starting point 2 can be used for e.g., ink (one filled glass capillary) or for an own felt tip pen.
2. Applied dots are air-dried for 3 min.
3. Open the developing chamber and fill in the solvent (see above).
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 20 min).

6. Open the separation chamber, remove TLC sheet and evaporate the developing solvent in the fume hood.

Chromatogram:

Legend:

(a) purple, (b) light red, (c) yellow, (d) green, (e) black, (f) redish color, (g) yellow orange

Remarks:

The results depend on the amount of sample applied. Make sure that the dots are not too big (too much sample) and the layer is not damaged. If too much of the dye mixture has been applied, a double development with the same solvent may improve the separation.

For the black felt tip pen an impressive separation pattern can be obtained, when the color is applied as a streak along the starting line (use a ruler or the spotting guide).

The solvent mixture suggested for this separation will not separate inks. These can be separated, if a mixture of *n*-butanol – glacial acetic acid – water (12:3:5, v/v/v) is used with the same adsorbent layer (development time approx. 25 min). Note how the chromatograms from the felt tip pens change in this solvent mixture.

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Description of TLC micro-set F1

This kit contains all chemicals required for the separation of:

- ✓ amino acids (test mixture, consisting of alanine, arginine, tryptophan and valine)
- ✓ amino acids in urine
- ✓ the heavy metal cations copper(II) and manganese(II)

For the separations with kit F1, the materials kit TLC micro-set M is prerequisite. (A detailed description is provided on page 69.)

Contents of TLC micro-set F1

- 1 manual (also as pdf download at www.mn-net.com/tlc)
- 50 glass capillaries 1 μ L
- 50 polyester sheets each of POLYGRAM® SIL G/UV₂₅₄ and CEL 300, 4 x 8 cm
- 100 mL each of *n*-butanol, ninhydrin spray reagent (in ethanol), acetone, 25 % ammonia solution, rubeanic acid spray reagent (in ethanol)*
- 50 mL each of 50 % acetic acid, 18 % hydrochloric acid*
- 8 mL each of the amino acid test mixture (see above), arginine reference solution (A) and tryptophan reference solution (B) (each in water)
- 8 mL each of the heavy metal cation test mixture (see above), Mn²⁺ reference solution (C) and Cu²⁺ reference solution (D) (each in water)

Ordering information

Designation	Pack of	REF
TLC micro-set F1*	1 kit	814200
Replacement parts for TLC micro-set F1		
Amino acid test mixture (components see above)	8 mL	814201
Collection of 4 individual components of the amino acid test mixture	4 x 8 mL	814202
Ninhydrin spray reagent (in ethanol)*	100 mL	814203
Heavy metal cation test mixture (Cu ²⁺ , Mn ²⁺)	8 mL	814204
Collection of 2 individual components (Cu ²⁺ , Mn ²⁺) of the cation test mixture	2 x 8 mL	814205
Rubeanic acid spray reagent (in ethanol)*	100 mL	814206
TLC polyester sheets POLYGRAM® SIL G/UV ₂₅₄ , 4 x 8 cm	4 x 50	814025
TLC polyester sheets POLYGRAM® CEL 300, 4 x 8 cm	4 x 50	814027
TLC micro-set M (materials kit)	1 kit	814100

* These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see safety instructions (see page 73–75) or MSDS (www.mn-net.com/msds). Further accessories see page 69.

Experiments with TLC micro-set F1


Separation of amino acids (test mixture and urine)

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

Amino acids are essential for the physiology of all living things. They are important for, e.g., cell growth, liver metabolism, blood coagulation, formation of haemoglobin, function of the nervous system and others. Finally – combined to larger units (peptides) – they are the basic units of proteins. Some can be synthesized by animals or humans (non-essential amino acids), some have to be supplied with the food (essential amino acids). Their separation and identification are thus important analytical tasks.

Natural amino acids can occur free (e.g., in blood serum, urine) or they can be obtained by cleavage of proteins, either by action of enzymes or by acid hydrolysis of substances containing proteins (e.g., hen's egg white). The human body secretes about 1.1 g of free amino acids per day.

Since amino acids are colorless, the spots have to be visualized after separation. This is achieved by spraying with the highly sensitive spray reagent ninhydrin. During heating, colored compounds are formed in a complex reaction (deamination, decarboxylation):

Equipment:

1 x POLYGRAM® SIL G/UV₂₅₄, 1 x developing chamber, 3–4 x glass capillaries, 1 x measuring cylinder, 1 x spotting guide, 1 x laboratory sprayer

Sample solution:

test mixture of amino acids
(aqueous solution of alanine, arginine, tryptophan, valine)

urine sample (optional)

Reference solutions A and B:

2 individual components (arginine (A) and tryptophan (B), each already dissolved in water)

Developing solvent:

n-butanol – glacial acetic acid – water (3:1:1, v/v/v)
prepared freshly by mixing 6 mL *n*-butanol with 4 mL 50% acetic acid (mixture alters its composition with time thus it cannot be prepared for stock)

Separation time:

approx. 60 min for a development height of 60 mm

Visualization:


spray reagent: ethanolic solution of ninhydrin;
subsequent heating in a drying oven
(approx. 5 min at 110 °C)

Procedure:

1. Use one glass capillary each, and apply as follows ½ of capillary volume using the spotting guide (see page 41: 2. Sample application):
starting point 1: reference solution A
(solution of arginine)
starting point 2: test mixture of amino acids
starting point 3: reference solution B
(solution of tryptophan)
starting point 4: one complete content of glass capillary – urine sample (optional)
2. Applied samples are dried with a hairdryer.
3. Open the developing chamber and fill in the solvent (see above) using the measuring cylinder.
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 60 min), then open the chamber in the fume hood, remove TLC sheet and mark the solvent front with a pencil.
6. Dry TLC sheet with the hairdryer – until acetic acid can no longer be smelled.
7. Fill the sprayer with 5–10 mL spray reagent (see above) and spray the TLC sheet, lying on a pad of filter paper, so that it is slightly moistened
8. Dry the TLC sheet with the hairdryer and spray it again.
9. Place the TLC sheet for approx. 5 min at 110 °C into the drying oven (orange-red to brown-red spots will appear).
10. Outline the spots with a pencil and determine the R_f values (see page 43: 5. Evaluation of thin layer chromatograms).

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Remarks:

It is recommended that the samples are applied in portions – with intermediate drying – to keep spots small and compact. This is very important because during the long developing time (approx. 60 min) zones may broaden by diffusion thus deteriorating the separation. The drying steps are essential for good results. The temperature of the drying cabinet must not exceed 130 °C.

It is interesting to run the same separation on a cellulose pre-coated TLC sheet (POLYGRAM® CEL 300). The color of the spots changes to violet, but spots are not as compact as on silica gel. The developing time is only 40 min, but valine and tryptophan cannot be separated. Furthermore, the sequence of the acids is changed as compared to silica gel.

With aluminium oxide useful separations cannot be obtained for this system.

Urine samples show a large number of substances, i.e. a broad spectrum of different amino acids. The result depends on the donor, the time of day and the type of food previously eaten. A defined correlation of single amino acids is not possible with this test arrangement. A double development can improve separation in some cases (do not forget thorough intermediate drying).

Experiments with TLC micro-set F1


Separation of heavy metal cations

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

In most cases in industry as well as in research and development, organic substances have to be separated. But thin layer chromatography is suited for the analyses of inorganic compounds, too.

In this example the heavy metals copper and manganese will be separated and identified. Since many heavy metals are very toxic, not only as elements (e.g., dust), but also as soluble salts, it is important to monitor their accumulation in nature as well as in human bodies. Their concentration plays an important role. Many of the heavy metals are essential to living organisms in traces, but higher concentrations result in numerous negative effects on health.

The heavy metal cations are visualized by spraying with rubeanic acid:

Several metal cations form colored inner complexes (neutral chelate complexes) with rubeanic acid. The complex formation depends on the hydrogen ion concentration. This is why in this experiment the TLC sheet has to be treated with ammonia vapor after spraying.

Equipment:

1 x POLYGRAM® CEL 300, 1 x developing chamber, 3 x glass capillaries, 1 x measuring cylinder, 1 x spotting guide, 1 x laboratory sprayer

Sample solution:

mixture of 2 heavy metal salts (cations Cu^{2+} , Mn^{2+}), already dissolved in water

Reference solutions C and D:

2 individual salts (cations Mn^{2+} (C) and Cu^{2+} (D)), each already dissolved in water

Developing solvent:

acetone – 18 % hydrochloric acid (4:1, v/v)
freshly prepared by mixing 10 mL acetone with 2.5 mL hydrochloric acid

Separation time:

up to a height of approx. 7 cm – possible fluctuation of separation time (approx. 10 min).

Visualization:


spray reagent: rubeanic acid (solution in ethanol) and treatment with vapor of 25 % ammonia solution

Procedure:

1. Use one filled glass capillary each and apply as follows using the spotting guide (see page 41: 2. Sample application):
starting point 1: reference solution C (Mn^{2+})
starting point 2: reference solution D (Cu^{2+})
starting point 3: sample solution (Cu^{2+} , Mn^{2+})
Note: To keep spots small apply the solutions in several portions and dry intermediately with hairdryer.
2. After applying dry again carefully.
3. Open the developing chamber and fill in the solvent (see above) using the measuring cylinder.
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation, then open the chamber in the fume hood, remove TLC sheet and mark the solvent front with a pencil.
6. Dry TLC sheet with a hairdryer.
7. Fill the sprayer with 5 mL spray reagent (see above) and spray the TLC sheet, lying on a pad of filter paper, until it is slightly moistened.
8. Dry the TLC sheet with a hairdryer and spray it again.
9. Place some milliliters of 25 % ammonia solution on a watch glass or in a Petri dish and hold TLC sheet (with the layer bottom-side) in the vapor of the solution.
10. Outline the visualized spots with a pencil.

Experiments with TLC micro-set F1

Chromatogram:

Legend:

(a) Cu^{2+} (black-brown), (b) Mn^{2+} (brown)

From the centers of the spots and the distance starting line – solvent front the characteristic R_f values can be calculated.

Remarks:

In this experiment, too, it is important to keep spots small and compact, by applying them in portions. Contrary to the copper spot the spot for manganese will only appear after treatment with ammonia. After about 30 min the color will disappear again. Because of partial separation of the solvent mixture a very thin secondary front can be seen in front of the copper spot. During separation the yellow $(\text{CuCl}_4)^{2-}$ complex zones can be observed.

Please note: It is important to keep the hydrochloric acid bottle closed for reproducible results! Even slight changes in the acid concentration will alter the R_f values.

TLC micro-set F2

Description of TLC micro-set F2

This kit contains all chemicals required for the analysis of:

- ✓ edible fats
- ✓ fats and cholesterol in blood

For the separations with kit F2, the materials kit TLC micro-set M is prerequisite. (A detailed description is provided on page 69.) Because of the directive 93/42/EWG clinical products have to be CE marked. For this reason, we do not add blood lancets and alcohol pads to the TLC micro-sets. Nevertheless, these can be obtained easily at any dispensing chemist's / pharmacy.

Contents of TLC micro-set F2

- 1 manual (also as pdf download at www.mn-net.com/tlc)
- 50 glass capillaries 1 µL
- 50 polyester sheets POLYGRAM® SIL G/UV₂₅₄, 4 x 8 cm
- 5 disposable pipettes 25 µL
- 5 sample vials N 11, 1.5 mL, with PE caps and seals
- 3 sample vials 30 mL (for butter, margarine and edible oil)
- 100 mL each of cyclohexane and molybdatophosphoric acid spray reagent (in ethanol)*
- 2 x 50 mL acetone with calibrated pipette (1 x enclosed)*
- 25 mL butan-2-one (ethyl methyl ketone)*
- 8 mL cholesterol reference solution (in acetone)*

Ordering information


Designation	Pack of	REF
TLC micro-set F2*	1 kit	814300
Replacement parts for TLC micro-set F2		
Cholesterol reference solution (in acetone)*	8 mL	814301
Molybdatophosphoric acid spray reagent (in ethanol)*	100 mL	814302
TLC polyester sheets POLYGRAM® SIL G/UV ₂₅₄ , 4 x 8 cm	4 x 50	814025
TLC micro-set M (materials kit)	1 kit	814100
* These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see safety instructions (see page 73–75) or MSDS (www.mn-net.com/msds). Further accessories see page 69.		

Experiments with TLC micro-set F2

Analysis of edible fats

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

Chemically speaking fats are glycerol esters of higher fatty acids (e.g., palmitic, stearic, oleic acids). Unsaturated fatty acids contain double bonds, and their glycerol esters are especially important for human nutrition. Besides these esters (the so-called triglycerides) natural and synthetic fats and oils contain other components. These include cholesterol, cholesterol esters, free fatty acids and phospholipids. Exposure to light and air result in the formation of hydroxy fatty acids or their triglycerides, the fat turns »rancid«.

In food, these substances are mainly found in butter, margarine and edible oils. Some of these lipids can also be found in blood; from a medical point of view the determination of cholesterol concentrations in blood is of great importance (see next experiment).

This experiment deals with commercial fats, which are separated into their main components. Depending on brands minor differences may occur.

Substances are made visible by spraying with molybdatophosphoric acid. With reducing compounds dark blue to greyish-black spots are formed (molybdenum blue).

Equipment:

1 x POLYGRAM® SIL G/UV₂₅₄, 1 x developing chamber, 4 x glass capillaries, 3 x sample vials, 30 mL, 1 x acetone vial with calibrated pipette, 1 x measuring cylinder, 1 x spotting guide, 1 x laboratory sprayer

Sample preparation:

Weigh 1 g each of butter, margarine and edible oil in the provided sample vials, add 10 mL acetone each and shake the well-closed vials. Some fats (e.g., butter) form only turbid liquids (emulsions). Nevertheless they can be used as samples without further treatment.

Reference solution:

solution of cholesterol in acetone (1 mg/mL); keep sample well closed!

Developing solvent:

cyclohexane – butan-2-one (10:1, v/v)
freshly mix (directly in the developing chamber)
10 mL cyclohexane with 1 mL butan-2-one

Separation time:

approx. 12 min for a height of 60–65 mm

Visualization:

with spray reagent molybdatophosphoric acid (dissolved in ethanol)


Procedure:

in advance: Prepare samples as described under sample preparation

1. Use one filled glass capillary each and apply as follows using the spotting guide (see page 41: 2. Sample application):
starting point 1: butter
starting point 2: edible oil
starting point 3: reference solution (cholesterol in acetone)
starting point 4: margarine
2. Let the applied spots dry (about 3 min).
3. Open the developing chamber and fill in the solvent (see above) using the measuring cylinder.
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 12 min), then open the chamber in the fume hood, remove TLC sheet, mark the solvent front with a pencil and evaporate the solvent in the fume hood (with a hairdryer).
6. Fill the sprayer with 5 mL spray reagent (see above) and spray the TLC sheet, lying on a pad of filter paper, until it is slightly moistured.
7. Dry the TLC sheet with a hairdryer and heat it for another 3–4 min in a drying oven at 120 °C (or with a hot air gun).
8. Outline the dark blue spots (molybdenum blue) on a light yellow background with a pencil.

Experiments with TLC micro-set F2

Chromatogram:

Legend:

(a) cholesterol esters, (b) triglycerides,
(c) cholesterol

Remarks:

The main components of fats are the triglycerides. Immediately above the cholesterol free fatty acids may appear. Hydroxy fatty acids or their triglycerides have lower R_f values, and phospholipids will remain at the starting line. The weak spots will disappear after some time, thus early marking is advisable.


Commercial brands of margarine are often promoted with the label »free of cholesterol«. But margarine can contain cholesterol, if butterfat (max. 2%) was processed.

Experiments with TLC micro-set F2

Analysis of fats and cholesterol in blood

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

Can you see blood? – Blood, too, contains lipids, which can be separated by TLC. In this experiment they will be separated and the cholesterol concentration will be determined semi-quantitatively.

Cholesterol is the oldest known steroid. It occurs free or as esters in human or animal organs and liquids. It is the main constituent of gall stones, from which it was isolated for the first time.

The role and significance of the blood cholesterol concentration to human health is not undisputed. It is certain, that it is an essential material produced by the organism in amounts of about 6–8 g per day. With a low-fat diet an additional 0.04–0.1 g, with fatty foods up to 1.4 g are consumed daily as free cholesterol or its esters. If the blood cholesterol concentration is too high, cholesterol can be deposited in veins causing arteriosclerosis and a resultant risk of heart attack. For this reason the determination of cholesterol in blood is diagnostically important.

This experiment is to identify lipid components in blood and to semi-quantitatively estimate the content of free cholesterol. From the known concentration of the reference solution (1 mg cholesterol/mL acetone) and the chromatogram of the blood sample spot sizes and intensities of the cholesterol spots are compared.

Equipment:

1 x POLYGRAM® SIL G/UV₂₅₄, 1 x developing chamber, 4 x glass capillaries, 1 x sample vials, 1,5 mL, 1 x disposable pipette, 1 x measuring cylinder, 1 x spotting guide, 1 x laboratory sprayer

Sample preparation:

With the calibrating pipette measure 0.5 mL acetone into a 1.5 mL sample vial and close until needed. Then, with an alcohol pad, clean the top of one of the donor's fingers. Press the tip with the thumb and prick with a sterile lancet. Because of the directive 93/42/EWG clinical products have to be CE marked. For this reason, we do not add blood lancets and alcohol pads to the TLC micro-sets.

Nevertheless, these can be obtained easily at any dispensing chemist's/pharmacy. Press further with the thumb to obtain a sufficiently large drop of blood, from which 25 µL is withdrawn with the disposable pipette (use a rubber cap if needed). Immediately release the blood into the sample vial with acetone, close well and shake. Allow the blood sample to deposit and keep well closed until further use. The supernatant clear solution (serum) contains the lipid components and serves as sample.

Reference solution:

Solution of 1 mg cholesterol in 1 mL acetone. Keep sample well closed!

Developing solvent:

cyclohexane – butan-2-one (10:1, v/v)
freshly mix (directly in the developing chamber)
10 mL cyclohexane with 1 mL butan-2-one

Separation time:

approx. 15 min for a height of 60–65 mm

Visualization:

with spray reagent molybdatophosphoric acid (dissolved in ethanol)

Procedure:


in advance: Prepare samples as described under sample preparation

1. Use one filled glass capillary each and apply as follows using the spotting guide (see page 41: 2. Sample application):
starting point 1: clear serum solution
starting point 2: twofold content of the glass capillary – clear serum solution
starting point 3: threefold content of the glass capillary – clear serum solution
starting point 4: reference solution
2. Let the applied spots dry (approx. 3 min).
3. Open the developing chamber and fill in the solvent (see above).
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 15 min), then open the chamber in the fume hood, remove TLC sheet, mark the solvent front with a pencil and evaporate the solvent in the fume hood (with a hairdryer).

Experiments with TLC micro-set F2

6. Fill the sprayer with 5 mL spray reagent (see above) and spray the TLC sheet, lying on a pad of filter paper, until it is slightly moistened.
7. Dry the TLC sheet with a hairdryer and heat it for another 3–4 min in a drying oven at 120 °C (or with a hot air gun).
8. Outline the substance spots with a pencil.

Chromatogram:

Legend:

- (a) cholesterol esters, (b) triglycerides,
(c) cholesterol

Remarks:

Blood contains not only cholesterol, but also triglycerides and cholesterol esters as lipid components (also see previous experiment).

The chromatograms will of course show individual differences. When 1 mL is applied, the concentration of triglycerides may be so small that they cannot be detected, as shown in the above figure.

Additionally the concentration of lipids depends on the age of the »blood donor« as well as on his nutrition, i.e. whether the donor is fasting or has eaten just a few hours before.

After comparing spot sizes and calculating the cholesterol concentration, the result can be further improved by repeating the experiment with different volumes of blood serum. Since the sample is dissolved in the volatile solvent acetone, it is possible to apply 4 µL, 5 µL or more of the serum solution.

Example for the calculation of free cholesterol contents in blood:

Usually the concentration is given as mg cholesterol per 100 mL blood. The dilution factor of 25 µL blood in 0.5 mL (500 µL) acetone equals 500/25, i.e. 20.

If the cholesterol spot of the 3 µL serum sample shows the same intensity and size as the spot of the 1 µL reference sample (with 1 mg/mL), then the sample had 1/3 of the concentration of the reference solution. This gives:

$$\frac{1}{3} \cdot \frac{500}{25} \cdot 100 = 667$$

i.e. 100 mL blood contained 667 mg free cholesterol.

Here is another interesting calculation for the experiment: Which absolute amount of cholesterol can still be accurately determined on the TLC sheet? This gives a good impression of the sensitivity of this method.

TLC micro-set F3

Description of TLC micro-set F3

This kit contains all chemicals required for:

- ✓ separation of analgetics (pain relievers)
- ✓ drug analysis as shown for cinchona bark

For the separations with kit F3, the materials kit TLC micro-set M is prerequisite.
(A detailed description is provided on page 69.)

Contents of TLC micro-set F3

- 1 manual (also as pdf download at www.mn-net.com/tlc)
- 50 glass capillaries 1 µL
- 50 polyester sheets POLYGRAM® SIL G/UV₂₅₄, 4 x 8 cm
- 5 each of Aspirin® and Thomapyrin® tablets
- 20 folded filters MN 615 1/4, 11 cm diameter
- 3 sample vials 8 mL (for Aspirin® sample, Thomapyrin® sample and cinchona bark extract)
- 5 g cinchona bark
- 100 mL each of ethanol, 2-propanol, toluene – diethyl ether (61:39, v/v) *,
spray reagent for caffeine (in aqueous acetone solution)*
and spray reagent according to Dragendorff-Munier (in water)
- 50 mL each of iron(III) chloride solution* and potassium hexacyanoferrate(III) solution (each in water)
- 30 mL ethyl acetate*
- 25 mL each of 12.5 % ammonia solution and diethylamine*
- 8 mL each of caffeine (A), paracetamol (B), quinine (C) reference solutions (each in ethanol)*

Ordering information

Designation	Pack of	REF
TLC micro-set F3*	1 kit	814400
Replacement parts for TLC micro-set F3		
Caffeine reference solution (in ethanol)*	8 mL	814407
Paracetamol reference solution (in ethanol)*	8 mL	814406
Quinine reference solution (in ethanol)*	8 mL	814405
Spray reagent for caffeine (in aqueous acetone solution)*	100 mL	814401
Spray reagent according to Dragendorff-Munier (in water)	100 mL	814402
Iron(III) chloride solution (in water)*	100 mL	814403
Potassium hexacyanoferrate(III) solution (in water)	100 mL	814404
TLC polyester sheets POLYGRAM® SIL G/UV ₂₅₄ , 4 x 8 cm	4 x 50	814025
TLC micro-set M (materials kit)	1 kit	814100

* These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see safety instructions (see page 73–75) or MSDS (www.mn-net.com/msds). Further accessories see page 69.

Experiments with TLC micro-set F3

Separation of analgetics (pain relievers)

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

Analgetics (pain relievers) are frequently used pharmaceuticals. One of the best known active components is acetylsalicylic acid, used for more than 100 years now.

Besides acetylsalicylic acid many other active compounds are used as analgetics. The total amount of analgetics produced is impossible to estimate. Many of these drugs are distributed as combination drugs (e.g., with caffeine).

For this experiment two brands have been chosen from the vast number of analgetics offered on the market: Aspirin® and Thomapyrin®. The detection of the active components is somewhat more complex than usual, because no spray reagent is known which visualizes all compounds at the same time. This is why two parallel chromatograms are prepared.

Equipment:

2 x POLYGRAM® SIL G/UV₂₅₄, 2 x developing chambers, 4 x glass capillaries, 2 x sample vials, 8 mL, 1 x beaker, 2 x folded filters, 1 x measuring cylinder, 1 x spotting guide, 1 x laboratory sprayer

Sample preparation:

- Pulverize one tablet of Aspirin® in a mortar and suspend it in 8 mL ethanol. Filter the suspension through a folded filter MN 615 1/4 into a beaker and wash with 8 mL ethanol. Transfer the clear filtrate into a sample vial as Aspirin® sample. Discard the remaining residue.
- Crush one tablet of Thomapyrin® in a mortar, suspend with 5 mL 2-propanol and filter through a folded filter into a beaker. Wash with 5 mL 2-propanol and transfer the clear filtrate into a sample vial as Thomapyrin® sample.

Reference solutions:

caffeine, dissolved in ethanol (A)
paracetamol, dissolved in ethanol (B)

Developing solvent:

ethyl acetate – ethanol (98:2, v/v),
prepared from 9.8 mL acetyl acetate and 0.2 mL ethanol.

Separation time:

approx. 10 min

Visualization:

Caffeine is visualized with a special spray reagent. For the detection of acetylsalicylic acid and paracetamol a mixture of two iron salt solutions (FeCl_3 and $\text{K}_3[\text{Fe}(\text{CN})_6]$) is used and mixed directly in the reagent sprayer just before use (see 5. Visualization).

Procedure:

in advance: Prepare samples as described under sample preparation

- Use one filled glass capillary each and apply as follows using the spotting guide (see page 41: 2. Sample application): Do this application in parallel on two TLC sheets:
starting point 1: Thomapyrin® sample
starting point 2: paracetamol solution (B)
starting point 3: Aspirin® sample
starting point 4: caffeine solution (A)
Note: To prevent mistakes – mark TLC sheet with I and II.
- Open the developing chamber and fill in the solvent using the measuring cylinder (see above).
- Place one TLC sheet each into each developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chambers (see image page 42).
- Wait for separation (approx. 10 min), then open the chambers in the fume hood, remove TLC sheets, mark the solvent fronts with a pencil and evaporate the solvent in the fume hood (with a hairdryer).

5. Visualization:

TLC sheet I (fume hood!)

Spray very thinly with special spray reagent for caffeine.

Dry shortly (cold air) – spray again and dry again with cold air.

→ mark brown caffeine spots with a pencil

TLC sheet II (drying oven at 110 °C; 5 min)


Fill 3 mL $\text{K}_3[\text{Fe}(\text{CN})_6]$ solution and subsequently 3 mL FeCl_3 solution into sprayer (disregard precipitations).

Moisten the heated TLC sheet with spray reagent in fume hood and dry it with a hairdryer

→ mark substance spots with a pencil
(blue: paracetamol; greyish brown to violet: acetylsalicylic acid)

Experiments with TLC micro-set F3

Chromatograms:

Legend:

(a) acetylsalicylic acid, (b) paracetamol, (c) caffeine

Remarks:

Whereas Aspirin[®] (starting point 3) contains only acetylsalicylic acid as active component, Thomapyrin[®] (starting point 1) consists of acetylsalicylic acid, paracetamol and caffeine (combined preparation). For the visualization of caffeine, do not dry with hot air, because in this case the stains will disappear immediately. In some cases a slight yellow coloration may be seen for paracetamol and acetylsalicylic acid due to adsorption of iodine.

For the visualization on sheet II, treatment in the drying cabinet should not last too long, because otherwise diffuse, blue spots will appear for acetylsalicylic acid. A smaller spot in front of the acetylsalicylic acid indicates salicylic acid (hydrolysis product in older samples). It is important for this test that the sprayer is clean before filling the reagents. If a blue color can be seen in the sprayer (»Berlin blue«), it has to be cleaned thoroughly with diluted sodium hydroxide solution. After excessive spraying, the whole sheet will turn blue after a short time. After a longer period of time the background will turn blue in any case due to the reduction of Fe^{3+} and the formation of »Berlin blue«. At this point it is important to emphasize again that all spray reagents be used under a well-functioning fume hood.

If a UV lamp is available, the whole separation can be seen on a single sheet under shortwave (254 nm) radiation.

Aspirin[®] Trademark of Bayer AG, Germany
Thomapyrin[®] Trademark of Boehringer Ingelheim, Germany

Analysis of drugs using the example of cinchona bark

Please study the general principles of operation (page 41–43) and the safety instructions (page 73–75) for TLC!

As already mentioned, thin layer chromatography has found many applications for the analyses of substance mixtures. For the analysis of drugs, too, it is a rapid and informative method.

In present day language the term »drugs« is often used as meaning drugs of abuse. In the original sense drugs are dried materials of plant origin with therapeutically active ingredients. They are used as pharmaceuticals either in substance or as extracts.

For this experiment the drug cinchona bark (*Cortex chinæ succirubrae*) has been chosen. It contains certain alkaloids which will be extracted and then separated by thin layer chromatography. Alkaloids are a group of nitrogen containing basic natural substances and pharmacologically highly effective. Depending on the content, one differentiates between main and secondary alkaloids. About 35 cinchona alkaloids are known, among them the 4 main alkaloids quinine, quinidine, cinchonine and cinchonidine.

Pharmaceutically the most important cinchona alkaloids are quinine and quinidine. Quinine is a remedy against malaria and acts antipyretic. Quinidine has a relaxing effect and is used as cardiac depressant (antiarrhythmic).

For visualization of separated substances the spray reagent according to Dragendorff-Munier is used as a general test for nitrogen containing compounds like alkaloids.

Equipment:

1 x POLYGRAM® SIL G/UV₂₅₄, 2 x developing chambers, 3 x glass capillaries, 1 x beaker, 1 x folded filter, 1 x sample vial, 8 mL, 1 x measuring cylinder, 1 x spotting guide, 1 x laboratory sprayer

Sample preparation:

Pulverize approx. 1 g cinchona bark in a mortar. Fill the crushed material into a beaker and add about 1 mL 12.5% ammonia solution and 5 mL methanol. Stirring or shaking occasionally, extract for about 10 minutes at 50–60 °C (water bath). Filter the suspension through a folded filter MN 615 1/4, wash with 2 mL methanol and transfer into the provided sample vial. In the clear, reddish brown filtrate precipitations may appear later, but they do not interfere with the separation.

Reference solution C:

quinine dissolved in ethanol (C)

Developing solvent:

toluene – diethylether – diethylamine
(by addition of 1 mL diethylamine to 9 mL of the mixture toluene – diethylether (61:39, v/v))

Separation time:

approx. 2 x 10 min with intermediate drying

Visualization:

Spray reagent according to Dragendorff-Munier, consisting mainly of basic bismuth nitrate, tartaric acid and potassium iodide. Minor precipitates are without significance.

Procedure:

in advance: Prepare samples as described under sample preparation


1. Use one glass capillary each and apply as follows using the spotting guide (see page 41: 2. Sample application):
starting point 1: content of glass capillary – extract
starting point 2: ½ content of glass capillary – quinine solution (C)
starting point 3: twofold content of glass capillary – extract

Note: To keep spots as small as possible, apply in several portions and dry intermediately with a hairdryer.

2. After applying dry again carefully.
3. Open the developing chamber and fill in the solvent using the measuring cylinder (see above).
4. Place the TLC sheet into the developing chamber, touching the wall of the chamber with the lower end of the TLC sheet. Close the developing chamber (see image page 42).
5. Wait for separation (approx. 10 min), then open the chamber in the fume hood, remove TLC sheet and mark the solvent front with a pencil; evaporate the solvent and place the TLC sheet once more into the developing chamber (same solvent) for 10 min.
6. Dry TLC sheet and fill the sprayer with 5 mL spray reagent (see above).
7. Briefly and evenly spray the TLC sheet, lying on a pad of filter paper, dry it and spray again briefly (very slightly).
8. Mark orange colored zones with a pencil.

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
Remarks:

The largest spots are the main alkaloids. Quinidine and cinchonidine are not separated. The smaller spots can be assigned to the numerous secondary alkaloids. Colored (reddish brown) substances, e.g., resins remain at the starting line.

Visualization can be already made after a single development.

If a UV lamp (254 nm) is available, the separation can be evaluated without use of a spray reagent.

TLC micro-set M

Description of TLC micro-set M

This kit is prerequisite for the separations with kits F1 to F3. In addition, it serves as basic equipment for the individual study of further thin layer chromatographic experiments.

Contents of TLC micro-set M (materials kit)

- 2 x 50 glass capillaries 1 μ L
- 2 spotting guides
- 1 rubber cap for capillaries
- 1 measuring cylinder 10 mL
- 1 beaker 25 mL
- 2 developing chambers
- 1 glass laboratory sprayer with rubber bulb
- 1 plastic syringe 1 mL
- 20 sheets filter paper MN 713, 15 x 21 cm
- 50 polyester sheets each of POLYGRAM[®], 4 x 8 cm
SIL G/UV₂₅₄, Alox N/UV₂₅₄ and CEL 300

Ordering information

Designation	Pack of	REF
TLC micro-set M (materials kit)	1 kit	814100
Replacement parts for TLC micro-set M		
Developing chambers for TLC micro-sets	4	814021
Glass laboratory sprayer with rubber bulb	1	814101
Glass capillaries 1 μ L	3 x 50	814022
Rubber caps for capillaries	2	814102
Plastic syringe, 1 mL content with graduation	1	814104
Spotting guides	2	814023
Measuring cylinders, glass, 10 mL content	2	814024
Filter paper MN 713, 15 x 21 cm	100	814103
TLC polyester sheets POLYGRAM [®] , 4 x 8 cm: 100 x SIL G/UV ₂₅₄ ; 50 x Alox N/UV ₂₅₄ ; 50 x CEL 300	1 set	814028

Further TLC products

As manufacturer of TLC plates and adsorbents MACHEREY-NAGEL offers a comprehensive product line for professional TLC. A small selection of important products is listed below. All other TLC products are described on the MN homepage www.mn-net.com.

● TLC polyester sheets POLYGRAM®

The same unbreakable and easy-to-handle polyester sheets as in TLC micro-sets are available in different sizes and coated with different adsorbents.

● TLC aluminium sheets ALUGRAM® Xtra

ALUGRAM® Xtra are state-of-the-art, easily cuttable aluminium sheets with silica layers, which show an outstanding wettability and are characterized by a high separation efficiency.

● TLC glass plates ADAMANT / HPTLC glass plates Nano-ADAMANT

For high demands in separation, as e.g., in pharmaceutical and food analysis, TLC glass plates ADAMANT and HPTLC glass plates Nano-ADAMANT are offered.

● TLC developing chambers and accessories

The advanced and professional user can individually select from the product range of TLC developing chambers and further accessories.

Ordering information

Designation	Pack of	REF
Further TLC products (selection)		
TLC polyester sheets POLYGRAM® SIL G/UV ₂₅₄ , 20 x 20 cm	25	805023
TLC aluminium sheets ALUGRAM® Xtra SIL G/UV ₂₅₄ , 20 x 20 cm	25	818333
TLC glass plates ADAMANT UV ₂₅₄ , 20 x 20 cm	25	821030
TLC glass plates ADAMANT UV ₂₅₄ , 5 x 10 cm	50	821010
HPTLC glass plates Nano-ADAMANT UV ₂₅₄ , 10 x 10 cm	25	821110
Simultaneous developing chamber for up to 5 TLC plates, up to 20 x 20 cm	1	814019
Simultaneous developing chamber for up to 2 TLC plates, up to 10 x 10 cm	1	814018
MN ALUGRAM® scissors, ground blade, black handle	1	818666
Chromatography paper MN 260, 7.5 cm x 17 cm (for chamber saturation)	100	814030

MN layers for TLC

Phase	Layer
Standard silica	
ADAMANT	silica 60, improved binder system, optimized particle size distribution
SIL G	silica 60, standard grade, particle size 5–17 µm
SILGUR	silica 60 with kieselguhr concentrating zone
Unmodified silica for HPTLC	
Nano-ADAMANT	nano silica 60, optimized binder system and particle size distribution
Nano-SIL	nano silica 60, standard grade, particle size 2–10 µm
Modified silica for HPTLC	
RP-18 W/UV ₂₅₄	nano silica with partial octadecyl modification, wettable with water
RP-2/UV ₂₅₄	silanized silica = dimethyl-modified silica 60
Nano-SIL CN	cyano-modified nano silica
Nano-SIL NH ₂	amino-modified nano silica
Nano-SIL DIOL	diol-modified nano silica
Aluminium oxide	
Alox N	aluminium oxide
Cellulose, unmodified and modified	
CEL 300	native fibrous cellulose MN 300
CEL 400	microcrystalline cellulose MN 400 (AVICEL®)
CEL 300 PEI	polyethyleneimine-impregnated cellulose ion exchanger
Layers for special separations	
POLYAMIDE-6	perlon = ε-polycaprolactame
CHIRALPLATE	RP silica with Cu ²⁺ ions and chiral reagent, for enantiomer separation
SIL G-25 HR	high purity silica 60 with gypsum, recommended for aflatoxin analysis
SIL G-25 Tenside	silica G with ammonium sulfate for separation of surfactants
Nano-SIL PAH	nano silica with special impregnation for PAH analysis

List of recommended textbooks for theory and practice of thin layer chromatography

- E. Hahn-Deinstrop Applied Thin-Layer Chromatography
2nd edition, Wiley-VCH, Weinheim, Germany, 2006
- P.E. Wall Thin-Layer Chromatography:
A Modern Practical Approach
Royal Society of Chemistry, London, UK, 2005
- J. Sherma, B. Fried Handbook of Thin-Layer Chromatography
3rd edition, CRC Press, Boca Raton, USA, 2003
- B. Fried, J. Sherma Thin-Layer Chromatography
4th edition, CRC Press, Boca Raton, USA, 1999
- B. Fried, J. Sherma Practical Thin-Layer Chromatography:
A Multidisciplinary Approach
CRC Press, Boca Raton, USA, 1996
- J. C. Touchstone Practice of Thin Layer Chromatography
Wiley-Interscience, New York, USA, 1992
- N. Grinberg Modern Thin-Layer Chromatography
Marcel Dekker, New York, USA, 1990
- F. Geiss Fundamentals of Thin Layer Chromatography /
Planar Chromatography
Dr. Alfred Hüthig Verlag, Heidelberg, Germany, 1987
- R. E. Kaiser (Editor) Planar Chromatography, Volume 1
Dr. Alfred Hüthig Verlag, Heidelberg, Germany, 1986
- J. G. Kirchner Thin Layer Chromatography
2nd edition, Wiley-Interscience, New York, USA, 1973
- J. C. Touchstone Quantitative Thin Layer Chromatography
Wiley-Interscience, New York, USA, 1973
- E. Stahl Thin Layer Chromatography
1st edition, Springer-Verlag, Berlin, Germany /
Academic Press, New York, USA, 1965
- G. B. Marini-Bettòlo Scientific Reports of the Istituto Superiore di Sanita:
Thin-Layer Chromatography
Elsevier Publishing Company, Amsterdam, Netherlands, 1964

Safety instructions

Safety instructions

All kit components with dangerous goods, respective hazard notes and safety instructions for TLC micro-sets from MACHEREY-NAGEL are listed in the following table.

In accordance with CLP regulation (Regulation on Classification, Labelling and Packaging of Substances and Mixtures; EU/1278/2008) inner packings have to be labeled only with the symbol and the product identifier. Less hazardous substances / mixtures with signal word WARNING and flammable substances / mixtures up to 125 mL or 125 g have not to be labeled with H (Hazard) and P (Precautionary) statements. This easing of labeling is NOT valid for sensitizing substances. Components in solutions with a content of <1 % are not specially declared, if they do not endanger in particular.

Corresponding MSDS (Material Safety Data Sheets) can be downloaded under www.mn-net.com/msds.

Dangerous good	GHS symbol + signal word	Hazard statements	Precautionary statements
TLC micro-set A, REF 814000			
ammonia solution 16 % / 2-propanol 35 % / water 49 % (25 % ammonia solution – 2-propanol (5:3, v/v))	 DANGER	226, 314, 335, 400	210, 260, 273, 280, 301+330+331, 303+361+353, 304+340, 305+351+338, 312, 403+233, 403+235
ethanol > 90 %	 DANGER	225	210, 233, 403+235
toluene	 DANGER	225, 304, 315, 336, 361d, 373	202, 210, 233, 301+310, 302+352, 308+313, 331, 332+313, 403+235, 405
toluene 69 % / cyclohexane 31 % (toluene – cyclohexane (2:1, v/v))	 DANGER	225, 304, 315, 336, 361d, 373, 410	202, 210, 233, 273, 301+310, 302+352, 308+313, 331, 332+313, 391, 403+235, 405
TLC micro-set F1, REF 814200			
acetic acid 50 %	 DANGER	314	260, 280, 301+330+331, 303+361+353, 304+340, 305+351+338
acetone	 DANGER	225, 319, 336, EUH066	210, 233, 280, 305+351+338, 337+313, 403+235
ammonia solution 25 %	 DANGER	314, 335, 400	260, 273, 280, 301+330+331, 303+361+353, 304+340, 305+351+338, 312, 403+233
<i>n</i> -butanol	 DANGER	226, 302, 315, 318, 335, 336	210, 233, 261, 280, 302+352, 304+340, 305+351+338, 310, 330, 332+313, 403+235
ethanol > 90 %	 DANGER	225	210, 233, 403+235
hydrochloric acid 18 %	 WARNING	315, 319, 335	261, 280, 302+352, 304+340, 305+351+338, 312, 332+313, 337+313, 403+233

Safety instructions

Dangerous good	GHS symbol + signal word	Hazard statements	Precautionary statements
TLC micro-set F2, REF 814300			
acetone	  DANGER	225, 319, 336, EUH066	210, 233, 280, 305+351+338, 337+313, 403+235
butan-2-one (ethyl methyl ketone)	  DANGER	225, 319, 336, EUH066	210, 233, 280, 305+351+338, 337+313, 370+378, 403+235
cyclohexane	   DANGER	225, 304, 315, 336, 410	210, 233, 273, 301+310, 302+352, 331, 332+313, 370+378, 391, 403+235, 405
ethanol > 90 %	 DANGER	225	210, 233, 403+235
molybdato phosphoric acid < 5 % + ethanol > 90 %	   DANGER	225, 315, 318	210, 233, 280, 302+352, 305+351+338, 332+313, 403+235
TLC micro-set F3, REF 814400			
acetone < 40 % + iodine < 2,5 % + iron(III) chloride < 5 %	  DANGER	225, 319, 336, 412, EUH066	210, 233, 273, 280, 305+351+338, 337+313, 403+235
ammonia solution 12,5 %	  DANGER	314, 335, 400	260, 273, 280, 301+330+331, 303+361+353, 304+340, 305+351+338, 312, 403+233
diethylamine	   DANGER	225, 302, 312, 314, 332, 335	210, 260, 280, 301+312, 301+330+331, 303+361+353, 304+340, 305+351+338, 403+233, 403+235
ethanol > 90 %	 DANGER	225	210, 233, 403+235
ethyl acetate	  DANGER	225, 319, 336, EUH066	210, 233, 280, 305+351+338, 337+313, 370+378, 403+235
iron(III) chloride < 5 %	 WARNING	319	280, 305+351+338, 337+313
2-propanol	  DANGER	225, 319, 336	210, 233, 280, 305+351+338, 337+313, 370+378, 403+235
toluene 66 % + diethyl ether 34 % (toluene – diethyl ether (61:39, v/v))	   DANGER	225, 302, 304, 315, 336, 361d, 373, EUH066	202, 210, 233, 301+310, 308+313, 330, 331, 332+313, 302+352, 403+235, 405

Hazard statements

H225	Highly flammable liquid and vapor.
H226	Flammable liquid and vapor.
H302	Harmful if swallowed.
H304	May be fatal if swallowed and enters airways.
H312	Harmful in contact with skin.

Safety instructions

H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H332	Harmful if inhaled.
H335	May cause respiratory irritation.
H336	May cause drowsiness or dizziness.
H361d	Suspected of damaging the unborn child.
H373	May cause damage to organs through prolonged or repeated exposure.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.
EUH066	Repeated exposure may cause skin dryness or cracking.

Precautionary statements

P202	Do not handle until all safety precautions have been read and understood.
P210	Keep away from heat / sparks / open flames / hot surfaces. – No smoking.
P233	Keep container tightly closed.
P260	Do not breathe vapors.
P261	Avoid breathing vapors.
P273	Avoid release to the environment.
P280	Wear protective gloves / eye protection.
P301+P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor / physician.
P301+P312	IF SWALLOWED: Call a POISON CENTER or doctor / physician if you feel unwell.
P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P303+P361+P353	IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower.
P304+P340	IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P313	IF exposed or concerned: Get medical advice / attention.
P310	Immediately call a POISON CENTER or doctor / physician.
P312	Call a POISON CENTER or doctor / physician if you feel unwell.
P330	Rinse mouth.
P331	Do NOT induce vomiting.
P332+P313	If skin irritation occurs: Get medical advice / attention.
P337+P313	If eye irritation persists: Get medical advice / attention.
P370+P378	In case of fire: Use water and all extinguishing agents for extinction.
P391	Collect spillage.
P403+P233	Store in a well-ventilated place. Keep container tightly closed.
P403+P235	Store in a well-ventilated place. Keep cool.
P405	Store locked up.
P501	Dispose of contents / container to regulated waste treatment.

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Water analysis

Besides TLC micro-sets, MACHEREY-NAGEL offers the *VISOCOLOR*[®] School reagent case. It is especially designed for schools and caters to the needs of both students and teachers. The case contains colorimetric and titrimetric tests for the determination of the most important water parameters.

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
VISOCOLOR[®] School reagent case

- Compact mini lab
- Easy test procedures
- Flexible applications

KATDE/EN200089 TLC micro-sets manual de/en 3/0/5/10.17 PD
Printed in Germany

Chromatography

Distributed By

Greyhound Chromatography and Allied Chemicals
6 Kelvin Park
Birkenhead
Merseyside, CH41 1LT

Tel: 0151 649 4000 Fax: 0151 649 4001

Email: info@greyhoundchrom.com

Web: <https://www.greyhoundchrom.com>

www.mn-net.com

MACHERY-NAGEL

TLC

Thin Layer Chromatography

State-of-the-art
TLC products


Quality
Efficiency
Selectivity

MACHERY-NAGEL

www.mn-net.com

MACHEREY-NAGEL – Thin Layer Chromatography for more than 5 decades

Why TLC?

- ✓ Fast and cost-saving separation technique
- ✓ Multiple sample application possible
- ✓ Developed plate serves as analytical documentation media
- ✓ Time consuming sample preparation steps can be omitted

MN ready-to-use layers for TLC and HPTLC

- ✓ Comprehensive range of plate sizes, surface chemistries and backings
- ✓ Pre-coated plates ready for immediate use
- ✓ Homogeneous, smooth and well adhering layers
- ✓ Available with UV indicator or non-impregnated
- ✓ Consistent high quality from batch-to-batch and from plate-to-plate

Benefits of TLC

TLC does not require complex, costly maintained instrumentation. The investment for performing successful TLC can be hundred times less than for HPLC. Since the separated compounds remain on the plate, they can be used for further experiments. Method development is simplified by TLC. The amount of solvents required for development is much less than with HPLC.


Standard analytical TLC plates and sheets

Thin layer chromatography can be used for both qualitative and quantitative analysis. Standard analytical TLC plates typically have adsorbent layers that are nominally between 0.20–0.25 mm in thickness.

Preparative TLC plates

Preparative TLC is used for purification and isolation of analytes from impurities. Preparative TLC layers (≥ 0.5 mm) are only on glass plates available.

In order to meet your individual application requirements three different types of backings are available.

TLC plates - glass backing


Glass plates are robust, heat proof and chemically resistant to all common mobile phases and visualization reagents.

POLYGRAM® TLC sheets - Polyester (PET) backing

Polyester sheets are easy to handle, lightweight and flexible. Developed **POLYGRAM®** sheets can also be stored for documentation in laboratory notebooks. Scissors cutting possible.

ALUGRAM® TLC sheets - aluminium backing

Aluminium sheets are easy to handle, lightweight and flexible. High performance silica on **ALUGRAM® Xtra** sheets provides outstanding wettability for precise colorization results, even with 100 % aqueous detection reagents. Moreover **ALUGRAM® Xtra** sheets are easy to cut with scissors. No flaking of silica occurs!

Physical properties of backing materials

Material	glass	polyester	aluminium
Thickness (approx.)	1.3 mm	0.2 mm	0.15 mm
Weight, packing and storage requirement	high	low	low
Torsional strength	ideal	low	relatively high
Temperature stability	high	max. 185°C	high
Susceptible to breakage	yes	no	no
Can be cut with scissors	no	yes	yes

Chemical resistance of support material

against solvents	high	high	high
against mineral acids and conc. ammonia	high	high	low

Stability of the binder system of NP plates in water

Suitability for aqueous detection reagents	depends on phase	very suitable	ALUGRAM®: limited suitability ALUGRAM® Xtra: well suited
--	------------------	---------------	---


Standard silica · TLC ready-to-use glass plates

MN offers SIL G and ADAMANT as silica coated glass plates.


Analytical TLC glass plates:

Silica 60, specific surface (BET) ~ 500 m²/g, mean pore size 60 Å, specific pore volume 0.75 mL/g, particle size 5–17 µm, thickness of layer 0.25 mm

2 different selectivities for the separation of nitroanilines, separations under identical conditions


SIL G

ADAMANT

Brilliant UV indicator on ADAMANT plates

ADAMANT

Competitor M

ADAMANT
without separation

Plate size [cm]	Plates per pack	SIL G*			ADAMANT	
		REF with UV 254 nm	REF with UV 254/366 nm	REF without UV indicator	REF with UV 254 nm	REF without UV indicator
2.5 x 7.5	100	809028.100				
5 x 10	200	809027.200			821010.200	
	50	809027		809017	821010	821040
5 x 20	100	809021	809121	809011	821015	
10 x 10	25	809020		809010	821020	821050
10 x 20	50	809022	809122	809012	821025	821070
20 x 20	25	809023	809123	809013	821030	821060

*Also available as preparative plates in 0.50, 1.00 and 2.00 mm thickness.

POLYGRAM® SIL G and SIL N-HR

Polyester sheets for TLC

Silica 60, specific surface (BET) $\sim 500 \text{ m}^2/\text{g}$, mean pore size 60 \AA , specific pore volume 0.75 mL/g , particle size $5\text{--}17 \text{ }\mu\text{m}$; standard grade

The binder system for **POLYGRAM®** sheets is also completely stable in purely aqueous eluents.

POLYGRAM® SIL G


Designation	Thickness of layer	Plate size [cm]	Fluorescent indicator	Plates per pack	REF
SIL G	0.20 mm	2.5 x 7.5	–	200	805902
SIL G	0.20 mm	4 x 8	–	50	805032
SIL G	0.20 mm	5 x 20	–	50	805012
SIL G	0.20 mm	20 x 20	–	25	805013
SIL G	0.20 mm	40 x 20	–	25	805014
SIL G UV ₂₅₄	0.20 mm	2.5 x 7.5	UV ₂₅₄	200	805901
SIL G UV ₂₅₄	0.20 mm	4 x 8	UV ₂₅₄	50	805021
SIL G UV ₂₅₄	0.20 mm	5 x 20	UV ₂₅₄	50	805022
SIL G UV ₂₅₄	0.20 mm	20 x 20	UV ₂₅₄	25	805023
SIL G UV ₂₅₄	0.20 mm	40 x 20	UV ₂₅₄	25	805024
SIL G UV ₂₅₄	0.20 mm	500 x 20	UV ₂₅₄	1 roll	805017

POLYGRAM® SIL N-HR

Different binder system compared to SIL G results in different separation characteristics

Special feature of **POLYGRAM® SIL N-HR**: **higher gypsum content**

SIL N-HR UV ₂₅₄	0.20 mm	5 x 20	UV ₂₅₄	50	804022
SIL N-HR UV ₂₅₄	0.20 mm	20 x 20	UV ₂₅₄	25	804023


ALUGRAM® Xtra SIL G

Standard silica layers on aluminium for TLC

- ✓ Silica 60, specific surface (BET) ~ 500 m²/g, mean pore size 60 Å, specific pore volume 0.75 mL/g, particle size 5–17 µm; standard grade
- ✓ Outstanding wettability for precise colorization results, even with 100 % aqueous detection reagents
- ✓ Excellent separation efficiency and reproducibility from lot to lot
- ✓ Easy and reliable cutting due to an optimized binder system, no flaking of silica
- ✓ Binder: highly polymeric product, which is stable in almost all organic solvents and resistant towards aggressive visualization reagents; also completely stable in purely aqueous eluents

For applications, visit www.mn-net.com/apps

Tailored to individual requirements

ALUGRAM® Xtra SIL G aluminium sheets

Designation	Thickness of layer	Plate size [cm]	Fluorescent indicator	Plates per pack	REF
SIL G	0.20 mm	20 x 20	–	25	818233
SIL G UV ₂₅₄	0.20 mm	4 x 8	UV ₂₅₄	50	818331
SIL G UV ₂₅₄	0.20 mm	20 x 20	UV ₂₅₄	25	818333

Nano TLC plates

Higher efficiency on smaller particles . . .

- ✓ Sharper separations
- ✓ Shorter developing times and migration distances
- ✓ Smaller sample volumes 0.01–0.1 μL
- ✓ Minimal diffusion
- ✓ Increased detection sensitivity

Analytical HPTLC glass plates:

Silica 60, specific surface (BET) $\sim 500 \text{ m}^2/\text{g}$, mean pore size 60 \AA , specific pore volume 0.75 mL/g , thickness of layer 0.20 mm , mean particle size range $2\text{--}10 \mu\text{m}$

Comparison of ADAMANT and Nano-ADAMANT plates

Separation of anthraquinone dyes

Layers: A: ADAMANT, B: Nano-ADAMANT

Sample: $1 \mu\text{L}$ about 0.1 %

Eluent: toluene-cyclohexane (4:3, v/v)

Migration time: A) 30 min, B) 15 min

Peaks:

1. Blue 3
2. Violet 2
3. Red
4. Green
5. Blue 1
6. Greenish Blue
7. Violet 1



Plate size [cm]	Plates per pack	Nano-SIL		Nano-ADAMANT	
		REF with UV 254 nm	REF without UV indicator	REF with UV 254 nm	REF without UV indicator
Glass plates					
5 x 5	100	811021	811011	821100	821130
10 x 10	25	811022	811012	821110	821140
10 x 20	50	811023	811013	821120	821150
ALUGRAM® Xtra					
5 x 20	50	818342	818240		
20 x 20	25	818343	818241		


Preadsorbent zone

After sample application in the kieselguhr layer the spots migrate to the kieselguhr/silica interface forming narrow bands. Separation then takes place in the silica layer.


Concentrating zone - SILGUR

- ✓ Concentrates sample spots on the plate
- ✓ Simplifies sample preparation and application


Channeled SILGUR plates

Channel-Plates with 19 channels help to prevent cross contamination by separating several samples. More samples can be separated on a plate, and spot areas can be more easily determined.

Designation	Thickness of layer	Plate size [cm]	Fluorescent indicator	Plates per pack	REF
Glass plates					
SILGUR-25	0.25 mm	10 x 20	–	50	810012
SILGUR-25	0.25 mm	20 x 20	–	25	810013
SILGUR-25 UV ₂₅₄	0.25 mm	10 x 20	UV ₂₅₄	50	810022
SILGUR-25 UV ₂₅₄	0.25 mm	20 x 20	UV ₂₅₄	25	810023
Glass plates – Channel Plates					
SILGUR-25-C UV ₂₅₄	0.25 mm	20 x 20	UV ₂₅₄	25	810123
ALUGRAM® Xtra aluminium sheets					
SILGUR	0.20 mm	10 x 20	–	20	818412
SILGUR	0.20 mm	20 x 20	–	25	818413
SILGUR UV ₂₅₄	0.20 mm	10 x 20	UV ₂₅₄	20	818422
SILGUR UV ₂₅₄	0.20 mm	20 x 20	UV ₂₅₄	25	818423
Glass plates					
Nano-SILGUR-20	0.20 mm	10 x 10	–	25	811032
Nano-SILGUR-20 UV ₂₅₄	0.20 mm	10 x 10	UV ₂₅₄	25	811042
ALUGRAM® Xtra aluminium sheets					
Nano-SILGUR	0.20 mm	10 x 10	–	25	818432
Nano-SILGUR UV ₂₅₄	0.20 mm	10 x 10	UV ₂₅₄	25	818442

TLC accessory

Description	REF
Simultaneous developing chamber for TLC	
20 x 20 cm, for up to 5 plates	814019
10 x 10 cm, for up to 2 plates	814018
MN ALUGRAM® scissors	818666

818666

814019

814018

HPLC

GC


SPE and Flash

Syringe filters

Vials and caps

... we Meet your Needs

www.mn-net.com

MACHERY-NAGEL


Distributed By


Greyhound Chromatography and Allied Chemicals, 6, Kelvin Park, Birkenhead, Merseyside, CH41 1LT
Tel: +44 (0) 151 649 4000 Email: info@greyhoundchrom.com Web: www.greyhoundchrom.com


Thin layer chromatography


Contents

Basics.....	266
Introductory kits.....	269
Summary of MN ready-to-use layers.....	272
Unmodified TLC silica layers.....	274
Silica layers with concentrating zone.....	278
Unmodified HPTLC silica layers.....	280
Modified silica layers.....	283
Further layers.....	288
Layers for special TLC separations.....	291
Chromatography papers.....	294
Accessories.....	295
Reagents.....	296
Adsorbents.....	297


Glass plates

ALUGRAM® Xtra aluminum sheets
ALUGRAM® aluminum sheets

POLYGRAM® polyester sheets

Thin layer chromatography (TLC) and high performance thin layer chromatography (HPTLC), also called planar chromatography, are, like all chromatographic techniques, based on a multi-stage distribution process involving

- Suitable adsorbents (the stationary phase) coated as a thin layer onto a suitable support (e.g., glass plate, polyester or aluminum sheet; also see page 272)
- Solvents or solvent mixtures (the mobile phase or eluent)
- Sample molecules

The principle of TLC is known for more than 100 years [11]. The real break-through as an analytical method, however, came about 50 years ago as a consequence of the pioneering work of Egon Stahl [12].

Today TLC has gained increasing importance as an analytical separation technique, which is probably due to effects of instrumentation and automation [13]. At the same time the applicability of thin layer chromatography was enhanced by development of new adsorbents and supports.

Today MACHEREY-NAGEL offers a versatile range of ready-to-use layers, which are the result of 50 years of continuous research and development.

Features of modern TLC / HPTLC

The success of thin layer chromatography as a highly efficient microanalytical separation method is based on a large number of advantageous properties:

- High sample throughput in a short time
- Suitable for screening tests
- Pilot procedure for HPLC and Flash chromatography
- After separation the analytical information can be stored for a longer period of time (the TLC ready-to-use layer acts as storage medium for data)
- Separated substances can be subjected to subsequent analytical procedures (e.g., IR, MS) at a later date
- Rapid and cost-efficient optimization of the separation due to easy change of mobile and stationary phase

Principle steps of a TLC separation

Sample preparation

For separation the sample must meet several requirements to obtain good results. Since the TLC plate is a disposable product, sample preparation in general is not as demanding as for other chromatographic methods. However, eventually several steps for sample pretreatment may be necessary. These include sampling, mechanical crushing, extraction steps, filtration and sometimes enrichment of interesting components or clean-up, i.e. removal of undesired impurities.

Our TLC micro-sets introduce some simple methods of sample pretreatment. The dyes or dye mixtures of the beginner's set do not require complicated procedures. The advanced sets require

the user to carry out some additional steps for preparing a sample, thus introducing the user to techniques often performed in industrial laboratories.


Thorough preparation of samples is an important prerequisite for the success of a TLC separation. For our range of products for more demanding sample pretreatment please see the chapter "SPE" from page 10.

Sample application

The most frequent technique is application with a glass capillary as spot or short streak.

Application as streak will yield better results especially for instrumental quantification. For both types of application some manual skill is required to obtain reproducible results. Substance zones which are too large from the beginning will cause poor separation since during chromatography they will become even larger and more diffuse.

A valuable aid for manual application especially of large volumes of very dilute samples is the concentrating zone (e.g., SILGUR-25 UV₂₅₄), which consists of a chromatographically inactive adsorbent (kieselguhr). The substances to be separated are concentrated to a small band in the concentrating zone and the separation starts at the beginning of the chromatographically active adsorbent silica.

Another method for sample concentration is a short pre-elution (few mm) with a solvent, in which all substances have a high R_f value.

If a quantitative evaluation with a TLC scanner is to follow the separation we recommend to use commercially available sample applicators for spotting. These range from simple spotting guides via nanoapplicators to completely automated spotting devices. Application as streak can be performed automatically by spraying of the sample without touching the layer of the TLC plate. Application as band over the whole width of the TLC plate is especially important for preparative TLC. After application allow the solvent of the samples to evaporate completely (about 10 min) or blow with cold or hot air. Development of a chromatogram should never start before the solvent of the applied samples is evaporated completely.


Developing a chromatogram – separation techniques

The most frequently used separation technique is ascending TLC in a trough chamber (standard method, linear development). Usually it is applied as single development. However, multiple development, with or without change of eluent (step technique) can improve separation results. For 2-dimensional development only 1 spot of the sample is applied in one edge of a plate. After chromatography in the first direction the plate is dried, turned by 90° and developed in the 2nd dimension with another eluent. Thus complicated mixtures give 2-dimensional chromatograms taking advantage of the different separating properties of two eluents.

For selection and optimization of the eluent numerous publications are available. A generally applicable standardized optimization method is described by H. Keuker et al. [14].

It is important to pay attention to the atmosphere in the developing chamber. If reproducible migration distances are required, saturation of the chamber atmosphere with eluent vapor is necessary. For this purpose the developing chamber is lined with well absorbing chromatography paper (e.g., MN 260) and charged with a correspondingly larger volume of eluent.

Evaluation of a thin layer chromatogram

Evaluation depends on the purpose of the chromatographic analysis. For qualitative determination often localization of substances is sufficient. This can be easily achieved by parallel runs with reference substances.


A parameter often used for qualitative evaluation is the R_f value (retention factor) or the 100-fold value hR_f . The R_f value is defined as follows:

$$R_f = \frac{\text{distance starting line} - \text{middle of spot}}{\text{distance starting line} - \text{solvent front}} = \frac{b}{a}$$

i.e. the R_f values are between 0 and 1, best between 0.1 and 0.8 (i.e. 10–80 for hR_f). If reproducible R_f values are to be obtained, it is essential that several parameters such as chamber saturation, composition of solvent mixtures, temperature etc. are strictly controlled.

Quantitative evaluation is possible by suitable calibration measurements. For this purpose either the area of a substance spot is measured or a photometric evaluation is performed directly on the layer. The latter procedure, however, requires a higher instrumental expense.

The following paragraphs describe the most frequently used methods for evaluation in TLC.


Qualitative detection

Qualitative evaluation is generally made directly on the TLC plate via characteristic R_f values of substances, i.e. the ratio of distance start – substance zone to distance start – solvent front and specific chemical reactions.

Visualization of separated substances

First of all it is necessary to recognize the position of a substance spot. Only in very few cases the sample is a dye which can be seen with the naked eye. Much more often for unspecific visualization substances can be viewed under UV light, since many substances show a UV absorption. If a fluorescent indicator is added to the layer, all substances absorbing in the respective region of wave length cause a quenching of the fluorescence, i.e. they appear as dark spots on the fluorescent layer. Customary fluorescent indicators are excited at 254 nm or (less frequently) at 366 nm with a mercury lamp. For our program of fluorescent indicators for TLC please see page 296.

Quenching of the fluorescence

Identification of separated substances is possible via the R_f value compared to the pure compound, which is often applied simultaneously on the same plate.

For a number of compounds their native fluorescence can be used for visualization, which is excited by UV light (mostly long-wave UV) (e.g., aflatoxins). This allows not only determination of the R_f value, but often enables a further qualitative assignment.

If these methods do not allow localization or characterization of a substance, post-chromatographic detection methods can be applied, chemical reactions on the TLC plate [15]. Quite un-specific reactions are iodine adsorption and the charring technique (spraying with sulfuric acid and heat treatment).

More reliable results are possible with specific reagents for spraying or dipping, which form colored or fluorescent compounds with the substances to be detected. Depending on the sensitivity of these reactions they are not only used for group or substance specific characterization (in addition to the R_f value) but also for quantification down to trace levels. As example take the ninhydrin reaction. Formation of a (usually red) zone with this detection method yields the information, that a certain group of substances, e.g., α -amino acids, are present. The R_f value allows further assignment to one or several single compounds.

For identification of a substance a combination of different detection methods can be useful. Thus almost all lipids can be converted to products with light green fluorescence by reaction with 2',7'-dichlorofluorescein. Adsorption of iodine vapor enables a differentiation between saturated and unsaturated lipids or lipids containing nitrogen. And finally the R_f value is a third means of identification.

Here are some general remarks concerning spraying: use all spray reagents under a fume hood. The developed, dried TLC plate or sheet is placed on a sheet of filter paper for spraying. Usually it is sufficient to fill the sprayer with about 5–10 mL solution. Spray from a distance of about 15 cm with the aid of a rubber ball or – if available – with pressurized air. It is always better to spray a layer twice very thinly and evenly (with intermediate drying), than to saturate the layer with excessive spray reagent. In the latter case spots tend to become diffuse. After visualization mark outlines of zones with a lead pencil, because some spots tend to fade after a while.

Especially for quantitative evaluation short dipping of the layer in the respective reagent solution is recommended. For this purpose automatic instruments are commercially available, which allow reproducible dipping.

When a substance is localized on the TLC plate (e.g., under UV), but not yet identified, TLC scanners allow recording of UV spectra of individual substance zones directly on the layer, or the zone is removed by scratching or cutting (for sheets), eluted and further analyzed, e.g., by FT-IR, RAMAN, NMR or mass spectroscopy.

Quantitative evaluation

Often TLC is considered to be only a semiquantitative analytical procedure. This is true for visual evaluation of spots, since the eye can only compare but not measure absolute values. If, however, a direct optical evaluation (“in situ” measurement) is performed on the TLC plate with a thin layer scanner, after measurement of calibration functions, exact quantitative results are possible. Commercial scanners offer many features such as evaluation in absorption and fluorescence, unattended programmed scanning of lanes, multi-wave length measurement, background correction, selectable base line for integration, recording of spectra,

evaluation of circular or anti-circular chromatograms with very high ease of operation. In addition to manual operation control by a computer is possible with respective data collection and storage. Usually wavelengths from 200 to 700 nm are available (visible and UV), e.g., all post-chromatographic (and of course all pre-chromatographic) visualization procedures are evaluated with the proper wavelength, which is determined with the instrument. Time requirements for all these possibilities are extremely low. Interlaboratory experiments with standard deviations of 2 % show how excellent results are obtainable [16].

TLC micro-sets introductory kits for science education

Beginner's set

- Features separations with simple developing solvents; samples are colored thus eliminating the need for visualization.
- All equipment needed is contained in the set.

TLC micro-set A for beginners

This kit contains all chemicals and accessories for the following separations:

- Separation of the fat-soluble (lipophilic)
Test dye mixture 1: butter yellow, indophenol, sudan blue II, sudan red G
- Separation of a mixture of anthraquinone dyes
Test dye mixture 2: blue 1, blue 3, green, green blue, red, violet 1, violet 2
- Separation of a mixture of food dyes
Test dye mixture 3: brilliant black BN (E151), fast red E, erythrosine (E127), yellow orange S (sunset yellow CFC, E110), naphthol red S, ponceau 4 R (E124), tartrazine (E102)
- Separation of dyes from felt tip pens

Advanced sets F1, F2 and F3

- Require some experience and skill from the user: some of the samples have to be pretreated before separation, and for identification of substances spray reagents have to be used

Contents of TLC micro-set A for beginners

- 1 manual
- 3 developing chambers
- 50 glass capillaries 1 μ L
- 1 spotting guide
- 2 felt tip pens
- 1 measuring cylinder 10 mL
- 50 polyester sheets 4 x 8 cm each of POLYGRAM®: SIL G/UV₂₅₄, Alox N/UV₂₅₄ and CEL 300
- 8 mL each of test dye mixture 1 (4 lipophilic dyes), test dyes sudan red G, and sudan blue II
- 8 mL each of test dye mixture 2 (7 anthraquinone dyes), test dyes blue 1 and violet 2
- 8 mL each of test dye mixture 3 (7 food dyes), test dyes yellow orange S, and brilliant black BN
- 100 mL each of toluene, toluene – cyclohexane (2:1, v/v), ethanol, 2.5 % sodium citrate solution, 25 % ammonia solution – 2-propanol (5:3, v/v)

Ordering information

Designation	Pack of	REF
TLC micro-set A for beginners*	1 kit	814000
Replacement parts for TLC micro-set A		
Test dye mixture 1*, solution of 4 lipophilic dyes in toluene (components see above)	8 mL	814001
Test dye mixture 2*, solution of 7 anthraquinone dyes in toluene – cyclohexane (2:1, v/v) (components see above)	8 mL	814002
Test dye mixture 3, aqueous solution of 7 food dyes (components see above)	8 mL	814003
Collection of 4 individual components of test dye mixture 1*	4 x 8 mL	814011
Collection of 7 individual components of test dye mixture 2*	7 x 8 mL	814012
Collection of 7 individual components of test dye mixture 3	7 x 8 mL	814013
Sodium citrate, 2.5 g in 100 mL bottle to fill up with distilled water	2.5 g	814029

* These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see SDS.

Information about the advanced sets F1, F2 and F3 can be found on page 270 and page 271.

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This kit contains all chemicals required for the separation of

- Amino acids (test mixture, consisting of alanine, arginine, tryptophan and valine)
- Amino acids in urine
- The heavy metal cations copper(II) and manganese(II)

Contents of TLC micro-set F1

1 manual, 50 glass capillaries 1 μ L
50 polyester sheets 4 x 8 cm each of POLYGRAM®:
SIL G/UV₂₅₄ and CEL 300
100 mL each of *n*-butanol, ninhydrin spray reagent (0.2 % in ethanol), acetone, 25 % ammonia solution, rubeanic acid spray reagent
50 mL each of 50 % acetic acid, 18 % hydrochloric acid
8 mL each of the amino acid test mixture (see left), tryptophan and arginine reference solutions
8 mL each of the heavy metal cation test mixture (see left), Cu²⁺ and Mn²⁺ reference solutions

TLC micro-set F2

This kit contains all chemicals required

- For analysis of edible fats
- For analysis of fats and cholesterol in blood

Contents of TLC micro-set F2

1 manual, 50 glass capillaries 1 μ L
50 polyester sheets 4 x 8 cm POLYGRAM®:
SIL G/UV₂₅₄
5 disposable pipettes 25 μ L
5 sample vials N 11 (1.5 mL) with PE caps and seals
3 sample vials 30 mL (for butter, margarine and edible oil)
100 mL each of cyclohexane and molybdato-phosphoric acid spray reagent
2 x 50 mL acetone with calibrated pipette
25 mL butan-2-one
8 mL cholesterol reference solution

TLC micro-set F3

This kit contains all chemicals required

- For separation of analgetics (pain relievers)
- For drug analysis as shown for cinchona bark

Contents of TLC micro-set F3

1 manual, 50 glass capillaries 1 μ L
50 polyester sheets 4 x 8 cm POLYGRAM®:
SIL G/UV₂₅₄
5 Aspirin® tablets, 5 Thomapyrin® tablets
20 folded filters MN 615 1/4, 11 cm diameter
3 sample vials 8 mL (for Aspirin® sample, Thomapyrin® sample, cinchona bark extract), 5 g cinchona bark
100 mL each of ethanol, 2-propanol, toluene – diethyl ether je 100 mL Ethanol, 2-Propanol, Toluol – Diethylether (61:39, v/v), spray reagent for caffeine and spray reagent according to Dragendorff-Munier
50 mL each of iron(III) chloride solution and potassium hexacyanoferrate(III) solution, 30 mL ethyl acetate
25 mL each of 12.5 % ammonia solution and diethylamine
8 mL each of caffeine, paracetamol, quinine reference solutions

All experiments with TLC micro-sets F1–F3 require the materials kit (see TLC micro-set M on page 271).

Ordering information

Designation	Pack of	REF
TLC micro-set F1*	1 kit	814200
Refill reagents for TLC micro-set F1		
Amino acid test mixtures (components see previous page)	8 mL	814201
Collection of 4 individual components of the amino acid test mixture	4 x 8 mL	814202
Cation test mixture (components see previous page)	8 mL	814204
Collection of 2 individual components of the cation test mixture (Cu ²⁺ , Mn ²⁺)	2 x 8 mL	814205
TLC micro-set F2*	1 kit	814300
Refill reagents for TLC micro-set F2		
Cholesterol reference solution*	8 mL	814301
TLC micro-set F3*	1 kit	814400
Refill reagents for TLC micro-set F3		
Quinine reference solution*	8 mL	814405
Paracetamol reference solution*	8 mL	814406
Caffeine reference solution*	8 mL	814407
Refill packs TLC sheets for all TLC micro-sets		
TLC polyester sheets POLYGRAM® SIL G/UV ₂₅₄ , 4 x 8 cm	4 x 50	814025
TLC polyester sheets POLYGRAM® Alox N/UV ₂₅₄ , 4 x 8 cm	4 x 50	814026
TLC polyester sheets POLYGRAM® CEL 300, 4 x 8 cm	4 x 50	814027
TLC polyester sheets POLYGRAM® 4 x 8 cm: 100 x SIL G/UV ₂₅₄ ; 50 x Alox N/UV ₂₅₄ ; 50 x CEL 300	1 kit	814028

* These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see SDS.

Accessories for TLC micro-sets can be found under TLC accessories on page 295.

Spray reagents can be found on page 296.

TLC micro-set M


This kit is prerequisite for the separations with kits F1 to F3. In addition, it serves as basic equipment for the individual study of further thin layer chromatographic experiments.

Contents of TLC micro-set M (materials kit)

2 x 50 glass capillaries 1 µL, 2 spotting guides
 1 rubber cap for capillaries
 1 measuring cylinder 10 mL
 1 beaker 25 mL
 2 developing chambers
 1 glass laboratory sprayer with rubber bulb
 1 plastic syringe 1 mL
 20 sheets filter paper MN 713 (15 x 21 cm)
 50 polyester sheets 4 x 8 cm each of POLYGRAM®:
 SIL G/UV₂₅₄, Alox N/UV₂₅₄ and CEL 300

Ordering information

Designation	Pack of	REF
TLC micro-set M (materials kit)	1 kit	814100

Advantages of MN plates and sheets for TLC

Continuous high quality

- Guaranteed by stringent production control including standardized lot tests, surface checks for roughness or cracks as well as hardness and adherence checks

Comprehensive range of phases for TLC / HPTLC

- There is no universal TLC plate which meets all possible types of analyses
- Our versatile range of TLC ready-to-use layers covers many different types of applications

Immediately ready for chromatographic separation

- Coatings or impregnations are not necessary

Homogeneous, smooth, well adhering layers

- An important criterion especially for reproducible quantitative evaluation

Electron microscope photograph of a cross section through a glass plate with silica layer (magnification x 500)

Adsorbents for MN plates and sheets for TLC

Classical adsorbents


- For ~ 80 % of all TLC separations silica 60 (mean pore diameter 60 Å = 6 nm) is used
- Other classical adsorbents are aluminum oxide, cellulose, kieselguhr, ion exchangers and polyamide

Special phases

- Modified silica, like C₁₈ (octadecyl-) cyano-, amino-, diol-, RP-2
- Special layers for specific separations, like PAH- or enantiomer separation

Particle size distribution and thickness of layer

- Are chosen to fit the given type of application (e.g., HPTLC, standard or preparative separations)
- Most MN ready-to-use layers are available with or without fluorescent indicator

Electron microscope photograph of a cross section through an aluminum sheet with silica layer (magnification x 500)

Supports for ready-to-use layers for TLC


	Glass plates G	POLYGRAM® P	ALUGRAM® A / ALUGRAM® Xtra Ax
Physical properties of support materials			
Material	glass	polyester	aluminum
Thickness (approx.)	1.3 mm	0.2 mm	0.15 mm
Weight, packaging and storage requirements	high	low	low
Torsional strength	ideal	low	relatively high
Temperature stability	high	max. 185 °C	high
Susceptible to breakage	yes	no	no
Can be cut with scissors	no	yes	yes
Chemical resistance of support materials			
Against solvents	high	high	high
Against mineral acids and conc. ammonia	high	high	low
Stability of the binder system of NP plates in water			
Suitability for aqueous detection reagents	depending on phase	very suitable	ALUGRAM®: limited suitability; ALUGRAM® Xtra: very suitable

Summary of MN ready-to-use layers

Summary				
Phase	Support*		Layer	Page
Standard silica particle size 5–17 µm				
ADAMANT	G		silica 60, improved binder system, optimized particle size distribution	274
SIL G	G	P A Ax	silica 60, standard grade	276
DURASIL	G		silica 60, special binder system	277
SILGUR	G		Ax silica 60 with kieselguhr concentrating zone	279
Unmodified silica for HPTLC particle size 2–10 µm				
Nano-SILGUR	G		Ax nano silica 60 with kieselguhr concentrating zone	279
Nano-ADAMANT	G		nano silica 60, improved binder system, optimized particle size distribution	281
Nano-SIL	G	A Ax	nano silica 60, standard grade	281
Nano-DURASIL	G		nano silica 60, special binder system	282
Modified silica for HPTLC particle size 2–10 µm				
Nano-SIL C18-50/ Nano-SIL C18-100	G		nano silica with partial or complete C ₁₈ modification	283
RP-18 W/UV ₂₅₄	G	A	nano silica with partial octadecyl modification, wettable with water	284
RP-2/UV ₂₅₄	G	A	silanized silica = dimethyl-modified nano silica 60	284
Nano-SIL CN	G	A	cyano-modified nano silica	285
Nano-SIL NH ₂	G	A	amino-modified nano silica	286
Nano-SIL DIOL	G		diol-modified nano silica	287
Aluminum oxide				
Alox-25 / Alox N	G	P A	aluminum oxide	288
Cellulose, unmodified and modified				
CEL 300	G	P A	native fibrous cellulose MN 300	289
CEL 400	G	P	microcrystalline cellulose MN 400 (AVICEL®)	289
CEL 300 PEI		P	polyethyleneimine-impregnated cellulose ion exchanger	290
CEL 300 AC		P	acetylated cellulose MN 300	290
POLYAMID-6				
POLYAMID-6		P	perlon = ε-polycaprolactame	290
Layers for special separations				
CHIRALPLATE	G		RP silica with Cu ²⁺ ions and chiral reagent, for enantiomer separation of amino acids	291
SIL N-HR		P	high purity silica 60, special binder system, higher gypsum content	291
SIL G-25 HR	G		high purity silica 60 with gypsum, recommended for aflatoxin analysis	292
SIL G-25 Tenside	G		silica G with ammonium sulfate for separation of surfactants	292
Nano-SIL PAH	G		nano silica with special impregnation for PAH analysis	292
IONEX-25 SA-Na		P	mixed layer of strongly acidic cation exchanger and silica	293
IONEX-25 SB-AC		P	mixed layer of strongly basic anion exchanger and silica	293
Alox / CEL-AC-Mix	G		mixed layer of aluminum oxide and acetylated cellulose	293
SILCEL-Mix	G		mixed layer of cellulose and silica	293

* G = Glass plates P = POLYGRAM® polyester sheets A = ALUGRAM® aluminum sheets Ax = ALUGRAM® Xtra aluminum sheets

ADAMANT ^G unmodified standard silica layers

★ Key features

- Outstanding hardness and abrasion resistance due to an optimized binder system
- Increased separation efficiency due to an optimized particle size distribution
- High suitability for trace analysis resulting from a UV indicator with increased brilliance and a lownoise background of the layer


🔧 Technical characteristics

- Silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 5–17 µm

Separation of steroids

MN Appl. No. 402930

Layers: ADAMANT UV₂₅₄, SIL G/UV₂₅₄
 Sample: 0.1 % solution in CHCl₃
 Eluent: chloroform – methanol (97:3, v/v)
 Migration distance: ADAMANT 50 mm in 10 min, SIL G 57 mm in 10 min
 Detection: UV


Substance	R _f ADAMANT	R _f SIL G
Cortisone	0.37	0.27
Corticosterone	0.43	0.30
Testosterone	0.50	0.39
Deoxycorticosterone	0.55	0.46
Progesterone	0.73	0.62

Separation of barbiturates

MN Appl. No. 402950

Layer: ADAMANT UV₂₅₄
 Sample volume: 1 µL
 Eluent: chloroform – acetone (95:5, v/v)
 Migration distance: 70 mm in 20 min
 Detection: UV

Substance	R _f
Thiamylal (0.5 %)	0.69
Thiopental (1.0 %)	0.65
Hexobarbital (5.0 %)	0.41
Pentobarbital (1.0 %)	0.26
Phenobarbital (1.0 %)	0.18

Ordering information

Plate size [cm]	2.5 x 7.5	5 x 10	5 x 10	5 x 20	10 x 10	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	100	50	200	100	25	50	25		

Glass plates

ADAMANT		821040	821040.200		821050		821060	0.25 mm	–	
ADAMANT UV ₂₅₄		821005	821010	821010.200	821015	821020	821025	821030	0.25 mm	UV ₂₅₄

ALUGRAM® Xtra SIL G Ax unmodified standard silica layers on aluminum

★ Key features

- Outstanding wettability for precise colorization results, even with 100 % aqueous detection reagents
- Excellent separation efficiency and reproducibility from lot to lot
- Easy and reliable cutting due to an optimized binder system, no flaking of silica

🔧 Technical characteristics


- Silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 5–17 µm
- Binder: highly polymeric product, which is stable in almost all organic solvents and resistant towards aggressive visualization reagents, also completely stable in purely aqueous eluents

Separation of nutmeg ingredients

MN Appl. No. 403590

Layer: ALUGRAM® Xtra SIL G UV₂₅₄
 Sample: shake 1.0 g freshly powdered drug for 3 min with 4 mL methanol and filter; apply 10 µL
 Eluent: toluene – ethyl acetate (95:5, v/v)
 Migration distance: 15 cm
 Detection: 254 nm: underivatized
 daylight and 366 nm: spray with 5 % ethanolic sulfuric acid, 1 % vanillic acid and heat to 105 °C

The chromatograms show the following zones with increasing R_f values: linalool (bluish grey), eugenol (yellowish brown), myristicin (reddish brown), and anethole (pink-violet). Other colored zones may appear.

Ordering information

Plate size [cm]	2.5 x 7.5	4 x 8	5 x 7.5	5 x 10	5 x 20	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	200	50	20	50	50	20	25		

ALUGRAM® Xtra aluminum sheets


SIL G			818230.20	818261	818232		818233	0.20 mm	–
SIL G/UV ₂₅₄	818329	818331	818330.20	818360	818332	818362	818333	0.20 mm	UV ₂₅₄

Further application examples can be found online in our application database at www.mn-net.com/apps

Unmodified TLC silica layers

SIL G G P A unmodified standard silica layers

Technical characteristics

- Silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 5–17 µm
- Thickness of layer for analytical plates 0.25 mm, for preparative plates 0.5 and 1 mm; for 2 mm preparative layers a slightly coarser material is used
- Indicators: manganese activated zinc silicate with green fluorescence for short-wave UV (254 nm); special inorganic fluorescent pigment with blue fluorescence for long-wave UV (366 nm)
- Binders: highly polymeric products, which are stable in almost all organic solvents and resistant towards aggressive visualization reagents; binder system for POLYGRAM® sheets is also completely stable in purely aqueous eluents

Ordering information

Glass plates

Plate size [cm]	2.5 x 7.5	5 x 10	5 x 10	5 x 20	10 x 10	10 x 20	20 x 20	Thickness of layer
Pack of [plates]	100	50	200	100	25	50	25	
SIL G-25		809017	809017.200	809011		809012	809013	0.25 mm
SIL G-25 UV ₂₅₄	809028.100	809027	809027.200	809021	809020	809022	809023	0.25 mm
SIL G-25 UV ₂₅₄₊₃₆₆				809121		809122	809123	0.25 mm

Glass plates

Pack of [plates]	(preparative TLC)		20	
SIL G-50			809051	0.50 mm
SIL G-50 UV ₂₅₄			809053	0.50 mm

Glass plates

Pack of [plates]	(preparative TLC)		15	
SIL G-100			809061	1.00 mm
SIL G-100 UV ₂₅₄			809063	1.00 mm

Glass plates

Pack of [plates]	(preparative TLC)		12	
SIL G-200			809073	2.00 mm
SIL G-200 UV ₂₅₄			809083	2.00 mm

POLYGRAM® polyester sheets

Plate size [cm]	2.5 x 7.5	4 x 8		5 x 20		20 x 20	40 x 20	
Pack of [plates]	200	50		50		25	25	
SIL G	805902	805032		805012		805013	805014	0.20 mm
SIL G/UV ₂₅₄	805901	805021		805022		805023	805024	0.20 mm
SIL G/UV ₂₅₄					roll 500 x 20 cm	805017		0.20 mm

ALUGRAM® aluminum sheets

Plate size [cm]	2.5 x 7.5	4 x 8	5 x 7.5	5 x 10	5 x 20	10 x 20	20 x 20	
Pack of [plates]	200	50	20	50	50	20	25	
SIL G			818030.20	818161	818032	818163	818033	0.20 mm
SIL G/UV ₂₅₄	818129	818131	818130.20	818160	818132	818162	818133	0.20 mm

Further application examples can be found online in our application database at www.mn-net.com/apps

Unmodified TLC silica layers

DURASIL ^G unmodified standard silica layers

Technical characteristics

- Silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 5–17 µm
- Hard, water-resistant and wettable layers due to a special binder system

Ordering information

Plate size [cm]	5 x 10	5 x 10	5 x 20	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	200	100	50	25		

Glass plates

DURASIL-25				812003	812004	0.25 mm	–
DURASIL-25 UV ₂₅₄	812005	812005.200	812006	812007	812008	0.25 mm	UV ₂₅₄


The most TLC layers are available as glass plate, polyester- or aluminum sheet (also see page 272 and 273).

MN TLC pre-coated layers – qualitative and individual tailored

Kieselguhr zone

- For rapid sample application
- Because kieselguhr is completely inert towards a large number of compounds, the samples always form a narrow band at the interface of the two adsorbents, irrespective of shape, size or position of the spots in the concentrating zone. Separation then takes place in the silica layer.

Silica layers with concentrating zone


SILGUR ^G ^{Ax} unmodified standard silica layers with concentrating zone

Technical characteristics

- Silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 5–17 µm
- Kieselguhr zone for rapid sample application (see page 278)
- Channel-plate with 19 channels help to prevent cross contamination by separating several samples
- More samples can be separated on a plate, and spot areas can be more easily determined

Ordering information

Plate size [cm]	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Glass plates				
Pack of [plates]	50	25		
SILGUR-25	810012	810013	0.25 mm	–
SILGUR-25 UV ₂₅₄	810022	810023	0.25 mm	UV ₂₅₄
Channel-Plates				
Pack of [plates]		25		
SILGUR-25-C UV ₂₅₄		810123	0.25 mm	UV ₂₅₄
ALUGRAM® Xtra aluminum sheets				
Pack of [plates]	20	25		
SILGUR	818412	818413	0.20 mm	–
SILGUR UV ₂₅₄	818422	818423	0.20 mm	UV ₂₅₄

Nano-SILGUR ^G ^{Ax} unmodified HPTLC silica layers with concentrating zone

Technical characteristics

- Nano silica 60, pore size 60 Å, specific surface (BET) ~ 500 m²/g, mean specific pore volume 0.75 mL/g, particle size 2–10 µm
- Kieselguhr zone for rapid sample application (see page 278)

Ordering information

Plate size [cm]	10 x 10	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		
Glass plates			
Nano-SILGUR-20	811032	0.20 mm	–
Nano-SILGUR-20 UV ₂₅₄	811042	0.20 mm	UV ₂₅₄
ALUGRAM® Xtra aluminum sheets			
Nano-SILGUR	818432	0.20 mm	–
Nano-SILGUR UV ₂₅₄	818442	0.20 mm	UV ₂₅₄


Sharper separation by nano silica

Nano silica for HPTLC


- Narrow fractionation of the silica particles allows theoretical plate heights, which are one order of magnitude smaller than on standard silica layers.

Advantages

- Shorter migration distances
- Lower amount of samples required
- Increased detection sensitivity with equal selectivity
- Less developing time

Nano silica



Standard silica

Comparison of ADAMANT and Nano-ADAMANT plates for separation of anthraquinone dyes

Layers: A) ADAMANT
B) Nano-ADAMANT
Sample: 1 μ L, about 0.1 %
Eluent: toluene – cyclohexane (4:3, v/v)
Migration time: A) 30 min, B) 15 min

Peaks:
1. Blue 3
2. Violet 2
3. Red
4. Green
5. Blue 1
6. Greenish blue
7. Violet 1

Unmodified HPTLC silica layers

Nano-ADAMANT ^G unmodified HPTLC silica layers

★ Key features

- Outstanding hardness and abrasion resistance due to an optimized binder system
- Increased separation efficiency due to an optimized particle size distribution
- High suitability for trace analyses resulting from a UV indicator with increased brilliance and a lownoise background of the layer

🔧 Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 2–10 µm

Ordering information

Plate size [cm]	10 x 10	10 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25	50		

Glass plates

Nano-ADAMANT	821140	821150	0.20 mm	–
Nano-ADAMANT UV ₂₅₄	821110	821120	0.20 mm	UV ₂₅₄

Nano-SIL ^{G Ax A} unmodified HPTLC silica layers

🔧 Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 2–10 µm
- Indicator: manganese activated zinc silicate with green fluorescence for short-wave UV (254 nm)
- Binder: highly polymeric product, which is stable in almost all organic solvents and resistant towards aggressive visualization reagents

Ordering information

Plate size [cm]	5 x 5	5 x 20	10 x 10	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	100	50	25	50	25		

Glass plates

Nano-SIL-20	811011	811012	811013	0.20 mm	–
Nano-SIL-20 UV ₂₅₄	811021	811022	811023	0.20 mm	UV ₂₅₄

ALUGRAM[®] Xtra aluminum sheets

Nano-SIL G	818240	818241	0.20 mm	–
Nano-SIL G/UV ₂₅₄	818342	818343	0.20 mm	UV ₂₅₄

ALUGRAM[®] aluminum sheets

Nano-SIL G	818141	0.20 mm	–
Nano-SIL G/UV ₂₅₄	818143	0.20 mm	UV ₂₅₄

Unmodified HPTLC silica layers

Nano-DURASIL ^G unmodified HPTLC silica layers

Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 2–10 µm
- Indicator: manganese activated zinc silicate with green fluorescence for short-wave UV (254 nm)
- Hard, water-resistant and wettable layers due to a special binder system
- Different selectivity compared to ADAMANT and SIL-G plates no reversed phase tendency, more polar than Nano-SIL

Ordering information

Plate size [cm]	10 x 10	10 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25	50		

Glass plates

Nano-DURASIL-20	812010	812011	0.20 mm	-
Nano-DURASIL-20 UV ₂₅₄	812013	812014	0.20 mm	UV ₂₅₄

MACHEREY-NAGEL CHROMABOND[®] SPE and Flash products

High-performance products for sample preparation

- Comprehensive range of RP- and normal phases as well as ion exchangers
- Polymer and silica based phases
- Phases for special applications like food or environmental analysis
- SPE polypropylene columns and cartridges, MULTI 96 plates and SPE accessories
- High throughput SPE
- Flash chromatography cartridges

More information from page 9 onwards as well as online at www.mn-net.com/chroma

Nano-SIL C18 ^G octadecyl-modified HPTLC silica layers

Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, pH stability 2–10, particle size 2–10 µm
- Indicator: acid-resistant product with a pale blue fluorescence for short-wave UV (254 nm), UV-absorbing substances appear as dark-blue to black spots on a light-blue background

Modification

- Partial (50 %) or complete (100 %) octadecyl modification, carbon content 7.5 and 14 %, respectively
- Order of polarity: silica > DIOL > NH₂ > CN > RP-2 > C18-50 > RP-18 W > C18-100

Recommended application

- Reversed phase separation mode with eluents from anhydrous solvents to mixtures with high concentrations of water (see table and figure below)
- Alkaloids, amino acids, preservatives, optical brighteners, barbiturates, polycyclic aromatic hydrocarbons (PAH), drugs, peptides, flavonoids, phenols, indole derivatives, steroids

Ordering information


Plate size [cm]	10 x 10	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

Glass plates

Nano-SIL C18-50	50 % silanized	81 1054	0.20 mm	–
Nano-SIL C18-50 UV ₂₅₄	50 % silanized	81 1064	0.20 mm	UV ₂₅₄
Nano-SIL C18-100	100 % silanized	81 1052	0.20 mm	–
Nano-SIL C18-100 UV ₂₅₄	100 % silanized	81 1062	0.20 mm	UV ₂₅₄

Eluent	v/v	Migration distances [mm/15 min]		
		C18-50	C18-100	RP-18 W
Methanol – H ₂ O	2:1	57	45	44
	1:1	52	21	40
	1:2	50	0	43
	1:3	40	0	45
	1:4	30	0	46
Acetonitrile – H ₂ O	0:1	0	0	54
	2:1	62	46	66
	1:1	52	30	54
	1:2	51	27	46
	1:3	48	15	44
Trichloromethane	1:9	20	0	42
		68	64	71

Migration of C18-50 and C18-100 silica layers as compared to RP-18 W plates

a) RP-18 W, b) Nano-SIL C18-50, c) Nano-SIL C18-100 all plates with UV indicator

— methanol – water; - - - acetonitrile – water

Elution properties of MN RP plates in mixtures of methanol – water and acetonitrile – water

Further application examples can be found online in our application database at www.mn-net.com/apps

Modified silica layers

RP-18 W/UV₂₅₄ G A octadecyl-modified HPTLC silica layers

🔧 Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 2–10 µm, for preparative plates (1 mm thickness of layer) standard silica 60, pH stability 2–10, particle size 5–17 µm
- Indicator: acid-resistant product with a pale blue fluorescence for short-wave UV (254 nm), UV-absorbing substances appear as dark-blue to black spots on a light-blue background

🔧 Modification

- Partial octadecyl (C₁₈) modification, wettable with water, carbon content 14 %
- Order of polarity: silica > DIOL > NH₂ > CN > RP-2 > C18-50 > RP-18 W > C18-100

✅ Recommended application

- NP or RP separation with eluents from anhydrous solvents to mixtures with high concentrations of water (see table and figure on previous page), relative polarity of the eluent determines the polarity of the layer
- Aminophenols, barbiturates, preservatives, nucleobases, polycyclic aromatic hydrocarbons, steroids, tetracyclines, plasticizers (phthalates)

Ordering information

Plate size [cm]	4 x 8	5 x 10	5 x 20	10 x 10	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
-----------------	-------	--------	--------	---------	---------	---------	--------------------	-----------------------

Glass plates

Pack of [plates]		50	25	50	25		
RP-18 W/UV ₂₅₄		811073	811075	811072	811071	0.25 mm	UV ₂₅₄
Pack of [plates] (preparative TLC)					15		
RP-18 W/UV ₂₅₄					811074	1.00 mm	UV ₂₅₄

ALUGRAM® aluminum sheets

Pack of [plates]	50	50	50	25	25		
RP-18 W/UV ₂₅₄	818144	818152	818145	818147	818146	0.15 mm	UV ₂₅₄

RP-2/UV₂₅₄ G A "silanized silica" = dimethyl-modified standard silica layers

🔧 Technical characteristics

- Silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, pH stability 2–10, particle size 5–17 µm
- Indicator: acid-resistant product with a pale blue fluorescence for short-wave UV (254 nm), UV-absorbing substances appear as dark-blue to black spots on a light-blue background

🔧 Modification

- Silanized silica with dimethyl modification, carbon content 4 %
- Order of polarity: silica > DIOL > NH₂ > CN > RP-2 > C18-50 > RP-18 W > C18-100

✅ Recommended application

- Normal phase or reversed phase separation modes with purely organic, organic - aqueous or purely aqueous eluents
- Active plant constituents, steroids

Ordering information

Plate size [cm]	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
-----------------	---------	---------	--------------------	-----------------------

Glass plates

RP-2/UV ₂₅₄	811081	811082	0.25 mm	UV ₂₅₄
------------------------	--------	--------	---------	-------------------

ALUGRAM® aluminum sheets

RP-2/UV ₂₅₄		818171	0.15 mm	UV ₂₅₄
------------------------	--	--------	---------	-------------------

Nano-SIL CN G A cyano-modified HPTLC silica layers

🔧 Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, pH stability 2–8, particle size 2–10 μm
- Indicator: acid-resistant product with a pale blue fluorescence for short-wave UV (254 nm), UV-absorbing substances appear as dark-blue to black spots on a light-blue background


🔧 Modification

- Cyanopropyl modification, carbon content 5.5 %
- Order of polarity: silica > DIOL > NH₂ > CN > RP-2 > C18-50 > RP-18 W > C18-100

✅ Recommended application

- NP or RP separation modes depending on the polarity of the developing solvent (see figure below)
- Steroid hormones, phenols, preservatives

R_f values of different steroids as a function of eluent composition

Layer: Nano-SIL CN/UV

Polarity of the eluent governs the type of separation mechanism:

Eluent system petroleum ether (PE) – acetone (NP mode)

the higher the concentration of PE, the stronger are the adsorptive interactions of the steroids with the stationary phase

Eluent system acetone – water (RP mode)

the sequence of elution of the steroids is reversed, the most nonpolar compounds are most strongly retained


Separation of preservatives

MN Appl. No. 401440

Layer: Nano-SIL CN/UV
 Sample volume: 400 nL
 Eluent: ethanol – water – glacial acetic acid (20:80:0.2) with 0.1 mol/L tetraethylammonium chloride
 Migration distance: 73 mm in 30 min
 Detection: TLC scanner, UV 254 nm

Peaks:

1. Propyl p-hydroxybenzoate
2. Ethyl p-hydroxybenzoate
3. Methyl p-hydroxybenzoate
4. Benzoic acid
5. Sorbic acid

Ordering information

Plate size [cm]	4 x 8	10 x 10	10 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	25	25		

Glass plates

Nano-SIL CN/UV	811115	811116	0.20 mm	UV ₂₅₄
----------------	--------	--------	---------	-------------------

ALUGRAM® aluminum sheets

Nano-SIL CN/UV	818184		0.15 mm	UV ₂₅₄
----------------	--------	--	---------	-------------------

Nano-SIL NH₂ G A amino-modified HPTLC silica layers

✔ Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, pH stability 2–8, particle size 2–10 µm
- Indicator: acid-resistant product with a pale blue fluorescence for short-wave UV (254 nm), UV-absorbing substances appear as dark-blue to black spots on a light-blue background


✔ Modification

- Aminopropyl modification, carbon content 3.5 %
- Order of polarity: silica > DIOL > NH₂ > CN > RP-2 > C18-50 > RP-18 W > C18-100
- Layer can be wetted equally well with pure water as with organic solvents

✔ Recommended application

- Vitamins, sugars, steroids, purine derivatives, xanthenes, phenols, nucleotides and pesticides

Influence of eluent composition on the separation of nucleotides

Layer: Nano-SIL NH₂/UV
 Eluent: MeOH – H₂O according to fig. + 0.18 mol/L NaCl
 Migration distance: 7 cm

c-AMP, AMP: adenosine monophosphate
 ADP: adenosine diphosphate
 ATP: adenosine triphosphate


Separation of sugars

MN Appl. No. 401590

Layer: Nano-SIL NH₂/UV
 Sample volume: 0.5 µL
 Eluent: ethyl acetate – pyridine – water – glacial acetic acid (60:30:10:5, v/v/v/v)
 Migration distance: 80 mm in 45 min, double development
 Detection: dry layer at 160 °C for 5 min, TLC scanner, UV 254 nm

Peaks (0.1 % each):

1. Lactose
2. Saccharose
3. Galactose
4. Glucose
5. Fructose
6. Arabinose
7. Xylose
8. Ribose

Ordering information

Plate size [cm]	4 x 8	10 x 10	10 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	25	25		

Glass plates

Nano-SIL NH ₂ /UV	811111	811112	0.20 mm	UV ₂₅₄
------------------------------	--------	--------	---------	-------------------

ALUGRAM® aluminum sheets

Nano-SIL NH ₂ /UV	818182		0.15 mm	UV ₂₅₄
------------------------------	--------	--	---------	-------------------

Further application examples can be found online in our application database at www.mn-net.com/apps

Nano-SIL DIOL G diol-modified HPTLC silica layers

Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, pH stability 2–8, particle size 2–10 µm
- Indicator: acid-resistant product with a pale blue fluorescence for short-wave UV (254 nm), UV-absorbing substances appear as dark-blue to black spots on a light-blue background

Modification

- Diol modification, carbon content 5.5 %
- Order of polarity: silica > DIOL > NH₂ > CN > RP-2 > C18-50 > RP-18 W > C18-100
- Layer can be wetted equally well with pure water as with organic solvents

Recommended application

- Steroids, pesticides and plant constituents
- For critical separations an alternative to silica
- Since it is less sensitive to the water content of the environment, leads to more reproducible results compared to silica


Separation of herbicides

MN Appl. No. 401950

Layer: Nano-SIL DIOL/UV
 Sample volume: 2 µL
 Eluent: petroleum ether (40–60 °C) – acetone (80:20, v/v)
 Migration distance: 70 mm
 Detection: TLC scanner, 230 nm

Peaks:
 (0.07 % each in methanol)

1. Metoxuron
2. Monuron
3. Metobromuron


Ordering information

Plate size [cm]	10 x 10	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

Glass plates

Nano-SIL DIOL/UV	811120	0.20 mm	UV ₂₅₄
------------------	--------	---------	-------------------

Alox **G P A** aluminum oxide layers**🔧** Technical characteristics

- Aluminum oxide, mean pore size 60 Å, specific surface (BET) ~ 200 m²/g
- Inert organic binder
- Indicator: manganese-activated zinc silicate

✓ Recommended application

- Terpenes, alkaloids, steroids, aliphatic and aromatic compounds
- We recommend to activate aluminum oxide layers before use by heating 10 minutes at 120 °C


Separation of bisadducts of fullerenes

MN Appl. No. 401930

F. Djojo, A. Hirsch, Chem. Eur. J. 4 (1998), 344–356

Layer: ALUGRAM® Alox N/UV₂₅₄
 Eluent: toluene – ethyl acetate (95:5, v/v)
 Detection: UV, 254 nm

Compound	R _f values
Bis[bis(4-phenyloxazolin)methane]fullerene 1	0.14
Bis[bis(4-phenyloxazolin)methane]fullerene 2	0.26


Separation of lipophilic dyes

MN Appl. No. 403010

Layer: Alox-25 UV₂₅₄
 Sample volume: 1000 nL
 Eluent: toluene – cyclohexane (2:1, v/v)
 Migration distance: 108 mm in 15 min
 Detection: TLC scanner, UV 254 nm

Peaks:

1. Indophenol
2. Sudan red G
3. Sudan blue II
4. Butter yellow

Ordering information

Plate size [cm]	4 x 8	5 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Glass plates					
Pack of [plates]		100	25		
Alox-25 UV ₂₅₄		807021	807023	0.25 mm	UV ₂₅₄
Pack of [plates] (preparative TLC)					
Alox-100 UV ₂₅₄			807033	1.00 mm	UV ₂₅₄
POLYGRAM® polyester sheets					
Pack of [plates]	50	50	25		
Alox N/UV ₂₅₄	802021	802022	802023	0.20 mm	UV ₂₅₄
ALUGRAM® aluminum sheets					
Pack of [plates]		50	25		
Alox N/UV ₂₅₄		818024	818023	0.20 mm	UV ₂₅₄

Further application examples can be found online in our application database at www.mn-net.com/apps

Cellulose MN 300 G P A native fibrous cellulose layers

Technical characteristics

- Fiber length (95 %) 2–20 µm, average degree of polymerization 400–500, specific surface acc. to Blaine 15 000 cm²/g, ≤ 20 ppm Fe, 6 ppm Cu, 7 ppm P; CH₂Cl₂- extract ≤ 0.25 %; residue on ignition at 850 °C ≤ 1500 ppm

Recommended application

- Partition chromatography of polar substances such as amino acids, carboxylic acids or carbohydrates

Ordering information

Plate size [cm]	4 x 8	5 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Glass plates					
Pack of [plates]			25		
CEL 300-10			808013	0.10 mm	–
CEL 300-10 UV ₂₅₄			808023	0.10 mm	UV ₂₅₄
CEL 300-25			808033	0.25 mm	–
CEL 300-25 UV ₂₅₄			808043	0.25 mm	UV ₂₅₄
Pack of [plates] (preparative TLC)			20		
CEL 300-50			808053	0.50 mm	–
CEL 300-50 UV ₂₅₄			808063	0.50 mm	UV ₂₅₄
POLYGRAM® polyester sheets					
Pack of [plates]	50	50	25		
CEL 300	801011		801013	0.10 mm	–
CEL 300 UV ₂₅₄		801022	801023	0.10 mm	UV ₂₅₄
ALUGRAM® aluminum sheets					
Pack of [plates]	50	50	25		
CEL 300	818155		818153	0.10 mm	–
CEL 300 UV ₂₅₄		818157	818156	0.10 mm	UV ₂₅₄

Cellulose MN 400 (AVICEL®) G P microcrystalline cellulose layers

Technical characteristics


- Prepared by hydrolysis of high purity cellulose with HCl, average degree of polymerization 40–200

Recommended application

- Carboxylic acids, lower alcohols, urea and purine derivatives

Ordering information

Plate size [cm]	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	25		
Glass plates				
CEL 400-10	808072	808073	0.10 mm	–
POLYGRAM® polyester sheets				
CEL 400		801113	0.10 mm	–
CEL 400 UV ₂₅₄		801123	0.10 mm	UV ₂₅₄

Cellulose MN 300 PEI P PEI-impregnated cellulose ion exchange layers

🔧 Technical characteristics

- Fibrous cellulose impregnated with polyethyleneimine

✅ Recommended application

- Analysis of nucleic acids, and of mutagenic substances with the ³²P postlabelling procedure

Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

POLYGRAM® polyester sheets

CEL 300 PEI	801053	0.10 mm	–
CEL 300 PEI/UV ₂₅₄	801063	0.10 mm	UV ₂₅₄

Cellulose MN 300 AC P acetylated cellulose layers

🔧 Technical characteristics

- Fibrous cellulose with 10 % content of acetylated cellulose for reversed phase chromatography

✅ Recommended application

- Reversed phase chromatography

Ordering information

Plate size [cm]	Acetyl content	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]		25		

POLYGRAM® polyester sheets

CEL 300 AC-10 %	10 %	801033	0.10 mm	–
-----------------	------	--------	---------	---

Polyamid-6 P ε-polycaprolactame layers

🔧 Technical characteristics

- Polyamide 6 = nylon 6 = perlon = ε-aminopolycaprolactame
- Separation mechanism based on hydrogen bonds to amide groups of the polymer matrix as well as on ionic, dipole and electron donor-acceptor interactions

✅ Recommended application

- Natural compounds, phenols, carboxylic acids, aromatic nitro compounds and especially amino acids

Ordering information

Plate size [cm]	5 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	25		

POLYGRAM® polyester sheets

POLYAMID-6	803012	803013	0.10 mm	–
POLYAMID-6 UV ₂₅₄	803022	803023	0.10 mm	UV ₂₅₄

Further application examples can be found online in our application database at www.mn-net.com/apps

CHIRALPLATE ^G special layer enantiomer separation

Technical characteristics

- Reversed phase nano silica impregnated with Cu²⁺ ions and a chiral selector (proline derivative)
- Separation based on ligand exchange, i.e. formation of ternary mixed-ligand complexes with the Cu(II) ions, differences in the stability of the diastereomeric complexes cause chromatographic separation

Recommended application

- Enantiomer separation of amino acids, *N*-methylamino acids, *N*-formylamino acids, α -alkylamino acids, thiazolidine derivatives, dipeptides, lactones, α -hydroxycarboxylic acids


Enantiomer separation of amino acids

MN Appl. No. 400520

Quantitative determination (remission location curves) of TLC-separated enantiomers of *tert.*-leucine:

Layer: CHIRALPLATE
 Eluent: methanol – water (10:80, v/v)
 Detection: dip in 0.3 % ninhydrin solution
 quantification with scanner, 520 nm

- a) L-*tert.*-leucine
 b) L-*tert.*-leucine + 0.1 % D-*tert.*-leucine
 c) L-*tert.*-leucine + 1 % D-*tert.*-leucine
 d) external reference sample

Ordering information

Plate size [cm]	5 x 20	10 x 10	10 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Glass plates						
Pack of [plates]	4					
CHIRALPLATE			811056		0.25 mm	UV ₂₅₄
Pack of [plates]	50	25	25	25		
CHIRALPLATE	811057	811059	811055	811058	0.25 mm	UV ₂₅₄


SIL N-HR ^P unmodified standard silica layers

Technical characteristics

- High purity silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 5–17 μm, different binder system compared to SIL G results in different separation characteristics
- A special feature of the POLYGRAM[®] SIL N-HR is a higher gypsum content

Ordering information

Plate size [cm]	5 x 20	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50	25		
POLYGRAM[®] polyester sheets				
SIL N-HR/UV ₂₅₄	804022	804023	0.20 mm	UV ₂₅₄

SIL G-25 HR G special layer for aflatoxin separation

🔧 Technical characteristics

- High purity silica 60 with gypsum and a very small quantity of a polymeric organic binder; softer than the standard silica layer, i.e. spots can be scratched and the layer absorbs faster

✅ Recommended application

- Aflatoxins

Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

Glass plates

SIL G-25 HR	809033	0.25 mm	–
SIL G-25 HR/UV ₂₅₄	809043	0.25 mm	UV ₂₅₄

SIL G-25 Tenside G special layer for separation of surfactants

🔧 Technical characteristics

- Silica G impregnated with ammonium sulfate

✅ Recommended application

- Detergents, alkanesulfonates, polyglycols

Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

Glass plates

SIL G-25 Tenside	810063	0.25 mm	–
------------------	--------	---------	---

Nano-SIL PAH G special HPTLC silica layer for PAH analysis

🔧 Technical characteristics

- Nano silica 60, mean pore size 60 Å, specific surface (BET) ~ 500 m²/g, specific pore volume 0.75 mL/g, particle size 2–10 µm
- Impregnated with caffeine, an electron acceptor for PAH analysis based on charge-transfer complexes

✅ Recommended application

- 6 PAHs according to German drinking water specifications (TVO) in accordance with German standard DIN 38407 part 7

Ordering information

Plate size [cm]	10 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	50		

Glass plates

Nano-SIL PAH	811051	0.20 mm	–
--------------	--------	---------	---

Further application examples can be found online in our application database at www.mn-net.com/apps

Layers for special TLC separations

IONEX ^P special mixed layers of silica with ion exchange resins

IONEX-25 SA-Na:

- Mixture of silica and a strongly acidic cation exchanger coated to polyester sheets

IONEX-25 SB-AC:

- Mixture of silica and a strongly basic anion exchanger coated to polyester sheets
- Both layers contain an inert organic binder

✓ Recommended application

- Amino acids, e.g., in protein and peptide hydrolyzates, in seeds and fodder, in biological fluids; for racemate separation in peptide syntheses, for the separation of nucleic acid hydrolyzates, aminosugars, amino acids, antibiotics, inorganic phosphates, cations and other compounds with ionic groups

Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

POLYGRAM[®] polyester sheets

IONEX-25 SA-Na	strongly acidic cation exchanger	806013	0.20 mm	–
IONEX-25 SB-AC	strongly basic anion exchanger	806023	0.20 mm	–

Mixed layers for TLC ^G

Alox/CEL-AC-Mix-25:

- Mixed layer of aluminum oxide G and acetylated cellulose, recommended for separation of PAH

SILCEL-Mix-25:

- Mixed layer of cellulose and silica, recommended for separation of preservatives and other antimicrobial compounds

Ordering information

Plate size [cm]	20 x 20	Thickness of layer	Fluorescent indicator
Pack of [plates]	25		

Glass plates

Alox/CEL-AC-Mix-25	810053	0.25 mm	–
SILCEL-Mix-25 UV ₂₅₄	810043	0.25 mm	UV ₂₅₄

Further application examples can be found online in our application database at www.mn-net.com/apps

Chromatography papers

Chromatography papers

- Paper chromatography is the oldest chromatographic technique separation due to partition of the analytes between special paper grades and the mobile phase, which penetrates the paper by capillary action ascending.
- Descending and circular techniques are possible

Please note

- Always treat chromatography papers with care
- Never touch them with fingers, because this will contaminate the surface
- Do not bend them sharply, because this will decrease the capillary action (preferably store them flat)

Direction

- Chromatography papers possess a preferred direction of the fibers with higher absorption properties (with our sheets 58 x 60 cm, the longer edge)
- We recommend to use them in the direction of higher absorption

Ordering information

Code	Weight [g/m ²]	Thickness [mm]	Description	Flow rate	Size [cm]	Pack of	REF
MN 214	140	0.28	smooth	90–100 mm/30 min	58 x 60	100 sheets	817001
MN 218	180	0.36	smooth	90–100 mm/30 min	58 x 60	100 sheets	817002
MN 260	90	0.20	smooth	120–130 mm/30 min	58 x 60	100 sheets	817003
MN 261	90	0.18	smooth	90–100 mm/30 min	58 x 60	100 sheets	817004
MN 827	270	0.70	soft carton	130–140 mm/10 min	58 x 60	100 sheets	817005
MN 866	650	1.70	soft carton	100–120 mm/10 min	38 x 38	100 sheets	817006
MN 866	650	1.70	soft carton	100–120 mm/10 min	80 x 80	100 sheets	817007
MN 214 ff	140	0.28	MN 214 defatted *	90–100 mm/30 min	56 x 58	100 sheets	817008

* This paper is extracted with organic solvents.

For further papers, filters and membranes, feel free to ask for our catalog "Filtration".

Accessories


• Beside ready-to-use layers for thin layer chromatography also accessories are required

• Selection of accessories for reliable separation in TLC

Ordering information

Designation	Pack of	REF
Simultaneous developing chamber for TLC, 20 x 20 cm	1	814019
Simultaneous developing chamber for TLC, 10 x 10 cm	1	814018
Developing chambers for TLC micro-sets	4	814021
Glass laboratory sprayer with rubber bulb	1	814101
Glass capillaries 1 µL	3 x 50	814022
Rubber caps for capillaries	2	814102
Plastic syringe, 1 mL content with graduation	1	814104
Spotting guides	2	814023
Measuring cylinders, glass, 10 mL content	2	814024
MN ALUGRAM® scissors, ground blade, black handle	1	818666
Filter paper MN 713, 15 x 21 cm	100	814103
Folded filters MN 615 1/4, 11 cm diameter	100	531011
Chromatography paper MN 260, 7.5 x 17 cm (for chamber saturation)	100	814030

Visualization reagents

• Small selection of frequently used spray reagents for post chromatographic detection reactions in TLC suited for spraying or dipping TLC plates

• A detailed description of many more detection procedures for TLC is available on request

Ordering information

Spray reagent	Solvent	Detection of	Pack of	REF
Aniline phthalate	2-propanol – ethanol (1:1)	reducing sugars, oxohalic acids	100 mL	814919
Bromocresol green	2-propanol	organic acids	100 mL	814920
Reagent for caffeine detection	water – acetone	caffeine	100 mL	814401
2',7'-Dichlorofluorescein	2-propanol	lipids (saturated, unsaturated)	100 mL	814921
4-(Dimethylamino)-benzaldehyde	2-propanol	terpenes, sugars, steroids	100 mL	814922
Reagent according to Dragendorff-Munier	water	alkaloids and other nitrogen compounds	100 mL	814402
Iron(III) chloride	water	phenolic compounds e.g., acetylsalicylic acid, paracetamol	100 mL	814403
Potassium hexacyanoferrate(III)	water		100 mL	814404
Molybdato-phosphoric acid	ethanol	lipids, sterols, steroids, reducing compounds	100 mL	814302
Ninhydrin	ethanol	amino acids, amines and amino sugars	100 mL	814203
Rhodamine B	ethanol	lipids	100 mL	814923
Rubeanic acid	ethanol	heavy metal cations	100 mL	814206

These products contain harmful substances which must be specially labeled as hazardous. For detailed information please see SDS.

Fluorescent indicators

UV indicators with efficient radiation for short-wave as well as long-wave UV ranges

• UV₂₅₄: manganese-activated zinc silicate with absorption maximum at 254 nm, green fluorescence, relatively susceptible towards acids: its fluorescence can be completely quenched by acidic solvents

• UV₃₆₆: inorganic fluorescent pigment with absorption maximum at 366 nm, blue fluorescence

Ordering information

	Composition	Absorption maximum	Color of fluorescence	Pack of 100 g
Fluorescent indicator UV ₂₅₄	manganese-activated zinc silicate	254 nm	green	816710.01
Fluorescent indicator UV ₃₆₆	inorganic fluorescent pigment	366 nm	blue	816720.01

Silica adsorbent for TLC

Pore size 60 Å, pore volume 0.75 mL/g, specific surface (BET) ~ 500 m²/g, pH 7 for a 10 % aqueous suspension

- Silica G: standard grade, particle size 2–20 µm, Fe < 0.02 %, Cl < 0.02 %, 13 % gypsum as binder
- Silica N: standard grade, particle size 2–20 µm, Fe < 0.02 %, Cl < 0.02 %, no binder
- Silica G-HR: high purity grade, particle size 3–20 µm, Fe < 0.002 %, Cl < 0.008 %, gypsum as binder
- Silica P: preparative grade, particle size 5–50 µm, Fe < 0.02 %, Cl < 0.02 %, organic binder
- Silica P with gypsum: preparative grade, particle size 5–50 µm, Fe < 0.02 %, Cl < 0.02 %, gypsum as binder

Ordering information

Designation	Fluorescent indicator	1 kg	5 kg
Silica G	–	816310.1	816310.5
Silica G/UV ₂₅₄	UV ₂₅₄	816320.1	816320.5
Silica N	–	816330.1	816330.5
Silica N/UV ₂₅₄	UV ₂₅₄	816340.1	816340.5
Silica G-HR	–	816410.1	816410.5
Silica P/UV ₂₅₄	UV ₂₅₄	816380.1	816380.5
Silica P/UV ₂₅₄ with gypsums	UV ₂₅₄	816400.1	816400.5

Polyamid adsorbent for TLC

Polyamide 6 = nylon 6 = perlon = ε-polycaprolactame

Ordering information

Designation	Fluorescent indicator	1 kg
Polyamid-DC 6	–	816610.1
Polyamid-DC 6 UV ₂₅₄	UV ₂₅₄	816620.1

Cellulose MN 301 native fibrous cellulose

- Standard grade, fiber length (95 %) 2–20 µm
- Average degree of polymerization 400–500, specific surface acc. to Blaine 15 000 cm²/g
- ≤ 20 ppm Fe, 6 ppm Cu, 7 ppm P, CH₂Cl₂ extract ≤ 0.25 %, residue on ignition at 850 °C ≤ 1500 ppm

Ordering information

Designation	1 kg	5 kg
Cellulose MN 301	816250.1	816250.5

Distributed By

Greyhound Chromatography and Allied Chemicals
 6 Kelvin Park
 Birkenhead
 Merseyside, CH41 1LT

Tel: 0151 649 4000 Fax: 0151 649 4001
 Email: info@greyhoundchrom.com
 Web: <https://www.greyhoundchrom.com>