



Spot the difference in the new SGE lineup!

SGE enhanced inlet liner range:

- Easy to choose.
- Easy to use.
- Confidence in your analysis.



Do you know how to select the right inlet liner for your analysis?

Choosing the right inlet liner and injection parameter can increase peak areas and **reduce detection limits by up to 300 %¹**.

Customer research conducted by SGE found that a significant number of GC users don't understand the importance of inlet liner selection, nor how it contributes to their analysis. The enhanced SGE inlet liner range aims to make it simple for all gas chromatographers to select the right liner.



Spot the difference in the new SGE lineup!

Easy to choose

- Color coded by geometry to simplify your selection.

Easy to use

- Contain o-rings so you're ready to go.

Confidence in your analysis

- Certified deactivation gives you confidence in your analysis.

Easy to choose

For ease of selection so you can optimize your results, SGE has color coded the liners by injection type. Color coding and sample types are listed in the table below. SGE has also created a liner selection tool that can be found on our website: www.sge.com/linertool



| Color | Injection Technique | Sample Types | Liner Geometry | How the Geometry Works |
|------------|---------------------|--|-----------------------|---|
| Dark Green | Splitless | <ul style="list-style-type: none"> • Trace level analyses. • Active compounds. | Taper / Gooseneck | <ul style="list-style-type: none"> • A bottom taper focuses sample onto the head of the column and minimizes sample contact with metal parts of the inlet. • Remember – the addition of quartz wool to your inlet liner promotes mixing of analytes, aids the vaporization of liquid samples, works as a trap to collect non-volatile residue in the sample (i.e. protects capillary column from 'dirty' samples). |
| Blue | Split | <ul style="list-style-type: none"> • General purpose. • Concentrated samples. • Dirty samples. | FocusLiner™ | <ul style="list-style-type: none"> • Ensures glass wool remains in the correct position in the liner. • Excellent reproducibility results from the wiping of the sample from the syringe needle and the prevention of droplet formation. • Minimizes high molecular weight discrimination. |
| Aqua | Splitless | <ul style="list-style-type: none"> • Trace level analyses. • Dirty samples. • Wide boiling point range. | Taper Focus | <ul style="list-style-type: none"> • Bottom taper focuses sample onto the head of the column and minimizes contact with metal parts of the inlet. • Ensures quartz wool remains in the correct position in the liner. • Excellent reproducibility results from the wiping of the sample from the syringe needle and the prevention of droplet formation. |
| Orange | Direct | <ul style="list-style-type: none"> • Trace level analyses. • Active compounds. | ConnectTite™ | <ul style="list-style-type: none"> • ConnectTite™ liners facilitate maximum transfer of sample to GC column and inhibit sample degradation due to hot metal components inside inlet. • Systems equipped with electronic pressure control require a hole in the liner body to maintain system gas flows. • ConnectTite™ liners that have a hole near the bottom are best suited for analyses where a tailing solvent peak could affect early eluting compounds. ConnectTite™ liners with a hole at the top of the liner will improve your analysis with aqueous injections or where compounds of interest elute away from the solvent peak. |
| Purple | Split Splitless | <ul style="list-style-type: none"> • General purpose. • Concentrated samples. • Dirty samples (only if quartz wool is present) • Gaseous samples (also Purge and Trap, Headspace). | Straight | <ul style="list-style-type: none"> • Straight liners facilitate higher split flows. • Narrow bore straight liners facilitate fast GC work. • Small injection volumes of less than 0.5 µL are best used with a narrow bore. • Narrow bore straight liners improve focussing of gaseous samples (purge & trap and headspace). |
| Yellow | Splitless LVI | <ul style="list-style-type: none"> • Trace level analyses. • Low boiling point compounds. • Active compounds. | Double Taper | <ul style="list-style-type: none"> • Bottom taper minimizes contact with metal parts of the inlet and focuses sample onto the head of the column. • Top taper aids in minimizing sample flashback. |
| Grey | PTV LVI | <ul style="list-style-type: none"> • Trace level analyses. • Large volume injections. | PTV/LVI | <ul style="list-style-type: none"> • PTV and LVI liners generally have sintered glass beads or powder to increase the surface area and trap non-volatile residue. • PTV liners use baffles or a wisp of quartz wool to aid in vaporization of samples and retain droplets during low temperature injections. • Side hole needles are recommended for these techniques to ensure effective distribution of sample within the liner. |

SGE has tested the inks used on our new inlet liners to confirm:

- The ink does not affect the chromatography of the sample.
- The ink does not deteriorate when in contact with strong solvents such as dichloromethane.
- The ink is not affected by prolonged exposure to temperatures above 400 °C.

Indication of no interference of inlet liner print or packaging with Polynuclear Aromatic Hydrocarbons (PAH) Analysis on BPX50.

| | |
|---------------------------------|------------------------|
| Inlet | 300 °C |
| Transfer Liner | 300 °C |
| Initial Temperature | Initial hold 1 minute |
| Rate 1 | 35 °C/minute |
| | Temperature 120 °C |
| | Hold 0.5 minutes |
| Rate 2 | 8 °C/minute |
| | Temperature 200 °C |
| | Hold 0 minutes |
| Rate 3 | 11 °C/minute |
| | Temperature 270 °C |
| | Hold 0 minutes |
| Rate 4 | 2 °C/minute |
| | Temperature 300 °C |
| | Hold 0 minutes |
| Rate 5 | 40 °C/minute |
| | Temperature 320 °C |
| | Hold 4 minutes |
| MS – Source Temperature | 260 °C |
| Scan | 50-600 amu |
| High Pressure Injection (45psi) | Splitless for 1 minute |
| Constant Velocity | 52 cm/sec |



1 µL injected of 20 ng/µL of:

- 1 naphthalene
- 2 2-methylnaphthalene
- 3 acenaphthene
- 4 acenaphthylene
- 5 fluorene
- 6 phenanthrene
- 7 anthracene
- 8 fluoranthene
- 9 pyrene
- 10 benzo(a)anthracene
- 11 chrysene
- 12 benzo(b)fluoranthene
- 13 benzo(j)fluoranthene
- 14 benzo(k)fluoranthene
- 15 benzo(a)pyrene
- 16 indeno(1,2,3-cd)pyrene
- 17 dibenzo(a,h)anthracene
- 18 benzo(g,h,i)perylene

Note:

- New enhanced SGE inlet liner was installed into the Shimadzu QP-2010.
- Inlet had a used, conditioned septum.
- Inlet liner was handled with gloves.
- Inlet liner was not conditioned.
- Analysis was performed as soon as the injector reached 300 °C.

Easy to use

SGE's enhanced inlet liners now come as a new, complete, packaged solution:

- 1, 5 and 25 packs.
- Complete with instrument appropriate o-rings or sealing rings.
- Each pack supplied with quality assurance test results.
- Color coded instrument stickers to identify which inlet liner type is in your instrument.
- SGE blister packs are now perforated enabling easy division of the 5 and 25 packs while maintaining blister integrity.
- 25 packs come in a re-usable container, with a range of attractive designs, that will be handy around the lab.



Confidence in your analysis


SGE standard inlet liner deactivation

SGE currently tests every batch of inlet liners for activity using the EPA 8081B method. This standard method ensures that each batch of inlet liners has less than 3 % Endrin breakdown. Now, SGE is validating this quality assurance by including a batch certificate with every pack.

When deactivation REALLY matters

Single pack deactivation certification options:

- Sometimes there is a need for fully traceable inlet liner certification. Customers who are ISO accredited or follow GLP will benefit from SGE's certified single packs. Add CERT when ordering.
- MS ready liners in single packs conditioned and MS tested post deactivation so they are ready for use straight out of the pack. Add MS when ordering.



To select the right inlet liner for your instrument, injection and sample - use our selection tool: www.sge.com/linertool

AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd
Toll Free: 1800 800 167
Tel: +61 (0) 3 9837 4200
Fax: +61 (0) 3 9874 5672
Email: support@sge.com

EUROPE

SGE Europe Ltd
European Head Office
Toll Free: 00800 2790 8999
Toll Free Fax: 00800 2626 2609
Tel: +44 1908 568 844
Fax: +44 1908 566 790
Tel France: +33 1 69 29 80 90
Fax France: +33 1 69 29 09 25
Tel Germany: +49 (0) 6155 / 60746 0
Fax Germany: +49 (0) 6155 / 60746 50
Email: europe@sge.com

CHINA

SGE Shanghai Representative Office
Tel: +86 21 6407 9382
Fax: +86 21 6407 9386
Email: china@sge.com

INDIA

SGE Laboratory Accessories Pvt Ltd
Tel: +91 22 24715896
Fax: +91 22 24716592
Email: sgeindia@vsnl.com

JAPAN

SGE Japan Inc
Tel: +81 45 222 2885
Fax: +81 45 222 2887
Email: japan@sge.com

MIDDLE EAST

SGE Gulf
Tel: +971 6 557 3341
Fax: +971 6 557 3541
Email: gulfsupport@sge.com

UNITED STATES OF AMERICA

SGE Incorporated
Toll Free: (800) 945 6154
Tel: +1 512 837 7190
Fax: +1 512 836 9159
Email: usa@sge.com

