

PRODUCT UPDATES FROM WELLINGTON LABORATORIES

WELLINGTON Reporter


April 13, 2014


NEW PRODUCTS


Mass-Labelled Telomer Sulfonates M2-4:2FTS & M2-8:2FTS

Polyfluorinated sulfonates (commonly referred to as telomer sulfonates) have been incorporated into aqueous film-forming foams (AFFFs) which are used to extinguish hydrocarbon-fueled fires. The effectiveness of the developed formulations for fighting Class B flammable liquid fires resulted in their widespread use and subsequent release into the environment.

Analytical laboratories developing methods for monitoring the levels of these compounds in environmental samples have reported matrix effects leading to ionization enhancement during LCMS analysis of water samples. For this reason, **Wellington** has prepared mass-labelled ($^{13}C_2$) reference standard solutions of 4:2FTS (**M2-4:2FTS**) and 8:2FTS (**M2-8:2FTS**) to complement our existing telomer sulfonate standards and aid researchers in the generation of accurate data.







Sodium 1H,1H,2H,2H-perfluoro-[1,2-¹³C₂]decane sulfonate

	Catalogue Number	Product (methanol)	Qty/Conc	
NEW	M2-4:2FTS	Sodium 1H,1H,2H,2H-perfluoro-[1,2-13C2]hexane sulfonate	1.2 ml	50 µg/ml
	M2-6:2FTS	Sodium 1H,1H,2H,2H-perfluoro-[1,2- ¹³ C ₂]octane sulfonate	1.2 ml	50 µg/ml
NEW	M2-8:2FTS	Sodium 1H,1H,2H,2H-perfluoro-[1,2-13C2]decane sulfonate	1.2 ml	50 µg/ml

The following native telomer sulfonates are also available...

Catalogue Number	Product (methanol)	Qty/Conc	
4:2FTS	Sodium 1H,1H,2H,2H-perfluorohexane sulfonate	1.2 ml	50 µg/ml
6:2FTS	Sodium 1H,1H,2H,2H-perfluorooctane sulfonate	1.2 ml	50 µg/ml
8:2FTS	Sodium 1H,1H,2H,2H-perfluorodecane sulfonate	1.2 ml	50 µg/ml

Distributed Throughout Europe and Middle East By-





Greyhound Chromatography & Allied Chemicals 6 Kelvin Park, Birkenhead, Merseyside, CH41 1LT, U.K. Tel: (+44)-0-151-649-4000 Fax: (+44)-0-151-649-4001 E-mail: info@ greyhoundchrom.com Web: www.greyhoundchrom.com