



Part No: APG90EB

# Compact pump without pressure sensor, with 10 ml/min ceramic pump head

# For HPLC and dosing applications

With its small footprint, the compact dual piston pump AZURA P2.1S can be used for a wide range of laboratory tasks. A pressure rating of up to 400 bar and chemical resistance to a wide range of eluents make it the perfect choice for LC and dosing applications. Pump heads with maximum flow rates of 10 and 50 ml/min are available.

Moreover, its exchangeable pump head allows adaptation of the pump for delivery of aggressive media and bioinert applications. The pump head can be heated or cooled with optional accessories.



AZURA P 2.1S pump

## **KEY FEATURES**



compact



Made in Germany

# **TECHNICAL DATA**

# **General Information**

Max. flow rate	10 ml/mir	า											
Flow rate range	0.001-10	ml/min											
Flow rate increment	0.001												
Maximum delivery pressure [psi]	5080 psi												
Maximum delivery pressure [bar]	350 bar												
Maximum delivery pressure [MPa]	35 MPa												
Pressure range													
		400											
		300											
	e [bar]												
	pressure [bar]	200											
	ď	100											
		0			0								
		(	0		2			e [ml/min]	6		8		0
				10 ml pu	mp head o	eramic	reco	mmended	r	ecommen	ided (conti	inuous)	
Pump head materials	Ceramic												
Maximum viscosity	100 ср												
Liquid temperature range	4-60 °C (39.2-140 °F)												
Gradient	Isocratic	Isocratic											
Leak management	no	no											
Wetted materials	GFP (graphite fiber reinforced PTFE), ruby, sapphire, aluminum oxide (Al <sub>2</sub> O <sub>3</sub> ), ceramic, FKM, PEEK, Zirconium oxide												

#### **Detailed information**

Best working conditions	0.1-8.0 ml/min
Continuous working conditions	0.1-4.0 ml/min
Flow rate accuracy	± 5 %
Flow rate accuracy conditions	measured at 5-50% of flow range using ethanol/water 10:90
Flow rate precision	< 0.5% RSD
Flow rate precision conditions	measured at 1 / 5 ml/min using ethanol/water 10:90
Pump type	Ultra-compact high pressure pump
Pulsation compensation	no
Active piston seal wash	no
Piston seal washing	yes
System protection	I <sub>min</sub> and I <sub>max</sub> adjustable
Pump head inlet (standard)	1/4"-28 UNF (flat bottom)
Pump head outlet (standard)	10−32 UNF (coned)

# Communication

Interfaces	display
Control	RS-232; analog; button on device
Analog inputs	0-10V

# General

Power supply	100-240 V; 50-60 Hz
Dimensions	121 x 129 x 220 mm
Weight	2.3 kg
Leak sensor	no
Ambient conditions	temperature range: 10–40 °C; 50–104 °F air humidity: below 90 % humidity (non condensing)

### **ACCESSORIES & SPARE PARTS**



Pump head 10 ml Pump head 10 ml, stainless steel

Article No.: AHB40



Accessory kit AZURA P 2.1S and P 4.1S included in shipkit

Article No.: FPG



Accessory kit AZURA small devices cable kit for all devices with footprint "S" included in shipkit

Article No.: FZA01



Pump head cooling and heating device

for 10 and 50 ml pump heads

Article No.: A2035-1



**DOCUMENTS** AND MORE INFORMATION

# AZURA Pump P 4.1S / P 2.1S User Manual EN (V6870)

Manual

English

PDF

2.0 MB

Distributed By:



Greyhound Chromatography and Allied Chemicals

6 Kelvin Park Birkenhead Merseyside CH41 1LT, UK

Tel: +44 (0)151 649 4000 Fax: +44 (0)151 649 4001

Email: info@greyhoundchrom.com Web: www.greyhoundchrom.com