

## WebSeal Microplate Sealing Mats DMSO Resistant

- Manufactured of pure silicone with extensive PTFE layer
- DMSO resistant for long term storage

### WebSeal System

PART NO.	DESCRIPTION	COLOUR	PK SIZE
WSM-2F	96 round well, domed base silicone/PTFE layer sealing mat	Blue/red	5/pk
WSM-2FBF	96 round well, flat base silicone/PTFE layer sealing mat	Blue/red	5/pk
WSM-3SF	96 square well, silicone/PTFE layer sealing mat	Blue/red	5/pk
WSM-5F	384 square well, silicone/PTFE layer sealing mat	Blue/red	5/pk
WSM-5RF	384 round well, silicone/PTFE layer sealing mat	Blue/red	5/pk

## WebSeal Mat Applicator

### WebSeal System

PART NO.	DESCRIPTION	COLOUR	PK SIZE
WSA-1	Hand held mat applicator	Black	each

## Minivap™

Sample concentrator, evaporating 500µL methanol in less than 10 minutes

### WebSeal System

PART NO.	DESCRIPTION	PK SIZE
CLS-229203	MiniVap Sample Concentrator	each

## UltraVap™

Programmable high speed sample concentrator with patented spiral needles, evaporating 500µL methanol in 6 minutes

### WebSeal System

PART NO.	DESCRIPTION	PK SIZE
CLS-229070	UltraVap 96 Well High Speed Sample Concentrator	each
CLS-229074	Replacement Spiral Needle Kit, with Fitting Tool	10/pk