

WebSeal Microplate Sealing System

- Manufactured from pure silicone with PTFE coating on underside
- Superior resealability after multiple injections
- Excellent chemical compatibility
- Dry heat autoclavable
- Eliminates cross-contamination well to well
- Available pre-slit for easy penetration and reduced vacuum formation

WebSeal Sealing Mats

Silicone/PTFE Coating

WebSeal System

| PART NO. | DESCRIPTION | COLOUR | PK SIZE |
|----------|---|---------|---------|
| WSM-2L | 96 round well, dome base silicone/PTFE coated sealing mat | Lt blue | 5/pk |
| WSM-2 | 96 round well, dome base silicone/PTFE coated sealing mat | Dk blue | 5/pk |
| WSM-2X | 96 round well, dome base, silicone/PTFE coated sealing mat pre-slit | Dk blue | 5/pk |
| WSM-2FB | 96 round well, flat base, silicone/PTFE coated sealing mat | Lt blue | 5/pk |
| WSM-2FBP | 96 round well, flat base, silicone/PTFE coated sealing mat pre-slit | Pink | 5/pk |
| WSM-2FBX | 96 round well, flat base, silicone/PTFE coated sealing mat pre-slit | Lt blue | 5/pk |
| WSM-3S | 96 square well, silicone/PTFE coated sealing mat | Lt blue | 5/pk |
| WSM-3SXY | 96 square well, silicone/PTFE coated sealing mat pre-slit | Yellow | 5/pk |
| WSM-3SX | 96 square well, silicone/PTFE coated sealing mat pre-slit | Lt blue | 5/pk |
| WSM-5 | 384 round well, silicone/PTFE coated sealing mat | Red | 5/pk |
| WSM-5X | 384 round well, silicone/PTFE coated sealing mat pre-slit | Lt blue | 5/pk |
| WSM-5R | 384 square well, silicone/PTFE coated sealing mat | Lt blue | 5/pk |
| WSM-5RX | 384 square well, silicone/PTFE coated sealing mat pre-slit | Lt blue | 5/pk |

WebSeal SYSTEM
 WebSeal SYSTEM
 WebSeal SYSTEM
 WebSeal SYSTEM

