

Organic Acids

Standards for Ion Chromatography

- Starting materials are available in organic acids or salts
- Analyte concentrations are reported for the anion component of the acid starting material
- Highest purity and stability
- Products included in our ISO 9001:2015, ISO/IEC 17025:2005 and ISO 17034:2016 scope of accreditation



Organic Acids					
Element	Source	Concentration	Matrix	Volume	Part #
Acetate	Sodium Acetate	100 ug/mL	H ₂ O	100 mL	IC-AC-100
Acetate	Sodium Acetate	100 ug/mL	H ₂ O	250 mL	IC-AC-250
Acetate	Sodium Acetate	100 ug/mL	H ₂ O	500 mL	IC-AC-500
Acetate	Sodium Acetate	1,000 ug/mL	H ₂ O	100 mL	IC-AC-M-100
Acetate	Sodium Acetate	1,000 ug/mL	H ₂ O	250 mL	IC-AC-M-250
Acetate	Sodium Acetate	1,000 ug/mL	H ₂ O	500 mL	IC-AC-M-500
Butyrate	Sodium Butyrate	100 ug/mL	H ₂ O	100 mL	IC-BTY-100
Butyrate	Sodium Butyrate	100 ug/mL	H ₂ O	250 mL	IC-BTY-250
Butyrate	Sodium Butyrate	100 ug/mL	H ₂ O	500 mL	IC-BTY-500
Butyrate	Sodium Butyrate	1,000 ug/mL	H ₂ O	100 mL	IC-BTY-M-100
Butyrate	Sodium Butyrate	1,000 ug/mL	H ₂ O	250 mL	IC-BTY-M-250
Butyrate	Sodium Butyrate	1,000 ug/mL	H ₂ O	500 mL	IC-BTY-M-500
Formate	Sodium Formate	100 ug/mL	H ₂ O	100 mL	IC-FM-100
Formate	Sodium Formate	100 ug/mL	H ₂ O	250 mL	IC-FM-250
Formate	Sodium Formate	100 ug/mL	H ₂ O	500 mL	IC-FM-500
Formate	Sodium Formate	1,000 ug/mL	H ₂ O	100 mL	IC-FM-M-100
Formate	Sodium Formate	1,000 ug/mL	H ₂ O	250 mL	IC-FM-M-250
Formate	Sodium Formate	1,000 ug/mL	H ₂ O	500 mL	IC-FM-M-500



Distributed By:
 Greyhound Chromatography and Allied Chemicals
 6 Kelvin Park
 Birkenhead, Merseyside, CH41 1LT, UK

Order: Tel: +44 (0) 151 649 4000
Email: sales@greyhoundchrom.com
Online: www.greyhoundchrom.com

Organic Acids (cont'd)

Element	Source	Concentration	Matrix	Volume	Part #
Glycolate	Sodium Glycolate	100 ug/mL	H ₂ O	100 mL	IC-GLY-100
Glycolate	Sodium Glycolate	100 ug/mL	H ₂ O	250 mL	IC-GLY-250
Glycolate	Sodium Glycolate	100 ug/mL	H ₂ O	500 mL	IC-GLY-500
Glycolate	Sodium Glycolate	1,000 ug/mL	H ₂ O	100 mL	IC-GLY-M-100
Glycolate	Sodium Glycolate	1,000 ug/mL	H ₂ O	250 mL	IC-GLY-M-250
Glycolate	Sodium Glycolate	1,000 ug/mL	H ₂ O	500 mL	IC-GLY-M-500
Lactate	Sodium Lactate	100 ug/mL	H ₂ O	100 mL	IC-LAC-100
Lactate	Sodium Lactate	100 ug/mL	H ₂ O	250 mL	IC-LAC-250
Lactate	Sodium Lactate	100 ug/mL	H ₂ O	500 mL	IC-LAC-500
Lactate	Sodium Lactate	1,000 ug/mL	H ₂ O	100 mL	IC-LAC-M-100
Lactate	Sodium Lactate	1,000 ug/mL	H ₂ O	250 mL	IC-LAC-M-250
Lactate	Sodium Lactate	1,000 ug/mL	H ₂ O	500 mL	IC-LAC-M-500
Oxalate	Sodium Oxalate	100 ug/mL	H ₂ O	100 mL	IC-OX-100
Oxalate	Sodium Oxalate	100 ug/mL	H ₂ O	250 mL	IC-OX-250
Oxalate	Sodium Oxalate	100 ug/mL	H ₂ O	500 mL	IC-OX-500
Oxalate	Sodium Oxalate	1,000 ug/mL	H ₂ O	100 mL	IC-OX-M-100
Oxalate	Sodium Oxalate	1,000 ug/mL	H ₂ O	250 mL	IC-OX-M-250
Oxalate	Sodium Oxalate	1,000 ug/mL	H ₂ O	500 mL	IC-OX-M-500
Propionate	Sodium Propionate	100 ug/mL	H ₂ O	100 mL	IC-PRO-100
Propionate	Sodium Propionate	100 ug/mL	H ₂ O	250 mL	IC-PRO-250
Propionate	Sodium Propionate	100 ug/mL	H ₂ O	500 mL	IC-PRO-500
Propionate	Sodium Propionate	1,000 ug/mL	H ₂ O	100 mL	IC-PRO-M-100
Propionate	Sodium Propionate	1,000 ug/mL	H ₂ O	250 mL	IC-PRO-M-250
Propionate	Sodium Propionate	1,000 ug/mL	H ₂ O	500 mL	IC-PRO-M-500
Valerate	Sodium Valerate	100 ug/mL	H ₂ O	100 mL	IC-VAL-100
Valerate	Sodium Valerate	100 ug/mL	H ₂ O	250 mL	IC-VAL-250
Valerate	Sodium Valerate	100 ug/mL	H ₂ O	500 mL	IC-VAL-500
Valerate	Sodium Valerate	1,000 ug/mL	H ₂ O	100 mL	IC-VAL-M-100
Valerate	Sodium Valerate	1,000 ug/mL	H ₂ O	250 mL	IC-VAL-M-250
Valerate	Sodium Valerate	1,000 ug/mL	H ₂ O	500 mL	IC-VAL-M-500



Distributed By:

Greyhound Chromatography and Allied Chemicals
 6 Kelvin Park
 Birkenhead, Merseyside, CH41 1LT, UK

Order: Tel: +44 (0) 151 649 4000

Email: sales@greyhoundchrom.com

Online: www.greyhoundchrom.com