



NEW PERFLUORINATED PRODUCTS **Native Telomer Acrylates and Acetates**


The production, use, and widespread distribution of fluorinated polymers has raised exposure concerns relating to both residual compounds (telomer alcohols, acrylates, and acetates) as well as potential degradation products (perfluoroalkylcarboxylic acids). Indeed, fluorotelomer acrylates have been found in indoor air samples presumably from residual off-gasing of fluorochemical products. In order to aid researchers in the analysis of these compounds, **Wellington** is adding native telomer acrylate (**8:2FTAcr** and **10:2FTAcr**) and native telomer acetate (**8:2FTOAc** and **10:2FTOAc**) reference standards to our line of perfluorinated compounds.




8:2FTAcr



10:2FTAcr




8:2FTOAc




10:2FTOAc

Catalogue Number	Product (isooctane)	Qty/Conc
8:2FTAcr	1H,1H,2H,2H-Perfluorodecyl acrylate	1.2 ml 50 µg/ml
10:2FTAcr	1H,1H,2H,2H-Perfluorododecyl acrylate (97%)	1.2 ml 48.5 µg/ml
8:2FTOAc	1H,1H,2H,2H-Perfluorodecyl acetate	1.2 ml 50 µg/ml
10:2FTOAc	1H,1H,2H,2H-Perfluorododecyl acetate	1.2 ml 50 µg/ml

Also, to supplement our existing inventory of perfluorinated compounds, we are introducing a chloro-perfluoro-1-octanesulfonate reference standard (**8Cl-PFOS**) and a branched perfluoroalkyl-carboxylic acid reference standard (**P37DMOA**).



8Cl-PFOS



P37DMOA

Catalogue Number	Product (methanol)	Qty/Conc
8Cl-PFOS	Sodium 8-chloroperfluoro-1-octanesulfonate	1.2 ml 50 µg/ml
P37DMOA	Perfluoro-3,7-dimethyloctanoic acid	1.2 ml 50 µg/ml

Distributed Throughout Europe and Middle East By-



Greyhound Chromatography & Allied Chemicals
6 Kelvin Park, Birkenhead,
Merseyside, CH41 1LT, U.K.
Tel: (+44)-0-151-649-4000 Fax: (+44)-0-151-649-4001
E-mail: info@greyhoundchrom.com
Web: www.greyhoundchrom.com

