



< [Vials & Syringes](#) < [Autosampler Syringes for Thermo Scientific/CTC/Leap GCs](#)

Autosampler Syringes for CTC: SGE

For syringe basics, including tips on handling & cleaning, needle termination codes, a needle gauge chart, needle point styles, and syringe termination information, [click here](#)

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

for CTC GCs



REQUEST ADDITIONAL INFORMATION

Description	Needle Length	Needle Gauge	Point Style	Needle term.	Vendor	Volume	Model #	qty.	Cat.#
Syringe	50mm	23	Cone	R	SGE	0.5µL	0.5BNR-C/F-0.63	ea.	24928
Syringe	50mm	26	Cone	R	SGE	0.5µL	0.5BNR-C/F-0.47	ea.	24926
Syringe	50mm	23	Cone	F	SGE	5µL	5F-C/F-5/0.63	ea.	24921
Syringe	50mm	26	Cone	F	SGE	5µL	5F-C/F-5/0.47	ea.	24920
Syringe	50mm	23	Cone	F	SGE	10µL	10F-C/F-0.63	ea.	24923
Syringe	50mm	23	Cone	R	SGE	10µL	10R-C/F-0.63	ea.	24932
Syringe	50mm	26	Cone	F	SGE	10µL	10F-C/F-5/0.47C	ea.	24922
Syringe	50mm	26	Cone	F	SGE	10µL	SK10F-C/F-5/0.47C	6-pk.	24925
Syringe	50mm	26	Cone	R	SGE	10µL	10R-C/F-5/0-0.47C	ea.	24930

Syringe	50mm	23	Cone	F	SGE	25µL	25F-C/T-GT-5/0.63C	ea.	22750
Syringe	50mm	26	Cone	F	SGE	25µL	25F-C/T-GT-5/0.47C	ea.	22750
Syringe	50mm	23	Cone	F	SGE	25µL*	25F-C/T-GT-5/0.63C	ea.	22751
Syringe	50mm	26	S-hole	R	SGE	25µL*	25R-C/T-GT-5/0.47H	ea.	22752
Syringe	50mm	23	Cone	F	SGE	100µL*	100F-C/T-GT-5/0.63C	ea.	22753
Syringe	50mm	26	Cone	R	SGE	100µL*	100R-C/T-GT-5/0.47C	ea.	22754

*Gas-tight syringe.

Click thumbnails to see enlarged view.



24923

Syringe, SGE, CE/CTC Autosampler, SGE 10-C/F-5/0.63 10ul/ASN-F/23/50mm/Cone



24922

Syringe, SGE, CE/CTC Autosampler, SGE 10-C/F-5/0.47C 10ul/ASN-F/26/50mm/Cone



22749

Syringe, SGE, 10F-C/T-GT-5/0.63C, Gas Tight, 10uL/23/50mm/cone

We accept Visa, MasterCard and American Express



Buy from an [Authorized Restek Distributor](#) to ensure top product quality, full warranty protection, and reliable technical support.

[HOME](#) | [NEWS](#) | [PRODUCTS](#) | [INFO. & SUPPORT](#) | [CUSTOMER SERVICE](#) | [CONTACT US](#)

Restek U.S., 110 Benner Circle, Bellefonte, PA 16823, 814-353-1300, 800-356-1688, fax: 814-353-1309