



GC Supplies - Thermo Scientific Instrument Quick Pick Guide

- Syringes
- Septa
- Inlet Liners
- Sealing Rings
- Ferrules
- SilTite™ FingerTite Ferrules
- Electron Multipliers



AUSTRALIA & PACIFIC REGION

CHINA

EUROPE

INDIA


JAPAN

MIDDLE EAST

USA

Autosampler Syringes for Thermo Scientific / CE Instruments / Fisons

All needles have a cone point style.



AUTOSAMPLER				NEEDLE SPECIFICATIONS						
TriPlus	AS3000	AS2000	AS200/800	Volume	Length (mm)	Gauge	Syringe Code	Syringe Part No.	Needle Part No.	Pack Size
Fixed Needle										
•				5 µL	50	26 (0.47)	5FX-5C	001105	-	-
•				10 µL	50	23 (0.63)	10F-C/T-GT-5/0.63C	002987	-	-
•		•		10 µL	80	23 (0.63)	10F-C/T-8/0.63C	002989	-	-
•		•		10 µL	80	26 (0.47)	10F-C/T-8/0.47C	002992	-	-
•	•	•	•	10 µL	50	25 (0.50)	10F-C/T-5/0.50C	002967	-	-
•	•	•	•	10 µL	50	23 (0.63)	10F-C/T-5/0.63C	002981	-	-
•	•	•	•	10 µL	50	26 (0.47)	10F-C/T-5/0.47C	002980	-	-
Removable Needle										
•	•	•		10 µL	50	26 (0.47)	10R-C/T-5/0.47C	002982	037785	2
•	•	•	•	10 µL	50	23 (0.63)	10R-C/T-5/0.63C	002984	037787	2
•		•	•	10 µL	80	26 (0.47)	10R-C/T-8/0.47C	002993	031535	3
•	•			0.5 µL	50	26 (0.47)	0.5BNR-C/T-5/0.47C	000490	033770	-
•	•			0.5 µL	50	23 (0.63)	0.5BNR-C/T-5/0.63C	000492	033772	-

Septa

Instrument	Diameter (mm)	Septum Type	Pack Size	Part No.
8000 Series, FOCUS, Trace/ULTRA GC™	17	TCS	50	0418491
	17	Auto-Sep T™	25	041886
	17	HT	25	0418990
	17	EC	25	041903

Inlet Liners

Description and Geometry Sketch	OD (mm)	ID (mm)	Length (mm)	Pack Size	Part no.
For TRACE™ and Focus™					
Splitless FocusLiner™ (for use with 70 mm Needle)	8	5	105	1	09204501
				5	092045
Splitless FocusLiner™ with Top-end Restriction (for use with 70 mm Needle)	8	5	105	1	09204601
				5	092046
				25	09204625
Split FocusLiner™ (for use with 50 mm Needle)	8	5	105	1	09204801
				5	092048
Splitless FocusLiner™ (for use with 50 mm Needle)	8	5	105	1	09204901
				5	092049
Splitless with Single Taper	8	3	105	1	09214101
				5	092141
				25	09214125
Splitless with Single Taper	8	5	105	1	09214401
				5	092144
				25	09214425
Splitless, Straight-through Liner	8	3	105	1	09214701
				5	092147
				25	09214725
SPME Liner	8	0.8	105	1	09214801
				5	092148
Split, Straight-through Liner	8	5	105	1	0921501
				5	092150
				25	09215025
Trace 2000 PTV Liner	2.75	1.75	120	1	09214201
				5	092142
				25	092154225
Sintered Glass, Large Volume Injection (LVI) Liner	2.75	0.78/1.8	120	1	09215501
				5	092155
				25	09215525
ConnectTite Liner Standard	8	5	105	1	09233401
				5	092334
ConnectTite Liner Top Hole	8	5	105	1	09233501
				5	092335
ConnectTite Liner Bottom Hole	8	5	105	1	09233601
				5	092336

- Taper / Gooseneck
- FocusLiner™
- Taper Focus
- ConnectTite
- Straight
- Double Taper
- PTV/LVI

Sealing Rings

Description	Usage	Pack Size	Part No.
Graphite Sealing Ring	8 mm ID. For 8000 & TRACE GCs.	10	0726004

Ferrules

Instrument	Column ID	Ferrule ID	Pack Size	Part No.
15% Graphite / 85% Vespel® Ferrules				
For All Injectors & Detectors (Not for 8000 Series, Focus, Trace/Ultra GC)	0.1-0.25 mm	0.4 mm	10	0726549
	0.32 mm	0.5 mm	10	0726557
	0.53 mm	0.8 mm	10	0726548
	1/8" OD Packed Columns	1/8"	10	072669
	1/4" OD Packed Columns	1/4"	10	072667
For All Focus, Trace/Ultra Injectors & Detectors at atmospheric pressure e.g. FID (Not for GC-MS)	0.1-0.25 mm	0.4 mm	10#	072696
	0.32 mm	0.5 mm	10#	072697
	0.53 mm	0.8 mm	10#	072698
Brass Nut for Focus, Trace/Ultra GC Injectors & Non-MS Detectors			2	1034085
100% Graphite Ferrules				
For Injectors & Detectors at atmospheric pressure e.g. FID (Not for GC-MS)	0.1-0.32 mm	0.5 mm	10	072619
	0.45-0.53 mm	0.8 mm	10	072614
	1/8" OD Packed Columns	1/8"	10	072622
	1/4" OD Packed Columns	1/4"	10	072621
SilTite™ Metal Ferrules				
For GC-MS Interface Connections (Starter Kit)	0.1-0.25 mm	0.4 mm	10*	073450
	0.32 mm	0.5 mm	10*	073451
	0.53 mm	0.8 mm	10*	073452
	1/32"	0.81 mm	10*	073453
Replacement SilTite™ Metal Ferrules				
For All GC-MS Interface Connections	0.1-0.25 mm	0.4 mm	10	073330
	0.32 mm	0.5 mm	10	073331
	0.53 mm	0.8 mm	10	073332
	1/32"	0.81 mm	10	073333
Replacement SilTite™ Nuts				
			5	073230

To be used in combination with brass nut (Part No. 1034085) * Includes 10 ferrules, 2 SilTite™ nuts.

SilTite™ FingerTite Ferrules

Description	Column ID	Ferrule ID (mm)	Pack Size	Part No.
Thermo (Available from SGE on 1/1/2011)				
SilTite™ FingerTite™ Thermo Injector/FID Starter Kit	0.1-0.25 mm	0.4	*	073614
SilTite™ FingerTite™ Thermo Injector/GC-MS ISQ Starter Kit	0.1-0.25 mm	0.4	*	073615
SilTite™ FingerTite™ Thermo Injector/EC-MS Starter Kit	0.1-0.25 mm	0.4	*	073616
SilTite™ FingerTite™ Ferrule 0.4 mm	0.1-0.25 mm	0.4	10	073630
SilTite™ FingerTite™ Ferrule 0.5 mm	0.32 mm	0.5	10	073631
SilTite™ FingerTite™ Ferrule Blanking	–	–	2	073633
SilTite™ FingerTite™ Female Nut	–	–	5	073636

* Each starter kit includes all the parts necessary to convert one GC system (one injector and one detector) to the SilTite™ FingerTite system suitable for installing 0.1 – 0.25 mm ID capillary columns. In addition there are 5 SilTite™ FingerTite nuts, one packet (10 ferrules) of 0.4 mm ID SilTite™ FingerTite ferrules and a ferrule install tool which allows you to seat the ferrule in the correct position on the capillary column.

Electron Multipliers

MS Description	Analyzer Type	Technique	Part No.
DSQ (Replacement multiplier)	Quadrupole	GC-MS	14633
Polaris-Q (Replacement multiplier only)	Ion Trap	GC-MS	14633
ITQ (Replacement multiplier only)	Ion Trap	GC-MS	14633
ISQ (Replacement multiplier only)	Quadrupole	GC-MS	14633

AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd
7 Argent Place, Ringwood
Victoria 3134 Australia
Toll Free: 1800 800 167
Tel: +61 (0) 3 9837 4200
Fax: +61 (0) 3 9874 5672
Email: support@sge.com

CHINA

SGE Shanghai Representative Office
Room 8221A, No. 808 Honqiao Road
Shanghai China
Tel: +86 21 6407 9382
Fax: +86 21 6407 9386
Email: china@sge.com

MIDDLE EAST

SGE Gulf
Y2 - 133, Sharjah Airport International
Free Zone (SAIF Zone)
P.O. Box 9377, Sharjah, UAE
Tel: +971 6 557 3341
Fax: +971 6 557 3541
Email: gulfsupport@sge.com

EUROPE

SGE Europe Ltd
European Head Office
1 Potters Lane
Kiln Farm Milton Keynes
MK11 3LA United Kingdom
Toll Free: 00800 2790 8999
Toll Free Fax: 00800 2626 2609
Tel: +44 1908 568 844
Fax: +44 1908 566 790
Tel France: +33 1 69 29 80 90
Fax France: +33 1 69 29 09 25
Tel Germany: +49 (0) 6155 / 60746 0
Fax Germany: +49 (0) 6155 / 60746 50
Email: europe@sge.com

INDIA

SGE Laboratory Accessories Pvt Ltd
430 Milan Industrial Estate
Abhudaya Nagar T. J. Road
Mumbai - 400 033 India
Tel: +91 22 24715896
Fax: +91 22 24716592
Email: sgeindia@vsnl.com

UNITED STATES OF AMERICA

SGE Incorporated
2007 Kramer Lane
Austin Texas 78758 USA
Toll Free: (800) 945 6154
Tel: +1 512 837 7190
Fax: +1 512 836 9159
Email: usa@sge.com

JAPAN

SGE Japan Inc
RK CUBE 3F 6-85, Ohta-machi
Naka-Ku, Yokohama-shi
Kanagawa 231-0011 Japan
Tel: +81 45 222 2885
Fax: +81 45 222 2887
Email: japan@sge.com



SGE's involvement in Gas Chromatography includes the support we give to our customers.

Technical Support

Contact us at techsupport@sge.com for technical support information.

Our Website

www.sge.com is a detailed source for information on chromatography, applications and SGE products. It also profiles the latest posters and technical materials generated by our technical experts based around the world.

50 Years Edition Product Selection Guide

Our product selection guide can be accessed at www.sge.com/selectionguide

SGE Approved Distribution Partner